visocolor® ECO Chlorine dioxide

en

Test kit for performing colorimetric tests on chlorine dioxide in drinking water, water reservoirs and disinfectant solutions

Method:

At a pH value of 5 to 6, chlorine dioxide reacts with *N,N*-diethyl-1,4-phenylenediamine (DPD) and forms a red-violet dye.

Measurement range:

0.2-3.8 mg/L CIO₂

Contents of test kit (*refill pack):

sufficient for 150 tests

16 mL CIO₂-1* 18 mL CIO₂-2* 25 mL CIO₂-3* 2 screw-plug measuring glasses

slide comparator

color chart

1 plastic syringe 5 mL 1 instructions for use*

Hazard warning:

Information regarding safety can be found on the box' label and in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Instructions for use:

also refer to the pictogram on the back of the color chart

Pour 5 mL water sample into both of the measuring glasses using the plastic syringe.
Place one of them on position A in the comparator.

- Only add the reagent to measuring glass B.

 2. Fill the second measuring glass with 2 drops of CIO₂-1.
- 3. Seal the glass and mix.
- 4. Open the glass after 2 min and add 3 drops of CIO2-2, seal the glass and mix
- 5. Add 3 drops of CIO₂-3, seal the glass and mix.
- Open the glass once again and place it on position B in the comparator. 6. Open the glass once again and place it on position 2 in the inspection hole on Slide the comparator until the colors match in the inspection hole on the reason of the recess on the
- top. Immediately check the measurement reading in the recess on comparator reed. Mid-values can be estimated. 8 After use, rinse out both measuring glasses thoroughly and seal them.
- The reagents can also be used for the **photometric evaluation** with photometers PF-12, PF-12^{Plus} and PF-3.

The method cannot be applied for the analysis of sea water.

Disposing of the samples:

Information regarding disposal can be found in the can download the SDS from www.mn-net.com/SDS. the safety data sheet. You

Interferences:

Free chlorine up to 5 mg/L is not determined with this procedure and thus, does not interfere.

Storage:

Store the test kit in a cool (< 25 °C) and dry place.