
visocolor[®] alpha

Gesamthärte

total Hardness / Dureté totale /
Durezza totale / Durezza total /
totaal Hardheid

REF 935042

100 Tests / 10 °d

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Disposing of the samples: Information regarding disposal can be found in the safety data sheet. You can download the SDS from www.mn-net.com/SDS.

Instructions for use: also refer to the pictogram on the inside

1. Rinse the test vessel several times with the water sample and fill to the ring mark (5 mL).
2. Hold the dropping bottle H₇ absolutely vertical and add the reagent drop by drop while smoothly shaking the test vessel until the colour turns from red to green. Count the number of drops. One drop corresponds to one degree of total water hardness (°d).

This technique can also be applied for analysing sea water dilution (1+29).

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Test kit for performing titrimetric tests on total hardness

Hazard warning: H₇ contains 2-propanol 15–20%.

Information regarding safety: you can download the SDS from www.mn-net.com/SDS.

Instructions for use: also refer to the pictogram on the inside

1. Rinse the test vessel several times with the water sample and fill to the ring mark (5 mL).
2. Hold the dropping bottle H₇ absolutely vertical and add the reagent drop by drop while smoothly shaking the test vessel until the colour turns from red to green. Count the number of drops. One drop corresponds to one degree of total water hardness (°d).

This technique can also be applied for analysing sea water dilution (1+29).

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