
Inhalt:

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
en

Contents:

- 1 electric air pump
- 1 lid with bore-hole
- 2 aerating bricks
- 4 Winkler oxygen flasks, vol. 250–300 mL

- 1 storage vessel, vol. 10 L
- 1 flexible tube, length 1 m
- 1 laboratory bottle, vol. 1 L

1. Set-up of BOD₅-accessories:

- Insert flexible tube through the bore-hole of the lid, and connect the end to the pump.
- Fill the storage vessel with water (see diluting water).
- Close the lid of the storage vessel firmly.
- Be sure to position the pump always above the water level in the vessel, otherwise water may flow back and damage the pump. **The pump should remain on a dust-free, clean surface.**
- Plug in to the mains.
- Switch the air pump off by unplugging it from the mains, before taking diluting water for BOD₅-measurements. Then open the storage vessel, and extract the required volume of diluting water.

2. Diluting water:

- Unchlorinated tap water is used to produce diluting water. The water must not contain more than 0.01 mg copper, chlorine or chloramine. Its pH-Value should be between 6 and 8 and must be re-corrected as necessary.
- The diluting water is **aerated at $(20 \pm 2) ^\circ\text{C}$ in the dark for at least 6 hours**. Make sure the concentration of dissolved oxygen in the diluting water is at least 8 mg/L O₂.
- The water must not be oversaturated with oxygen. Therefore, leave to stand in an open container for about 1 hour **before** measurement. Use this solution within 24 hours and discard the remains

The added laboratory flask in the BOD₅-Accessories Set will be used for the preparation and mixture of the sample dilutions, the Winkler oxygen flasks for the incubation of prepared sample dilutions within the determination of BOD₅.

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