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

NANOCOLOR® BOD₅-TT Accessories Set

Contents:

- | | |
|----------------------|--------------------------------|
| 1 electric air pump | 1 storage vessel, vol. 1000 mL |
| 1 lid with bore-hole | 1 flexible tube, length 1 m |
| 2 aerating bricks | 2 reaction vessels, vol. 40 mL |

1. Set-up of BOD₅-TT accessories

- Insert flexible tube through the bore-hole of the lid, and connect the end to the pump.
- Fill the storage vessel with water (see diluting water).
- Close the lid of the storage vessel firmly.
- Be sure to position the pump always above the water level in the vessel, otherwise water may flow back and damage the pump.
The pump should remain on a dust-free, clean surface.
- Plug in to the mains.
- Switch the air pump off by unplugging it from the mains, before taking water for dilution for BOD₅ measurements. Then open the storage vessel, and extract the required volume of diluting water.

2. Diluting water

- Unchlorinated tap water is used to produce diluting water. The water must not contain more than 0.01 mg/L copper, chlorine or chloramine. The pH-value should be between 6 and 8 and must be recorrected if necessary.
- The diluting water is **aerated at (20 ± 2) °C in the dark** for at **least 1 hour**. Make sure the concentration of dissolved oxygen in the diluting water is at least 8 mg/L O₂.
- The water must not be oversaturated with oxygen. Therefore, leave to stand in an open container for about 1 hour before measurement. Use this solution within 24 hours and discard the remains.

3. Inoculating water and inoculated diluting water

If the sample does not contain enough suitable microorganisms, then so-called **inoculating water** is added to the diluting water, which can be obtained in the following ways:

- Biomass from a sewage plant from preliminary sedimentation channel
- Biomass from a sewage plant from drainage channel with a maximum of 300 mg/L COD or 100 mg/L TOC
- Surface water which contains municipal sewage

Producing inoculated diluting water:

Add 5 mL inoculating water to 500 mL aerated diluting water immediately before use. Mix well and store at +20 °C.

At the end of the day, discard not-used inoculated diluting water. Oxygen consumption in the inoculated diluting water, which forms the control mixture, must not exceed 1.5 mg/L O₂ after 5 days at +20 °C.

Examples for samples, which require the use of inoculated diluting water:

- Industrial, biologically inactive sewage
- Surface water without any municipal sewage
- Frozen samples!
- Control test with glucose-glutamic acid standard solution in accordance with DDIN EN 1899-1-H51 (**NANOCONTROL BOD₅**, REF 92582)!

The added reaction vessels in the BOD₅-TT Accessories Set will be used for the preparation of the control and sample dilutions. For further informations see instructions for use of the test kit BOD₅-TT - Test 8-25 (REF 985825).

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