



HI931

Automatic Potentiometric Titrator

Dear Customer,

Thank you for choosing a Hanna Instruments[®] product.

This manual has been written for HI931 Automatic Potentiometric Titrator with software version 1.03 and higher.

Please read this instruction manual carefully before using this instrument as it provides the necessary information for correct use of this instrument and a precise idea of its versatility.

If you need additional technical information, do not hesitate to e-mail us at tech@hannainst.com.

Visit www.hannainst.com for more information about Hanna Instruments and our products.

*All rights are reserved. Reproduction in whole or in part is prohibited without the copyright owner's written consent, Hanna Instruments Inc., Woonsocket, Rhode Island, 02895, USA.
Hanna Instruments reserves the right to modify the design, construction, or appearance of its products without advance notice.*

INTRODUCTION

The **HI931** is an automatic potentiometric titrator with high accuracy, great flexibility and repeatability.

The titrator is designed to perform a variety of potentiometric titrations, allowing the user to obtain both good results and high-speed analysis.

The main attributes of the **HI931** titrator are:

- Small footprint, requires minimal bench space
- Casing made with strong, chemically resistant plastic
- Flexible electrode holder supports up to 3 electrodes, 4 dispensing tubes, 1 temperature sensor and removable stirrer
- Electrode holder positions electrodes in the center of beaker, allowing for smaller sample sizes
- Support for 100 titration methods
- User-customizable reports
- Integrated research grade pH/mV/ISE meter
- Clearly displayed warning and error messages

This manual provides information regarding installation and functionality of the titrator and refined operation suggestions.

Before using the titrator, it is recommended you become familiar with its various features and functionality.

This manual is divided into four parts:

PART 1: Quick Start Guide

Helps the user quickly setup and operate **HI931** Automatic Potentiometric titrator. It covers basic connections, user interface and how to run a titration.

PART 2: Instruction Manual

Provides a comprehensive description of the operating principles, user interface, general options, methods, titration mode, optimization, maintenance etc.

PART 3: Applications

Contains complete instructions for commonly-used analyses. Additional methods and method packs are available, contact your local Hanna Instruments office for more details.

PART 4: Titration Theory

Outlines the principles of operation of the titrator. It covers the chemistry of titrations, titration types and result calculations

TABLE OF CONTENTS

PART 1. QUICK START GUIDE

1.1. Safety Measures	1-1	1.9. The First Titration	1-5
1.2. Abbreviations	1-1	1.9.1. Required Solutions.....	1-5
1.3. Titrator Connections	1-2	1.9.2. Priming the Burette.....	1-5
1.3.1. Front View.....	1-2	1.9.3. Method Selection	1-5
1.3.2. Rear View.....	1-2	1.9.4. Setting Method Parameters.....	1-6
1.4. User Interface	1-3	1.9.5. Setting Up Titration Report.....	1-6
1.4.1. Keypad	1-3	1.9.6. Preparing the Sample.....	1-6
1.4.2. Display	1-3	1.9.7. Performing a Titration.....	1-7
1.5. Language	1-4	1.9.8. Titration Screen.....	1-7
1.6. Contextual Help	1-4	1.9.9. Titration Graph.....	1-7
1.7. Methods	1-4	1.9.10. Titration Termination	1-7
1.7.1. Standard Methods.....	1-4	1.9.11. Results.....	1-8
1.7.2. User-Defined Methods.....	1-4	1.9.12. Viewing the Last Titration Data.....	1-8
1.8. How to Calibrate a pH Electrode	1-4	1.9.13. Printing the Titration Report	1-8
1.8.1. Preparation	1-4	1.9.14. Saving Data to USB Storage Device	1-9
1.8.2. Calibration Procedure	1-5	1.9.15. Titration Report.....	1-9

PART 2. INSTRUCTION MANUAL

2.1. Setup	2-1	2.3.9. Language.....	2-20
2.1.1. Unpacking.....	2-1	2.3.10. Total Volume Alert.....	2-21
2.1.2. Safety Measures	2-2	2.3.11. Titrant Age Reminder.....	2-21
2.1.3. HI931 Titrator Technical Specifications.....	2-2	2.3.12. USB Link with PC.....	2-22
2.1.4. Installation.....	2-4	2.3.13. Setup Balance Interface.....	2-22
2.2. User Interface	2-9	2.3.14. Printer Mode.....	2-25
2.2.1. Start Up.....	2-9	2.3.15. Reset to Default Settings.....	2-26
2.2.2. Keypad	2-9	2.3.16. Optimize Memory Space	2-26
2.2.3. Display	2-11	2.3.17. Update Software.....	2-26
2.2.4. Menu Navigation	2-12	2.4. Titration Methods	2-27
2.3. General Options	2-13	2.4.1. Selecting Methods.....	2-27
2.3.1. Save to USB.....	2-13	2.4.2. Standard Methods.....	2-28
2.3.2. Restore from USB.....	2-14	2.4.3. User-Defined Methods.....	2-29
2.3.3. Administration.....	2-15	2.4.4. Viewing / Modifying Method	2-30
2.3.4. Temperature	2-16	2.4.5. Method Options	2-31
2.3.5. Date & Time Setting	2-18	2.4.6. Printing.....	2-54
2.3.6. Display Settings.....	2-19	2.5. Titration Mode	2-54
2.3.7. Beeper	2-19	2.5.1. Running a Titration	2-54
2.3.8. Stirrer	2-20	2.5.2. Stopping a Titration	2-55

2.6. pH Mode	2-56	2.8.3. ISE Calibration.....	2-79
2.6.1. Display	2-56	2.8.4. Logging.....	2-80
2.6.2. pH Setup.....	2-57	2.9. Auxiliary Functions	2-80
2.6.3. pH Calibration	2-63	2.9.1. Burette.....	2-80
2.6.4. Logging.....	2-65	2.9.2. Stirrer	2-82
2.7. mV Mode	2-65	2.9.3. Results.....	2-83
2.7.1. Display	2-66	2.10. Maintenance & Peripherals	2-86
2.7.2. mV Setup.....	2-66	2.10.1. Burette Maintenance	2-86
2.7.3. Relative mV Calibration	2-68	2.10.2. Peripherals	2-89
2.7.4. Logging.....	2-69	2.11. Accessories	2-92
2.8. ISE Mode	2-69	2.11.1. Solutions.....	2-92
2.8.1. Display	2-70	2.11.2. Sensors.....	2-95
2.8.2. ISE Setup.....	2-70	2.11.3. Titrator Components	2-97

PART 3. APPLICATIONS

HI0001EN – 0.1N Sodium Hydroxide Titrant Concentration	3-1	HI1005EN – Acidity of Water.....	3-13
HI0002EN – 0.1N Hydrochloric Acid Titrant Concentration.....	3-3	HI1007EN – Chloride in Water.....	3-15
HI0003EN – 0.1M Sodium Thiosulfate Titrant Concentration	3-5	HI1008EN – Neutralization with Sulfuric Acid	3-17
HI0010EN – 0.1M Ferrous Ammonium Sulfate Titrant Concentration	3-7	HI1009EN – Neutralization with Sodium Hydroxide.....	3-19
HI0200EN – 0.02M Silver Nitrate Titrant Concentration	3-9	HI1011EN – Troubleshooting 1.....	3-21
HI1004EN – Alkalinity of Water.....	3-11	HI1012EN – Troubleshooting 2.....	3-23

PART 4. TITRATION THEORY

4.1. Titration Theory	4-1	4.4.3. Sources of Error.....	4-10
4.1.1. Introduction.....	4-1	4.5. Calculations	4-11
4.1.2. Uses of Titrations	4-1	4.5.1. Sample Calculation by Mass.....	4-11
4.1.3. Advantages & Disadvantages.....	4-1	4.5.2. Sample Calculation by Volume	4-11
4.2. Types of Titrations	4-2	4.5.3. Standardize Titrant by Mass.....	4-12
4.2.1. Titrations According to the Measurement Method	4-2	4.5.4. Standardize Titrant by Volume	4-12
4.2.2. Titrations According to the Reaction Type.....	4-3	4.5.5. Blank Titration.....	4-12
4.2.3. Titrations According to the Titration Sequence.....	4-8	4.5.6. Multiple Endpoint Titration.....	4-13
4.3. Titration Procedure	4-9	4.6. Glossary	4-14
4.3.1. Manual Titration	4-9	4.7. List of Figures	4-17
4.3.2. Automatic Titration	4-9	Certification	4-18
4.4. Titration Results	4-10	Recommendations for Users	4-18
4.4.1. Accuracy.....	4-10	Warranty	4-18
4.4.2. Repeatability.....	4-10		

PART 1. QUICK START GUIDE

1.1. SAFETY MEASURES

The following safety measures must be followed:


1. Never connect or disconnect the pump assembly or other peripheral with the titrator turned on.
2. Verify that the burette and the attached tubing are assembled correctly.
3. Always check that the titrant bottle and the titration beaker are placed on a flat, stable surface.
4. Always wipe up spills and splashes immediately.
5. Avoid the following environmental working conditions
 - Severe vibrations
 - Direct sunlight
 - Atmospheric relative humidity above 95% non-condensing
 - Environment temperatures below 10 °C and above 40 °C
 - Explosion hazards
6. Have the titrator serviced by qualified service personnel only.

1.2. ABBREVIATIONS

ABS	Acrylonitrile Butadiene Styrene	PTFE	Polytetrafluoroethylene
GLP	Good Laboratory Practice	PVDF	Polyvinylidene fluoride
PEI	Polyetherimide	RPM	Revolutions per minute
eq / kg	Equivalentents per kilogram	mmol / kg	Millimoles per kilogram
eq / L	Equivalentents per liter	mmol / L	Millimoles per liter
g / 100 mL	Grams per 100 milliliters	M (mol / L)	Molarity (moles per liter)
g / L	Grams per liter	mol / kg	Moles per kilogram
μg / L	Micrograms per liter	mol / L	Moles per liter
meq / kg	Milliequivalentents per kilogram	N (eq / L)	Normality (equivalentents per liter)
meq / L	Milliequivalentents per liter	ppb (μg / kg)	Parts per billion (micrograms per kilogram)
mg / 100 mL	Milligrams per 100 milliliters	ppb (μg / L)	Parts per billion (micrograms per liter)
mg / g	Milligrams per gram	ppm (mg / kg)	Parts per million (milligrams per kilogram)
mg / kg	Milligrams per kilogram	ppm (mg / L)	Parts per million (milligrams per liter)
mg / L	Milligrams per liter	ppt (g / kg)	Parts per thousand (grams per kilogram)
mmol / g	Millimoles per gram	ppt (g / L)	Parts per thousand (grams per liter)
% (g / 100 g)	Percent by weight (grams per 100 grams)		
%w / v	Percent weight by volume		


1.3. TITRATOR CONNECTIONS

1.3.1. FRONT VIEW



1. Titrant bottle
2. Aspiration tube (Titrant inlet)
3. Pump assembly
4. Burette assembly
5. Dispensing tube (Titrant outlet)
6. Support rod
7. Stirrer
8. Temperature sensor
9. pH electrode
10. Stirrer holder
11. Dispensing tip
12. Titration beaker

1.3.2. REAR VIEW




1. Burette blank support
2. On/Off button
3. Extension interface connector
4. Balance interface connector
5. Printer connector
6. PC keyboard connector
7. Power adapter connector
8. PC interface
9. Analog board

2026 Titrator image shown above.

1.4. USER INTERFACE


1.4.1. KEYPAD

The titrator's keypad has 27 keys grouped in five categories, as follows:



1.4.2. DISPLAY

The titrator has a 5.7" graphical backlit color display.



- | | |
|------------------------------------|-------------------------------|
| 1. Time and Date | 7. Results units |
| 2. Method name | 8. Stirrer information |
| 3. Result | 9. End point titration volume |
| 4. Reminders or Warnings | 10. mV or pH reading |
| 5. Temperature compensation status | 11. Virtual option keys |
| 6. Temperature reading | |


The user interface contains several screens. In each screen, many information fields are present at the same time. The information is displayed in an easy-to-read manner.

Virtual option keys describe the function performed when the corresponding soft key is pressed.

1.5. LANGUAGE

To change the language, press **General Options** from the main screen. Highlight *Language* option. Using the **▲** and **▼** keys, select the language from the options listed in the **Set Language** screen and press **Select**.

Restart the titrator in order to apply the new language setting.



1.6. CONTEXTUAL HELP

Information about the titrator can be easily accessed by pressing **?**. The contextual help can be accessed at any time and it provides useful information about the current screen.

1.7. METHODS

The HI931 titrator can store up to 100 methods (standard and user-defined).

1.7.1. STANDARD METHODS

Each titrator is supplied with a package of standard methods. Standard method packs are developed at Hanna Instruments® to meet analysis requirements of specific industries (e.g., water treatment, wine, dairy, etc.).

1.7.2. USER-DEFINED METHODS

User defined methods allow the user to create and save their own methods. Each new method is based on an existing method which is altered to suit a specific application.

1.8. HOW TO CALIBRATE A pH ELECTRODE

To enter pH calibration mode, press **Mode**, then **pH**, then **pH Calibr.**.

1.8.1. PREPARATION

Pour small quantities of the buffer solutions into clean beakers. If possible, use plastic beakers to minimize any EMC interferences. For accurate calibration and to minimize cross-contamination, use two beakers for each buffer solution: one for rinsing the electrode and one for calibration. If you are measuring in the acidic range, use pH 7.01 or 6.86 as the first buffer and pH 4.01/3.00 or 1.68 as the second buffer. If you are measuring in the alkaline range, use pH 7.01 or 6.86 as the first buffer and pH 10.01/9.18 or 12.45 as the second buffer.



For extended range measurements (acidic and alkaline), perform a five-point calibration by selecting five buffers across the entire pH range.

1.8.2. CALIBRATION PROCEDURE





During calibration, the user has a choice of 8 standard buffers: (pH 1.68, 3.00, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45) and up to 5 custom buffers.

For accurate measurements it is recommended to perform a five-point calibration. However, at least a two-point calibration is suggested. For pH titrations, the selected buffers should bracket your endpoint (e.g. if your endpoint value is at 8.5, use 7.01 or 6.86 and 9.18 or 10.01 for calibration).

To begin calibration:

1. Press . If the instrument was calibrated before, previous calibration can be cleared by pressing .

Note: It is very important to clear calibration history when a new electrode is used.

2. Immerse the pH electrode and the temperature probe approximately 4 cm (1.5") into a buffer solution and stir gently.
3. If necessary, select the pH calibration buffer value with  or .
4. Once the reading has stabilized, press  to update the calibration. The calibration buffer will be added to the Calibrated Buffers section.
5. Rinse the pH electrode and the temperature probe, then immerse them into the next buffer solution and follow the above procedure or press  to exit the calibration.




1.9. THE FIRST TITRATION

1.9.1. REQUIRED SOLUTIONS

- Titrant - 500 mL of 0.1 M (mol/L) Sodium Hydroxide (NaOH) in a titrant bottle.
- Sample - 0.1 mol/L Hydrochloric Acid (HCl).
- Distilled or deionized water.

Note: Analytical grade reagents and water should be used for accurate results.

1.9.2. PRIMING THE BURETTE

1. Insert the aspiration tube in the titrant bottle and the dispensing tube in a waste beaker.
2. From the main screen press .
3. Highlight the *Prime Burette* option and then press .
4. Enter the number of burette rinses. At least 3 rinses are recommended.
5. Press  to start. The message "Executing..." will be displayed.

Note: Make sure you have continuous liquid flow inside the burette. For accurate results, the aspiration tube, the dispensing tube and the syringe must be free of air bubbles.

1.9.3. METHOD SELECTION

For this analysis we will use the **HI1009 Neutralization w/ NaOH** standard method.

To select this method:

1. Press  from the Idle screen.
2. Use the  and  keys to highlight **HI1009 Neutralization w/ NaOH** method.
3. Press .

1.9.4. SETTING METHOD PARAMETERS

To display the method parameters, press Method Options.

The **View/Modify Method** screen will be displayed.

Only certain parameters can be changed.

For this titration, the NaOH titrant concentration and the size of the HCl sample need to be entered.

1. Highlight *Titrant Conc.* option, then press Select. The Titrant Concentration screen will be displayed.
2. Enter the correct value, then press Accept.
3. Highlight *Analyte Size* option, then press Select.
4. Enter the volume of the sample (e.g.: 5 mL), then press Accept.
5. Press Escape, highlight *Save Method* option and then press Select.

Titrant Concentration

Enter the titrant concentration.

0.10123 M (mol/L)

ACCEPT	Escape	Delete Digit		Exponent
--------	--------	--------------	--	----------

1.9.5. SETTING UP TITRATION REPORT

Users can select the information that is stored for each titration.

To obtain proper information at the end of the titration, perform the following operations:

1. From the main screen, press results and the **Data Parameters** screen will be displayed.
2. Highlight *Setup Titration Report* option and press Select.
3. Mark the fields to be included with the * symbol using the ▲ and ▼ keys, and press Select to toggle the selection.
4. Press Save Report and then press Escape to return to the main screen.

1.9.6. PREPARING THE SAMPLE

1. Add 50 to 65 mL of distilled / deionized water to the titration beaker.
2. Use a pipette or burette to add 5.0 mL of the sample (0.1M Hydrochloric Acid (HCl)) into the same beaker.
3. Slide the stirrer assembly up.
4. Place the beaker under the stirrer assembly.
5. Lower the stirrer assembly until the electrodes are submersed and the stirrer is close to the bottom of the beaker.
6. Adjust the level of the sample solution with distilled / deionized water so that the pH electrode bulb is completely immersed in the sample solution and the reference junction of the electrode is 5-6 mm below the surface.


1.9.7. PERFORMING A TITRATION

From the main screen, press **start stop**. You will be prompted to enter the analyte size. Enter 5 mL and press **enter**. The titrator will start the analysis.

At the end of the titration, the message "Titration Completed" will appear on the display with the final concentration of the analyte in the sample and the equivalence endpoint volume.


1.9.8. TITRATION SCREEN

During a titration, the following screen is displayed:



1.9.9. TITRATION GRAPH

After a few doses are dispensed, **View Curve** will become active. Press **View Curve** to display the real-time titration graph. The curves displayed are plots of the pH and the 1st derivative versus Titrant Volume. See [PART 2. INSTRUCTION MANUAL](#) for more information. The two graphs are scaled to fit in the same screen window. Press **Select Trace** to change the y-axis scale to either the pH values or the 1st derivative values.



1.9.10. TITRATION TERMINATION


The titration is terminated when the conditions of the Termination Criteria have been met.

The titration is normally terminated when the first equivalence endpoint is detected according to the selected algorithm. To ensure the correct detection and interpolation of the equivalence endpoint, the titrator will dispense a few additional doses after the endpoint was reached. The titration result can be displayed either in the main screen or in the Graph of Titration Data screen.

When the titration has ended, the titrator will display the equivalence endpoint volume and the final concentration of the analyte together with the "Titration Completed" message.


To view the titration graph and/or results, press **results**.

When the titration ends, an "x" will mark the endpoint on the pH versus titrant volume curve in the Graph of Titration Data screen. The value of the endpoint volume is also displayed next to the endpoint



1.9.11. RESULTS

The results obtained from titration are stored in a report file that can be displayed, transferred to a USB storage device or a PC, or printed.



1.9.12. VIEWING THE LAST TITRATION DATA

To view the last titration report:

1. From the main screen, press **results**. The **Data Parameters** screen will be displayed.
2. From the **Data Parameters** screen highlight *Review Last Report* option and press **Select**. The **Review Result** screen will be displayed.
3. Use the **Page Up** and **Page Down** keys to display information related to the last titration performed.

1.9.13. PRINTING THE TITRATION REPORT

Connect a DOS / Windows-compatible parallel printer directly to the DB 25 pin connector located on the back of the titrator.

Note: Prior to connecting the printer, ensure that the titrator and the printer have been turned off.

From the **Review Report** screen, press **Print Report**. During the information transfer to the printer, the message "Printing" will be displayed on the screen.

Press **Escape** to return to the **Data Parameters** screen.


Press **Escape** again to return to the main screen.

1.9.14. SAVING DATA TO USB STORAGE DEVICE

Note: The USB Storage Device has to be formatted FAT or FAT32.

This feature allows saving the results of titrations logging sessions on a USB storage device.

1. From the main screen, press General Options, the **General Options** screen will be displayed.
2. Highlight *Save Files to USB Storage Device* option using the \triangle and ∇ keys.
3. Insert the USB storage device into the USB socket.
4. Press Select, the **List of Files on Titrator** screen will be displayed.
5. Use the \triangleleft and \triangleright keys to select the report files.



6. Press Copy All to transfer all available reports to the USB storage device, or highlight the name of the report file to be transferred and press Copy File. Transferring a report file will automatically transfer the corresponding log file and titration graph (*.BMP file if applicable).
7. Press Escape to return to the **General Options** screen.
8. Press Escape again to return to the main screen.

1.9.15. TITRATION REPORT

While scrolling with the Page Up and Page Down keys, the fields below can be seen on the titrator display or printed. The same information is available on the saved report file (Ti_00011.rpt in this example, with all report fields selected).

HI931 - Titration Report

Method Name: Neutralization w/ NaOH
 Time & Date: 15:01 Jun 13, 2018
 Report ID: Ti_00011

Calibration Data

Buffer	Potential	Efficiency	Temp.
Time and Date			
4.010pH	169.3mV	98.8%	24.0°C A
			11:44 Jun 13, 2018
7.010pH	-5.8mV	98.7%	23.9°C A
			11:42 Jun 13, 2018
10.010pH	-180.7mV	98.7%	24.0°C A
			11:46 Jun 13, 2018

GLP & Meter Information

Sample Name:
Company Name:
Operator Name:
Electrode Name:
Field 1:
Field 2:
Field 3:
Titrator Software Version: v1.00
Base Board Software Version: v1.00
Pump 1 Software Version: v1.00
Stirrer 1 Software Version: v1.00
Titrator Serial Number: TT180525011
Analog Board1 Serial Number: AB180525005
Pump 1 Serial Number: DP180525004
Stirrer 1 Serial Number: OS180524001
Analog 1 Calibration Date: May 25, 2018

Method Parameters

Name: Neutralization w/ NaOH
Method Revision: 3.0
Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
Pump Configuration:
 Titrant pump: Pump 1
Dosing Type: Dynamic
 Min Vol: 0.050 mL
 Max Vol: 0.500 mL
 delta E: 20.000 mV
End Point Mode: pH 1EQ point, 1st Der
Recognition Options
 Threshold: 50 mV/mL
 Range: NO
 Filtered Derivatives: NO
Pre-Titration Volume: 0.000 mL
Pre-Titration Stir Time: 0 sec
Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 15 sec
Electrode Type: pH
Blank Option: No Blank
Calculations: Sample Calc. by Volume
Dilution Option: Disabled
Titrant Name: 0.1N NaOH
Titrant Conc.: 0.1000 N (eq/L)
Analyte Size: 10.0000 mL
Analyte Entry: Fixed
Maximum Titrant Volume: 20.000 mL
Potential Range: -2000.0 to 2000.0 mV
Volume/Flow Rate: 25 mL / 50.0 mL/min
Signal Averaging: 1 Reading
Significant Figures: XXXXX

N (eq/L) --> meq/L

V eq 1000meq

--*-----

L eq

mL L

--*-----

1000mL

V = volume dispensed in liters.

0.100 eq/L -> titrant conc.

10.000 mL -> sample volume

Nr	Volume[mL]	mV	pH	Graphic	Temp. [°C]	Time
0	0.000	274.4	2.219	0.0	24.9	A 00:00:00
1	0.050	274.4	2.220	1.0	25.0	A 00:00:07
2	0.100	274.4	2.220	0.0	25.0	A 00:00:10
3	0.200	274.3	2.222	-0.8	25.0	A 00:00:12
4	0.400	274.0	2.227	-1.6	25.0	A 00:00:15
5	0.800	273.2	2.241	-2.0	25.0	A 00:00:18
6	1.300	271.5	2.271	-3.4	25.0	A 00:00:24
7	1.800	269.5	2.304	-3.9	25.1	A 00:00:30
8	2.300	267.2	2.344	-4.7	25.1	A 00:00:37
9	2.800	264.4	2.393	-5.7	25.1	A 00:00:43
10	3.300	260.8	2.455	-7.2	25.1	A 00:00:50
11	3.800	256.1	2.535	-9.3	25.1	A 00:00:58
12	4.300	250.3	2.635	-11.7	25.1	A 00:01:05
13	4.800	241.9	2.779	-16.8	25.1	A 00:01:14
14	5.300	228.3	3.011	-27.2	25.1	A 00:01:23
15	5.800	193.0	3.614	-70.5	25.1	A 00:01:31
16	6.077	21.0	6.556	-620.0	25.1	A 00:01:48
17	6.128	-38.2	7.568	-1183.2	25.1	A 00:02:03
18	6.177	-123.6	9.031	-1708.0	25.1	A 00:02:19
19	6.227	-157.7	9.616	-682.8	25.1	A 00:02:28
20	6.278	-174.5	9.903	-335.8	25.1	A 00:02:35
21	6.339	-187.8	10.130	-215.9	25.1	A 00:02:42

Titration Results

Method Name: Neutralization w/ NaOH
 Time & Date: 15:01 Jun 13, 2018
 Analyte Size: 10.0000 mL
 End Point Volume: 6.144 mL
 pH Equivalence Point: 8.063
 Result: 61.444 meq/L
 Initial & Final pH: 2.219 to 10.130
 Titration Duration: 2:42 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

PART 2. INSTRUCTION MANUAL

2.1. SETUP

2.1.1. UNPACKING

Remove the titrator and accessories from the packaging and examine it carefully. For further assistance, please contact your local Hanna Instruments® office or email us at tech@hannainst.com.

Each HI931 potentiometric titrator is supplied with:

- Titrator
- Pump assembly
- Burette assembly
 - » Burette with 25 mL syringe
 - » Aspiration tube with fitting and protection tube
 - » Dispensing tube with dispensing tip, protection tube and tube guide
 - » Tube locks
 - » Tool for burette cap removal
 - » Light spectrum protection screen
- Electrodes holder and stirrer
 - » Stirrer holder
 - » Overhead stirrer
 - » Propellers (3 pcs.)
 - » Support rod
- Burette blank support
- Pump and burette locking screws with plastic head
- Temperature sensor
- Shorting cap
- Power adapter
- USB cable
- USB memory stick
- HI900 PC application (installation kit on USB memory stick)
- Quality certificate
- Instruction manual

If any of the items are missing or damaged, please contact your local Hanna Instruments office or email us at tech@hannainst.com.

See [2.11.3. Titrator Components](#) section for component pictures.

Note: Save all packing materials until you are sure that the instrument functions correctly. Any damaged or defective items must be returned in their original packing materials together with the supplied accessories.

2.1.2. SAFETY MEASURES

The following safety measures must be followed:

1. Never connect or disconnect the pump assembly or other peripheral with the titrator turned on.
2. Verify that the burette and the attached tubing are assembled correctly.
3. Always check that the titrant bottle and the titration beaker are placed on a flat, stable surface.
4. Always wipe up spills and splashes immediately.
5. Avoid the following environmental working conditions
 - Severe vibrations
 - Direct sunlight
 - Atmospheric relative humidity above 95% non-condensing
 - Environment temperatures below 10 °C and above 40 °C
 - Explosion hazards
6. Have the titrator serviced by qualified service personnel only.


2.1.3. HI931 TITRATOR TECHNICAL SPECIFICATIONS

Analysis Type	Standard titration (Standardization, Fixed pH / mV, Equivalence point pH / mV)		
Endpoint Mode	Fixed mV		
	Fixed pH		
	mV equivalence point (1 st or 2 nd derivative)		
	pH equivalence point (1 st or 2 nd derivative)		
Burette	Size	5 mL / 10 mL / 25 mL / 50 mL	
	Resolution	0.001 mL	
	Flow Rate	0.3 mL to 2 x burette volume per minute	
	Accuracy	±0.005 mL (5 mL burette)	
		±0.010 mL (10 mL burette)	
±0.025 mL (25 mL burette)			
±0.050 mL (50 mL burette)			
Stirrer	Range	200 to 2500 RPM	
	Resolution	100 RPM	
mV	Range	–2000.0 to 2000.0 mV	
	Resolution	0.1 mV	
	Accuracy	±0.1 mV	
	Calibration	Single point, offset	
pH	Range	–2.000 to 20.000 pH	
	Resolution	0.1 / 0.01 / 0.001 pH	
	Accuracy	±0.001 pH	
	Calibration	Up to five points with standard or custom buffers	


ISE	Range	1×10^{-6} to 9.999×10^{10}
	Resolution	1 / 0.1 / 0.01
	Accuracy	± 0.001 pH
	Calibration	Up to five points
Temperature	Range	–5.0 to 105 °C 23.0 to 221.0 °F 268.2 to 378.2 K
	Resolution	0.1 °C / 0.1 °F / 0.1 K
	Accuracy	± 0.1 °C / ± 0.2 °F / ± 0.1 K
Data Storage	Methods	Up to 100 titration methods (standard and user-defined)
	Reports	Up to 100 titration and pH / mV / ISE reports
Connections	Measurement	1 × BNC socket (pH, ORP, ISE half-cell and ISE combination electrodes) 1 × 4 mm banana socket (reference electrode) 1 × RCA socket (temperature sensor) 1 × 6-pin connector (stirrer)
	Peripheral	1 × 6-pin mini DIN (external PC keyboard) 1 × DB-25 socket (printer) 1 × USB standard B (PC connection) 1 × DB-9 socket (analytical balance) 1 × USB standard A (USB flash drive)
Additional Specifications	Electrode Holder	4 × multi-purpose slots (titrant tubes) 3 × 12-mm electrode slots 1 × temperature sensor slot 1 × overhead stirrer slot
	Display	5.7" graphical color display with backlight
	Power Supply	100 - 240 VAC, 50 / 60 Hz
	Power Draw	0.5 amps
	Enclosure Material	PC-ABS, Steel
	Keypad	Polyester
	Dimensions	315 × 205 × 375 mm (12.4 × 8.1 × 14.8 ")
	Weight	Approximately 4.3 kg (9.5 lbs.) with 1 pump, stirrer and sensors
	Operating Environment	10 to 40 °C (50 to 104 °F); up to 95 % RH
	Storage Environment	–20 to 70 °C (–4 to 158 °F); up to 95 % RH

2.1.4. INSTALLATION

2.1.4.1. Titrator Front View




2.1.4.2. Titrator Rear View



2026 Titrator image shown above.

2.1.4.3. Titrator Right-Side View




2.1.4.4. Titrator Assembly

Note: Assembly operations must be completed before connecting the titrator to the power supply!

Assembling Support Rod

To insert support rod into the titrator case:

1. Remove protective cap from titrator case
2. Insert the support rods into the titrator case
3. Turn the support rod clockwise to secure it in place




Attaching Stirrer & Electrode

To attach the stirrer to the titrator, follow these steps:

1. Place the electrode holder (B) in the stirrer support housing (A). The stirrer support housing can be inverted if necessary.
2. Slide the electrode holder into the support rod and set the desired height using the holder button.
3. Insert electrode (C), temperature sensor (D) and stirrer (E) into the dedicated holes in the electrode holder. Push them until they are in stable position.

Note: For small sample sizes, use the electrode adapter (F) in the center of the holder.




Connecting the Pump

To connect the pump, follow these steps:

1. Retrieve the pump cable from inside the bay. The pump 1 connector is located in the left bay and pump 2 connector in the right bay.
2. Connect the cable to the pump as shown below. The pump connector is located on the bottom of the pump.
3. Lower the pump into the titrator, then slide it towards the front of the titrator case until it is firmly latched.
4. Secure the pump with the locking screw.

This procedure can be repeated to connect a second pump.



Attaching Burette Blank Support


To attach the support, follow these steps:

1. Insert and lower the burette blank support into the titrator's bay.
2. Slide it towards the front of the titrator case until it is firmly latched.
3. Secure the burette blank support with the locking screw.


Attaching the Burette

To attach the burette to the pump, follow these steps:

1. Make sure that the mark from the valve actuating cap and from the burette body are aligned.




2. Slide the burette into the support on the burette pump. Ensure correct coupling between the syringe plunger (A) and the pump piston (B).



2.1.4.5. Electrical Connections

1. Connect the electrode to the BNC connector (A).
2. Connect the temperature sensor to the RCA connector (C).
3. Connect the stirrer to the MINI-DIN connector (D).
4. Connect the power adapter cable to the power input connector (F).



2026 Titrator image shown above.

Function


Type of Connector

A. Connection for pH, ORP, ISE half-cell and ISE combination electrodes	BNC socket
B. Reference electrode	Ø 4 mm banana socket
C. Temperature sensor	RCA socket
D. Stirrer	6-pin connector
E. USB interface	USB standard B
F. Power input connector (24 VDC)	DC power jack connector
G. External PC keyboard	6-pin mini DIN (Standard PS2)
H. Printer	DB-25 socket
I. Extension	5-pin connector
J. Balance interface	DB-9 socket (RS-232)
K. Power switch	

2.2. USER INTERFACE

2.2.1. START UP

Once the instrument is assembled and installed, follow the steps below to start the titrator:


1. Connect the titrator to a power outlet with the supplied power adapter.
2. Turn on the titrator using the power button located on the back of the instrument.
3. Wait until the titrator finishes the initialization process.
4. Press  when prompted or wait a few seconds for titrator to start.



Note: All the performed initialization processes must be successfully completed. If one fails, restart the titrator. If the problem persists, contact your nearest Hanna Instruments® Service Center.






2.2.2. KEYPAD

The titrator's keypad is grouped into five categories, as follows:



2.2.2.1. Function Keys

If one of these keys is pressed, the associated function is immediately performed. Some of the keys are active only in specific screens:

	Starts or stops a titration process
	Turns the selected stirrer on and off
	Reserved
	Access the data parameters menu (reports, GLP, meter information, report setup)
	Displays contextual help

2.2.2.2. Option Keys

These keys are assigned to the virtual keys on the display. Their functions are listed in the boxes above the buttons and vary depending on the displayed screen.




An underlined virtual key can also be activated by pressing .

2.2.2.3. Arrow Keys

These keys have the following functions:

- Move the on-screen cursor.
- Increase or decrease the stirrer speed and other settings.
- Select a character (alphanumeric screen only).
- Navigate through menu options.

2.2.2.4. Numeric Keys

-  Used for numeric entries.
-  Toggles between positive and negative values.
-  Used for decimal point.

2.2.2.5. Enter Key


This key has the following functions:

- Accepts alphanumeric data entry.
- Executes the default (underlined) virtual option key.

2.2.3. DISPLAY

The titrator has a large color graphical display. The main screen is shown below with short explanations of the screen segments.

2.2.3.1. The Main Screen




1. Time and Date
2. Method name
3. Result
4. Reminders or Warnings
5. Temperature compensation status
6. Temperature reading
7. Results units
8. Stirrer information
9. End point titration volume
10. mV or pH reading
11. Virtual option keys

The user interface contains several screens. For each titrator function, several screens may be used.

Method name	Displays the name of the selected method.
Time and Date	Displays the current date and time.
Temperature reading	Displays the measured temperature.
ATC	Automatic temperature compensation
Manual	Manual temperature compensation
Manual	Temperature probe is not connected, manual temperature compensation
Stirrer information	The selected stirrer, actual and set stirrer speed is displayed in RPM. When stirrer is off, the stirrer information is not displayed.
Endpoint volume	Displays the volume delivered to reach the titration endpoint. When no titration has been performed, the displayed volume is "0.000 mL".
Result	Displays the titration result.
mV or pH reading	Displays the current reading. The reading will be in mV or pH.
mV	Indicates actual potential reading.
rel mV	Indicates relative potential reading.
pH	Indicates actual pH value.
Titration status	Displays the status of the selected titration.
	No Result is displayed when a titration has not been performed.
Reminders	Indicates when a task needs to be performed and displays errors.
Pump 1 selected	Displays the active pump.

2.2.4. MENU NAVIGATION

2.2.4.1. Selecting an Option



Press the option key below the virtual key. For example, to access the **Method Options** screen, press the option key below it.

2.2.4.2. Selecting a Menu Item




Use the Δ and ∇ arrow keys to move the cursor. When the menu is larger than the display, a scroll bar is active on the right side.

The \uparrow and \downarrow keys can be used to scroll through the pages.

To activate the selected menu item, press **enter** or **Select**.

2.2.4.3. Entering Text



Use **Delete Letter** to erase previous text. Use the arrow keys to highlight the letter then press **enter**. Use the same procedure to enter the whole name.


For editing, use the **Cursor Left** and **Cursor Right** keys.

When editing is complete, press **Accept**.

The method name will be updated and displayed in the name field of the **View / Modify Method** screen.

When all the desired parameters have been set, press **Escape**.

2.2.4.4. Saving Modifications




The **Saving Method** screen allows the user to save the modifications. To exit without saving, press **Escape** or highlight *Exit Without Saving Method* option and then press **Select**. To save the modifications, highlight *Save Method* option and then press **Select**.

Note: To access the contextual help menu, press ? at any time. Help is related to the displayed screen. Press Escape or ? to return to the previous screen.

2.3. GENERAL OPTIONS


The **General Options** screen gives access to options that are not directly related to the titration process or pH / mV / ISE measurement. To access this screen, press General Options from the main screen.



2.3.1. SAVE TO USB

This option allows the user to save files from the titrator to a USB storage device.

Note: The USB Storage Device has to be formatted FAT or FAT32.



On the titrator, the available file types are:

Standard method HIXXXYY.MTD (e.g.: HI0001EN.MTD, HI1004EN.MTD)

User-defined method USERXXX.MTD (e.g.: USER0001.MTD)

Report Ti_XXXX.RPT, mV_XXXX.RPT, pH_XXXX.RPT, ISEXXXX.RPT, mVrXXXX.RPT (e.g.: Ti_00001.RPT, mV_00001.RPT, pH_00001.RPT, ISE00001.RPT, mVr00001.RPT)

Insert the USB storage device into the USB port on the right side of the titrator.

Use the ← and → keys to select the file type. The number of files and the file names will be displayed.

Use the ↑ and ↓ keys to scroll through the list.

The option keys allow the following operations:

- Escape Returns to the **General Options** screen.
- Copy File Copies highlighted file from the titrator to USB storage device.
- Copy All Copies all displayed files from the titrator to USB storage device.
- Delete File Deletes the highlighted file.
- Delete All Deletes all displayed files.


Note: The saved files will be stored on the USB key in the **HI931** folder, as follows:

Methods USB Drive\HI931\Methods*.mtd

Reports USB Drive\HI931\Reports*.rpt

2.3.2. RESTORE FROM USB

This screen allows the user to transfer files from the USB storage device to the titrator.



The file types that can be transferred are:

Standard method HIXXXYY.MTD (e.g.: HI0001EN.MTD, HI1004EN.MTD)

User-defined method USERXXX.MTD (e.g.: USER0001.MTD)

Report Ti_XXXX.RPT, mV_XXXX.RPT, pH_XXXX.RPT, ISEXXXX.RPT, mVrXXXX.RPT (e.g.: Ti_00001.RPT, mV_00001.RPT, pH_00001.RPT, ISE00001.RPT, mVr00001.RPT)

Insert the USB storage device into the USB port on the right side of the titrator.

Use the < and > keys to select the file type. The number of files and the file names will be displayed.

Use the ^ and v keys to scroll through the list.

The option keys allow the following operations:

- Escape Returns to the **General Options** screen.
- Copy File Copies the highlighted file from the USB storage to the titrator.
- Copy All Copies all displayed files from the USB storage to the titrator.
- Delete File Deletes the highlighted files from the USB storage device.
- Delete All Deletes all displayed files from the USB storage device.

Note: In order to restore files from USB Key, please ensure that the methods and / or reports you wish to transfer to the titrator are in the correct folder:

Methods USB Drive\HI931\Methods*.mtd

Reports USB Drive\HI931\Reports*.rpt

2.3.3. ADMINISTRATION

A 4-digit numeric PIN can be set to prevent unauthorized changes from being made.

When the user enters administration mode and a pin has not been set, the user will be prompted to enter a new PIN.

Titrator Administration				
Administrator PIN has not been set. Enter a 4-digit PIN to enable Administrator function.				
Enter PIN: <input type="password" value="----"/>				
Confirm PIN: <input type="password" value="----"/>				
Your PIN must be 4-digits long.				
Next	Escape	Delete Digit		

Once a PIN has been set, the titrator can be locked. When a titrator is locked, the users cannot modify methods or delete reports. Basic functions are still available (review reports, save to USB, etc.).

Titrator Administration				
Titrator is UNLOCKED.				
<div style="border: 1px solid black; padding: 5px; width: fit-content; margin: auto;"> Lock Titrator Enter PIN: <input type="password" value="----"/> </div>				
Accept	Escape	Delete Digit		

To return to administration mode, the titrator can be unlocked by entering the PIN.

Titrator Administration				
Titrator is LOCKED.				
Unlock Titrator	Escape			Recovery PIN

If the PIN is lost or forgotten, press recovery pin and contact technical support to supply the required information.

Recovery PIN				
For recovery PIN, please contact your vendor. When requesting PIN please provide following information:				
Titrator Serial Number: 12345678				
Code: 0078				
Recovery PIN: [REDACTED]				
Accept	Escape	Delete Digit		

2.3.4. TEMPERATURE

The temperature menu allows access to all three menu options related to temperature: source, setting and units.

Temperature Menu							
Select temperature option to be modified.							
<table border="1"> <tbody> <tr> <td>Temperature Source</td> </tr> <tr> <td>Manual Temperature Setting</td> </tr> <tr> <td>Temperature Units</td> </tr> </tbody> </table>					Temperature Source	Manual Temperature Setting	Temperature Units
Temperature Source							
Manual Temperature Setting							
Temperature Units							
Select	Escape						

2.3.4.1. Temperature Source


Option: Automatic Temperature or Manual Temperature

Select the temperature source used for temperature compensation.

When Automatic Temperature is selected, "ATC" is displayed on the main screen and the temperature is read by the temperature probe.

When Manual Temperature is selected, "Manual" is displayed on the main screen and a preset temperature value is used for temperature compensation.


Note: The selected temperature source will be indicated in the report files: A for Automatic and M for Manual.



2.3.4.2. Manual Temperature Setting

Option: -5.0 to 105.0 °C (23.0 to 221.0 °F, 268.2 to 378.2 K)

If the temperature probe is not connected, the user can manually set the temperature used by the titrator for compensation.



2.3.4.3. Temperature Units

Option: °C, °F, K


The temperature ranges are as displayed in the **Temperature Units** screen.

Temperature Menu													
Select temperature option to be modified.													
Temperature Source Manual Temperature Setting Temperature Units													
<table border="1"> <tbody> <tr> <td>Celsius</td> <td>-5.0 to 105.0</td> <td>°C</td> </tr> <tr> <td>Fahrenheit</td> <td>23.0 to 221.0</td> <td>°F</td> </tr> <tr> <td>Kelvin</td> <td>268.2 to 378.2</td> <td>K</td> </tr> </tbody> </table>					Celsius	-5.0 to 105.0	°C	Fahrenheit	23.0 to 221.0	°F	Kelvin	268.2 to 378.2	K
Celsius	-5.0 to 105.0	°C											
Fahrenheit	23.0 to 221.0	°F											
Kelvin	268.2 to 378.2	K											
Select	Escape												


2.3.5. DATE & TIME SETTING

This screen allows the user to set the date and time.

Use the  and  keys or the numeric keys to modify the date and time.

 Moves move the cursor to the next field.

 or  Changes the time format.

Date and Time Setting				
Enter the date.				
	2	10	2018	
day	month	year		
Enter the time.				
20	41	41		
hour	minute	second		
Press <Next> to move to the next entry.				
Accept	Escape	Delete Digit	Next	AM/PM

2.3.6. DISPLAY SETTINGS


This screen allows the user to customize the display settings.

 Increases the backlight time-saver interval.

 Decreases the backlight time-saver interval.

The backlight intensity can be adjusted using the  and  keys.

There are 8 levels of backlight intensity, ranging from 0 to 7.



The displayed color palette allows for selection of appropriate backlight intensity.

The backlight time-saver option protects the display during standby periods, when no keys have been pressed for a set amount of time, the backlight will turn off.


If the backlight is off, press any key to reactivate the backlight.

The range for backlight time-saver interval is between 1 and 60 minutes. To disable the backlight time-saver, increase the time to the maximum allowed, the Off indication will be displayed.

2.3.7. BEEPER

Option: On or Off


If enabled (on) an audible alert will sound after a titration is completed, when an invalid key is pressed or when a critical error occurs during titration.



2.3.8. STIRRER

Option: Enabled or Disabled


The stirrer can be disabled in individual titration method, if necessary.



2.3.9. LANGUAGE

Using the \triangle and ∇ keys, select the language from the options listed and press Select.

Restart the titrator in order to apply the new language setting.



2.3.10. TOTAL VOLUME ALERT

Option: Off, 0 to 10000 mL

This screen allows a programmable reminder to appear when the titrant reservoir is below 100 mL. The titrant volume will decrease as the titrant is used.

After the new titrant volume has been entered in the **Total Volume Alert** screen, a warning message appears on the main screen reminding the user to re-standardize the newly added titrant.

Off Disables this option.

Total Volume Alert				
Enter the amount of titrant available to the titration/reagent system from its reservoir. The mLs will decrease as the titrant/reagent is depleted.				
1000 mL				
A reminder will appear when less than 100 mLs of titrant volume is left.				
ACCEPT	Escape	Delete Digit		Off

2.3.11. TITRANT AGE REMINDER

Option: Off, 0 to 31 days


A programmable reminder will appear when it is time to verify the titrant concentration or to change the titrant.

Off Disables this option.

Titrant Age Reminder				
Enter the number of days to pass since the last Titr. Vol. updating or the last Start pressing, whereafter the reminder appears.				
30 days				
The range is from 0 to 31 days.				
Start	Escape	Delete Digit		Off

2.3.12. USB LINK WITH PC

In order to use this feature, the USB cable needs to connect the titrator with the PC. Make sure that HI900 PC application is running on the PC.



“Active / Inactive” message shows the status of the USB link with the PC.


“Active” means that the titrator is using the USB communication with the PC and not with another device.

“Ready” shows that the titrator is able to communicate with the PC.

During transfer of information between the PC and the titrator, press “Transmit” and the status is displayed.

2.3.13. SETUP BALANCE INTERFACE

This screen allows the user to setup an analytical balance for automatic acquisition of sample mass prior to titration or standardization.



The balance is connected to the titrator via RS 232 interface.

- Enable Balance Enables the selected balance.
- Disable Balance Disables the selected balance (automatic weight acquisition will not be available).
- Escape Returns to the **General Options** screen.
- New Balance Adds a new balance to the list.
- Edit Customizes the serial communication parameters. The **Balance Configuration** screen will open.
- Delete Deletes the highlighted balance.

Note: At least one balance must be in the list.

Be sure that the balance configuration settings match the settings of your balance. It may be necessary to change settings on your balance or titrator. Users should consult their balance instruction manual.

Verify connection with the balance is working properly by pressing the  key.

Balance Configuration																
Select the option to be modified.																
<table border="1"> <thead> <tr> <th>Balance Name</th> <th>Default</th> </tr> </thead> <tbody> <tr> <td>Baud Rate</td> <td>9600</td> </tr> <tr> <td>Data Bits</td> <td>8 Bits</td> </tr> <tr> <td>Parity</td> <td>No Parity</td> </tr> <tr> <td>Stop Bit</td> <td>1 bit</td> </tr> <tr> <td>Request Command</td> <td>B</td> </tr> </tbody> </table>					Balance Name	Default	Baud Rate	9600	Data Bits	8 Bits	Parity	No Parity	Stop Bit	1 bit	Request Command	B
Balance Name	Default															
Baud Rate	9600															
Data Bits	8 Bits															
Parity	No Parity															
Stop Bit	1 bit															
Request Command	B															
Select	Escape		Test Balance													

2.3.13.1. Balance Name

Option: Up to 24 characters

Assign a name for your customized balance.

Balance Name																																																																																																																																					
Select the highlighted letter by using the arrow keys then press "Enter". Select the empty field for a space. Press Accept to save the entered text.																																																																																																																																					
<table border="1"> <tbody> <tr><td>A</td><td>B</td><td>C</td><td>D</td><td>E</td><td>F</td><td>G</td><td>H</td><td>I</td><td>J</td><td>K</td><td>L</td></tr> <tr><td>M</td><td>N</td><td>O</td><td>P</td><td>Q</td><td>R</td><td>S</td><td>T</td><td>U</td><td>V</td><td>W</td><td>X</td><td>Y</td></tr> <tr><td>Z</td><td>a</td><td>b</td><td>c</td><td>d</td><td>e</td><td>f</td><td>g</td><td>h</td><td>i</td><td>j</td><td>k</td><td>l</td></tr> <tr><td>m</td><td>n</td><td>o</td><td>p</td><td>q</td><td>r</td><td>s</td><td>t</td><td>u</td><td>v</td><td>w</td><td>x</td><td>y</td></tr> <tr><td>z</td><td>À</td><td>Á</td><td>Â</td><td>Ã</td><td>Ä</td><td>Å</td><td>Æ</td><td>Ç</td><td>È</td><td>É</td><td>Ê</td><td>Ë</td></tr> <tr><td>Ì</td><td>Í</td><td>Î</td><td>Ï</td><td>Ñ</td><td>Ò</td><td>Ó</td><td>Ô</td><td>Õ</td><td>Ö</td><td>×</td><td>÷</td><td>ä</td></tr> <tr><td>å</td><td>ç</td><td>è</td><td>é</td><td>ê</td><td>ë</td><td>ì</td><td>í</td><td>î</td><td>ï</td><td>ñ</td><td>ò</td><td>ó</td></tr> <tr><td>ô</td><td>õ</td><td>ö</td><td>÷</td><td>×</td><td>÷</td><td>×</td><td>÷</td><td>×</td><td>÷</td><td>×</td><td>÷</td><td>×</td></tr> <tr><td>0</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td><td>9</td><td>×</td><td>÷</td><td>,</td></tr> <tr><td>? !</td><td>()</td><td>[]</td><td>< ></td><td>=</td><td>/</td><td>+</td><td>-</td><td></td><td></td><td></td><td></td><td></td></tr> </tbody> </table>					A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	a	b	c	d	e	f	g	h	i	j	k	l	m	n	o	p	q	r	s	t	u	v	w	x	y	z	À	Á	Â	Ã	Ä	Å	Æ	Ç	È	É	Ê	Ë	Ì	Í	Î	Ï	Ñ	Ò	Ó	Ô	Õ	Ö	×	÷	ä	å	ç	è	é	ê	ë	ì	í	î	ï	ñ	ò	ó	ô	õ	ö	÷	×	÷	×	÷	×	÷	×	÷	×	0	1	2	3	4	5	6	7	8	9	×	÷	,	? !	()	[]	< >	=	/	+	-					
A	B	C	D	E	F	G	H	I	J	K	L																																																																																																																										
M	N	O	P	Q	R	S	T	U	V	W	X	Y																																																																																																																									
Z	a	b	c	d	e	f	g	h	i	j	k	l																																																																																																																									
m	n	o	p	q	r	s	t	u	v	w	x	y																																																																																																																									
z	À	Á	Â	Ã	Ä	Å	Æ	Ç	È	É	Ê	Ë																																																																																																																									
Ì	Í	Î	Ï	Ñ	Ò	Ó	Ô	Õ	Ö	×	÷	ä																																																																																																																									
å	ç	è	é	ê	ë	ì	í	î	ï	ñ	ò	ó																																																																																																																									
ô	õ	ö	÷	×	÷	×	÷	×	÷	×	÷	×																																																																																																																									
0	1	2	3	4	5	6	7	8	9	×	÷	,																																																																																																																									
? !	()	[]	< >	=	/	+	-																																																																																																																														
Lab Balance																																																																																																																																					
Accept	Escape	Delete Letter	Cursor Left	Cursor Right																																																																																																																																	

2.3.13.2. Baud Rate

Option: 4800, 9600, 19200, 38400


Set the serial communication baud rate.

Balance Configuration																
Select the option to be modified.																
<table border="1"> <tbody> <tr> <td>Balance Name:</td> <td>Lab Balance</td> </tr> <tr> <td>Baud Rate:</td> <td>9600</td> </tr> <tr> <td>Data Bit:</td> <td></td> </tr> <tr> <td>Parity:</td> <td>N</td> </tr> <tr> <td>Stop Bit:</td> <td>4800</td> </tr> <tr> <td>Request Command:</td> <td>9600</td> </tr> </tbody> </table>					Balance Name:	Lab Balance	Baud Rate:	9600	Data Bit:		Parity:	N	Stop Bit:	4800	Request Command:	9600
Balance Name:	Lab Balance															
Baud Rate:	9600															
Data Bit:																
Parity:	N															
Stop Bit:	4800															
Request Command:	9600															
Select	Escape		Test Balance													

2.3.13.3. Data Bits

Option: 5, 6, 7, 8 bits


Set the number of data bits.



2.3.13.4. Parity

Option: No Parity, Even, Odd


Set the parity of data packet.



2.3.13.5. Stop Bit

Option: 1 bit or 2 bits


Set the number of stop bits.



2.3.13.6. Edit Request Stop


Option: Up to 10 characters

Type the syntax for weight request command.



2.3.14. PRINTER MODE

Option: Ansi, Ascii, Text



Ansi Use this mode when the printer is set as Ansi. When in this mode, all available accented characters and symbols will be printed.

Ascii Use this mode when the printer is set as Ascii. When in this mode, only some of the available accented characters and symbols will be printed.

Text This mode is recommended when the user doesn't need to print accented characters.

2.3.15. RESET TO DEFAULT SETTINGS

Note: This will delete all user-defined methods and restore all manufacturer settings such as titrator configuration, standard method parameters, etc.

Confirmation of Reset				
Are you sure you want to reset the titrator to manufacturer settings?				
This will delete the calibration data, all user methods, balances and reports.				
Reset	Escape			

2.3.16. OPTIMIZE MEMORY SPACE

This screen allows the user to run a memory defragmentation utility to increase access speed to memory storage. Press Accept and then restart the titrator. Do not disconnect the power supply during this operation.

Optimize Memory Space				
This option is used in order to clean up the memory space.				
Please ensure the power is not disconnected during this operation.				
Accept	Escape			

2.3.17. UPDATE SOFTWARE

This screen allows the user to update the titrator software from a USB storage device containing a software setup kit.

Update Software				
Current version:		HI931 v1.00		
New version:		HI931 v1.01		
Are you sure you want to update the current software with the new version?				
Accept	Escape	Refresh		

To update the software:

1. Copy the "Setup931" folder to a USB storage device.
2. Insert the USB storage device into the USB port.
3. Go to **General Options**, then **Update Software**. The titrator will display the current and new software versions.
4. Press . When prompted, remove the USB storage device and restart the titrator.

2.4. TITRATION METHODS


All parameters required to complete an analysis are grouped into a method.

The titrator is supplied with a pack of standard methods, these methods have been developed by Hanna Instruments® and can be used to create user-defined methods.


Standard and user-defined methods can be upgraded, saved or deleted by connecting the titrator to a PC using the HI900 PC application or a USB flash drive.


2.4.1. SELECTING METHODS

To select a method, press  from the main screen. A list of available methods will be displayed.



In the **Analysis Methods** screen, you can view the list of all available methods (standard and user-defined methods).

To select a method, highlight the method then press . The name of the selected method will be displayed on the main screen.



2.4.2. STANDARD METHODS

The standard methods are developed for the most common types of analysis and can be used as templates to create new user-defined methods.

Only specific method parameters can be modified by the use. See [2.4.5. Method Options](#) for more information.

2.4.2.1. Upgrading Standard Methods

To upgrade the titrator with new standard methods, follow the steps below:

From USB storage device:

1. Insert the USB storage device into the USB port, located on the right side of the titrator.
2. Press **General Options** from the main screen.
3. Using **▲** and **▼** keys, highlight *Restore from USB Storage Device* option and choose **Select**.
4. Using **<** and **>** keys, navigate through file types to find “standard method files”.
5. Press the **Copy File** or **Copy All** key to upgrade the titrator with the standard methods.
6. Press **Escape** to return to **General Options** screen.

From PC:

You can upgrade the titrator with standard methods from a PC using the **HI900** PC application. See [2.3.12. USB Link with PC](#) for more information.

2.4.2.2. Deleting Standard Methods

Standard methods can be removed from the titrator by following one of the procedures below.

From General Options Screen:


1. Using the **▲** and **▼** keys, highlight *Save to USB* option and press **Select**.
2. Using the **<** and **>** keys, navigate through the file types menu to find the list of “standard method files”.
3. Press the **Delete** or **Delete All** keys to remove unnecessary standard methods.
4. Press **Escape** to return to the **General Options** screen.

From PC:

The not required standard methods can be removed from the titrator using the **HI900** PC application. See [2.3.12. USB Link with PC](#) for more information.

2.4.2.3. Restoring the Standard Methods to the Manufacturer Settings

You can restore the standard methods to the default settings by highlighting a standard method and pressing **Reset to Default**.









2.4.3. USER-DEFINED METHODS

User-defined methods are created by users, by modifying a standard method or previously created user-defined method. All method parameters can be modified to suit user-specific requirements.

2.4.3.1. Creating User-Defined Methods





To create a new user-defined method, start from a standard or previously generated user-defined method and follow these steps:


1. Press  from the main screen.
2. Using the  and  keys, highlight an existing method from the method list.
3. Press . A new user-defined method will be generated.
4. Press  to activate the new method.



Note: The titrator can hold 100 methods (standard and user-defined). When the limit is reached, a warning message is displayed.

2.4.3.2. Deleting User-Defined Methods

1. To remove a user-defined method, press  from the main screen.
2. Highlight the user-defined method that you want to delete and press , a confirmation screen will appear.
3. Press  again to confirm, or press  to cancel the operation.



2.4.4. VIEWING / MODIFYING METHOD

To modify the method parameters, press Method Options from the main screen. A list of all the parameters for the selected method will be displayed. Press the ^ and v keys to highlight the option you want to modify and choose Select.

View/Modify Method																																
Id: HI0001EN Modified: 12:04 Sep 12, 2018																																
Select the option to be modified.																																
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border: 1px solid black;">Name:</td> <td style="border: 1px solid black;">0.1N Sodium Hydroxide</td> </tr> <tr> <td style="border: 1px solid black;">Method Revision:</td> <td style="border: 1px solid black;">3.0</td> </tr> <tr> <td style="border: 1px solid black;">Stirrer Configuration</td> <td style="border: 1px solid black;"></td> </tr> <tr> <td style="border: 1px solid black;">Titrant pump:</td> <td style="border: 1px solid black;">Pump 1</td> </tr> <tr> <td style="border: 1px solid black;">Dosing Type:</td> <td style="border: 1px solid black;">Dynamic</td> </tr> <tr> <td style="border: 1px solid black;">End Point Mode:</td> <td style="border: 1px solid black;">pH 1EQ point,1st Der</td> </tr> <tr> <td colspan="2" style="border: 1px solid black;">Recognition Options</td> </tr> <tr> <td style="border: 1px solid black;">Pre-Titration Volume:</td> <td style="border: 1px solid black;">5.000 mL</td> </tr> <tr> <td style="border: 1px solid black;">Pre-Titration Stir Time:</td> <td style="border: 1px solid black;">60 sec</td> </tr> <tr> <td style="border: 1px solid black;">Measurement Mode:</td> <td style="border: 1px solid black;">Signal Stability</td> </tr> <tr> <td style="border: 1px solid black;">Electrode Type:</td> <td style="border: 1px solid black;">pH</td> </tr> <tr> <td style="border: 1px solid black;">Blank Option:</td> <td style="border: 1px solid black;">No Blank</td> </tr> <tr> <td style="border: 1px solid black;">Calculations:</td> <td style="border: 1px solid black;">Stdz. Titrant by Weight</td> </tr> <tr> <td style="border: 1px solid black;">Dilution Option:</td> <td style="border: 1px solid black;">Disabled</td> </tr> </table>					Name:	0.1N Sodium Hydroxide	Method Revision:	3.0	Stirrer Configuration		Titrant pump:	Pump 1	Dosing Type:	Dynamic	End Point Mode:	pH 1EQ point,1st Der	Recognition Options		Pre-Titration Volume:	5.000 mL	Pre-Titration Stir Time:	60 sec	Measurement Mode:	Signal Stability	Electrode Type:	pH	Blank Option:	No Blank	Calculations:	Stdz. Titrant by Weight	Dilution Option:	Disabled
Name:	0.1N Sodium Hydroxide																															
Method Revision:	3.0																															
Stirrer Configuration																																
Titrant pump:	Pump 1																															
Dosing Type:	Dynamic																															
End Point Mode:	pH 1EQ point,1st Der																															
Recognition Options																																
Pre-Titration Volume:	5.000 mL																															
Pre-Titration Stir Time:	60 sec																															
Measurement Mode:	Signal Stability																															
Electrode Type:	pH																															
Blank Option:	No Blank																															
Calculations:	Stdz. Titrant by Weight																															
Dilution Option:	Disabled																															
Escape	Print Method	Page Up	Page Down																													

To exit the **View / Modify Method** screen, press the Escape key and highlight *Save Method* or *Exit Without Saving Method*.

Saving Method						
Select a menu option.						
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border: 1px solid black;">Save Method</td> </tr> <tr> <td style="border: 1px solid black;">Exit Without Saving Method</td> </tr> </table>					Save Method	Exit Without Saving Method
Save Method						
Exit Without Saving Method						
"Escape" - exits without saving method.						
Select	Escape					

Select Saves modifications.

Escape Discards the changes.

2.4.5. METHOD OPTIONS

Note: Not all method options can be changed for standard methods.

2.4.5.1. Name

Option: Up to 24 characters

Method Name				
Select the highlighted letter by using the arrow keys then press "Enter". Select the empty field for a space. Press Accept to save the entire name.				
<pre> A B C D E F G H I J K L M N O P Q R S T U V W X Y Z a b c d e f g h i j k l m n o p q r s t u v w x y Z Å Ä Å Å Ä Ç È É Ê Ë Ì Í Î Ï Ð Ñ Ò Ó Ô Õ Ö Ø Ù Ú Û Ü Ý Þ à á â ã ä å ç è é ê ë ì í î ï ð ó ô õ ö ø ù ú û ü ý ÿ ÿ \ / \$ % ^ & # : 0 1 2 3 4 5 6 7 8 9 * , . ? ! () [] < > = / + - </pre>				
copy of 0.1N Sodium Hydr				
Accept	Escape	Delete Letter	Cursor Left	Cursor Right

2.4.5.2. Method Revision

Option: Up to 3 characters

Method Revision				
Select the highlighted letter by using the arrow keys then press <Enter>. Select the empty field for a space. The revision string format is "X.X".				
<pre> B C D E F G H I J K L M N O P Q R S T U V W X Y Z a b c d e f g h i j k l m n o p q r s t u v w x y z Å Ä Å Å Ç È É Ê Ë Ì Í Î Ï Ð Ñ Ò Ó Ô Õ Ö Ø Ù Ú Û Ü Ý Þ à á â ã ä å ç è é ê ë ì í î ï ð ó ô õ ö ø ù ú û ü ý ÿ ÿ \ / \$ % ^ & # : 0 1 2 3 4 5 6 7 8 9 * , . ? ! () [] < > = / + - * / \ _ & ^ ' : </pre>				
1.0				
Accept	Escape	Delete Letter	Cursor Left	Cursor Right


2.4.5.3. Stirrer Configuration

Use the arrow keys to select the menu option.

Stirrer Configuration				
Select a menu option.				
<pre> Stirrer: Stirrer 1 Stirring Speed: 1400 RPM </pre>				
Select	Escape			


Stirrer

Option: Stirrer 1 or Disabled



Stirrer Speed


Option: 200 to 2500 RPM



The stirrer will remain on for as long as the method is active. When the stirrer is running, the speed can be adjusted at any time by using the Δ and ∇ keys.

2.4.5.4. Pump Configuration

Option: Pump 1, Pump 2 (if installed)



2.4.5.5. Dosing Type

Option: Linear Dosing or Dynamic Dosing

Dosing Type


Select the dosing type.

Linear Dosing
Dynamic Dosing

Select
Escape

Linear Dosing

Linear dosing dispenses a pre-defined volume of titrant with every addition.



Linear dosing is recommended for titrations with a slower reaction rate, difficult nonaqueous titrations, and specific applications.

Note: For steep and normal titration curves, smaller volume increments are recommended, to obtain many points around the equivalence point. For flat titration curves, larger volume increments are recommended for equivalence point detection.

To set the dosing volume, select Linear Dosing and enter the optimum dose.

Dosing volume ranges are:

5 mL burette	0.001 to 4.750 mL
10 mL burette	0.001 to 9.500 mL
25 mL burette	0.005 to 23.750 mL
50 mL burette	0.005 to 47.500 mL

Dynamic Dosing


The titrator determines the titrant dose by trying to maintain a certain potential change (*delta E*) with each addition.

After a titrant dose, if the potential change is lower than the set *delta E*, the next dose will be progressively increased until *max Vol* is attained. If the potential change is still lower than the set value, the titration will continue with *max Vol* doses.

After a titrant dose, if the potential change is higher than the set *delta E*, the next dose will be progressively decreased until *min Vol* is attained. If the potential change is still higher than the set value, the titration will continue with *min Vol* doses.

The titrant is added in volumes that depend on the proximity of the endpoint as shown in the graph below.

Dynamic dosing allows for larger doses far from the endpoint, reducing the total titration time. Closer to the endpoint, smaller doses are made, providing more data and improved accuracy.



Dynamic Dosing

Enter min Vol, max Vol and delta E.

0.030 mL - min Vol

0.500 mL - max Vol

4.500 mV - delta E

Press Next to move to the next entry.

Accept	Escape	Delete Digit	Next
--------	--------	--------------	------

The following parameters must be set:

min Vol The smallest dose to be dispensed during a titration.

The *min Vol* must be greater than or equal to:

5 mL & 10 mL burette 0.001 mL

25 mL & 50 mL burette 0.005 mL

max Vol The largest dose to be dispensed during a titration.

The *max Vol* must be less than or equal to 4.000 mL.

delta E Sets the fixed potential jump that has to be achieved after each titrant dose.

The allowed range is between 0.1 and 99.999 mV.

Recommendations for dosing parameters:

For steep and normal titration curves the recommended settings are:

delta E 3.5 to 9 mV

min Vol 0.010 to 0.025 mL (25 mL burette)

max Vol 0.075 to 0.250 mL (25 mL burette)

For flat titration curves the recommended settings are:

delta E 10 to 15 mV

min Vol 0.050 to 0.150 mL (25 mL burette)

max Vol 0.400 to 0.600 mL (25 mL burette)

To achieve the highest levels of accuracy and reproducibility, it is recommended that 20 to 80% of the nominal burette volume used for each titration is consumed. If lower volumes of titrant are required, a smaller burette can be used.

2.4.5.6. Endpoint Mode

Option: Equivalence Endpoint (pH or mV) or Fixed Endpoint (pH or mV)

Titration End Point Mode				
Select the end point detection.				
Equivalence End Point (pH)				
Equivalence End Point (mV)				
Fixed End Point (pH)				
Fixed End Point (mV)				
Select	Escape			

Fixed Endpoint (pH or mV)

Fixed Endpoint (pH)

Option: -2.000 to 20.000 pH

The titration is terminated when the preset pH value has been exceeded. The endpoint volume is a calculated value based on the dispensed volume when pH is under the preset value and the dispensed volume when pH exceeds the preset value.

Preset pH End Point				
Enter the end point pH value.				
8.600 pH				
The range is from -2.000 to 20.000 pH.				
Accept	Escape	Delete Digit		

Fixed Endpoint (mV)

Option: -2000.0 to 2000.0 mV

The endpoint detection algorithm is the same as for pH, but the threshold value is expressed in mV.


Preset mV End Point

Enter the end point mV value.

0.0 mV

The range is from -2000.0 to 2000.0 mV.

Accept	Escape	Delete Digit	
--------	--------	-----------------	--



Equivalence Endpoint (pH or mV)

The titration is terminated when the equivalence point is detected (the point where the added quantity of titrant equals the quantity of analyte present in the sample).

Titration End Point Mode

Select the end point detection.

Equivalence End Point (pH)

Equivalence End Point (mV)

Fixed End Point (pH)

Fixed End Point (mV)

Select	Escape		
--------	--------	--	--

Endpoint Determination

Option: 1st derivative or 2nd derivative

End Point Determination

Select the end point determination.

1st derivative

2nd derivative

Select	Escape		
--------	--------	--	--


The equivalence point detection algorithm requires three additional titrant doses to be dispensed after the equivalence point is reached.

The reported endpoint volume is a calculated value based on a number of points around the equivalence point.

The potentiometric titration curve is the response in mV potential, or pH, between the indication of the electrode versus the volume of titrant added.


The inflection point of the titration curve is assumed to be the equivalence point of the chemical reaction.

For non-symmetric titration curves, the theoretical error can be reduced by using the dynamic dosing.



1st Derivative


When 1st derivative is used to recognize the equivalence point, the titration curve inflection point (EQP) is the point where the 1st derivative reaches its maximum value.



The detection algorithm looks for the maximum value of the 1st derivative. The 1st derivative must be greater than the threshold value at the maximum point. See [2.4.5.7. Recognition Options \(Equivalence Endpoint Only\)](#) section for more information.

2nd Derivative

When 2nd derivative is used to recognize the equivalence point, the titration curve inflection point (EQP) is the point where the second derivative crosses zero.



The detection algorithm looks for the point where the second derivative changes sign.

The checked point, or 1st derivative, must be greater than the threshold value. See [2.4.5.7. Recognition Options \(Equivalence Endpoint Only\)](#) section for more information.

2.4.5.7. Recognition Options (Equivalence Endpoint Only)

The **Recognition Options** screen is a set of parameters used to avoid false detection of the equivalence point due to the chemical system (titrant / sample species and concentrations) and / or electrode response.

Recognition Options	
Select the options for equivalence point recognition.	
Threshold	500 mV/mL
Range	NO
Filtered Derivatives	NO
Select	Escape

Threshold


Option: 1 to 9999 mV / mL

This parameter must be set by the user according to the analysis.

The threshold represents the absolute value of the 1st derivative, expressed in mV / mL, which the detection algorithm does not search for the equivalence point.

Threshold				
Enter the threshold for equivalence point detection.				
EQ 1 Threshold: 500 mV/mL				
Recommended value is between: 1 and 450 mV/mL for FLAT Curve, 450 and 1800 mV/mL for NORMAL Curve, 1800 and 9999 mV/mL for STEEP Curve.				
Accept	Escape	Delete Digit		Next Threshold

The recommended value is 40% of the absolute value of the 1st derivative.



Depending on the titration curve profile, the following guide can be used:

- Flat** 1 to 450
- Normal** 50 to 1800
- Steep** 1800 to 9999


Range

Option: -2.000 to 20.000 pH or -2000.0 to 2000.0 mV


Range is an optional feature for equivalence point recognition.

Select Yes in the Range Options screen to enable.

The titrator will only look for an equivalence point between the set values.



The Limit 2 value must not be equal to the Limit 1 value.



Filtered Derivatives

Option: Yes or No

This option adds a filtering procedure in the 1st and 2nd derivative computation algorithm that reduces the influence of pH or mV noise.

Select Yes in the Filtered Derivative Option to enable.

Filtered Derivatives Option

Select option for filtered derivatives.

NO

YES

"NO" - without filtered derivatives.
"YES" - with filtered derivatives.

Select	Escape		
--------	--------	--	--

Noise can be due to:

- Chemical system properties (sample, titrant, solvent), such as slow chemical reactions or unbuffered samples such as wastewater, tap water, wine
- Electrode response
- Incorrect method parameters settings such as *Signal Stability*, *Stirring Speed*, etc.
- Insufficient titrant additions

Note: A shift in the endpoint volume by 1 or 2 doses may be seen due to filtering.

2.4.5.8. Pre-Titration Volume

During a titration, the equivalence point is reached after many titrant doses. These doses take up extra time while having no relevance for equivalence point detection.

Pre-titration volume adds a large initial dose to jump directly to the proximity of the equivalence point.

This first dose occurs after the pre-titration stir time is completed.

The ranges for pre-titration volumes are shown below:

5 mL burette 0.001 to 4.750 mL

10 mL burette 0.001 to 9.500 mL

25 mL burette 0.005 to 23.750 mL

50 mL burette 0.005 to 47.500 mL

Pre-Titration Volume				
Enter the initial titrant volume to be dispensed.				
9.000 mL				
Press Help to view the valid ranges for the pre-titration volume.				
Accept	Escape	Delete Digit		

To disable a pre-titration volume, enter 0.000 mL.

Note: A pre-titration volume is highly recommended whenever possible. Fewer doses will considerably shorten the overall titration duration.

2.4.5.9. Pre-Titration Stir Time

Option: 0 to 180 seconds

When enabled, the sample is mixed for a set period of time before any titrant is added. This allows the sample to become homogeneous.


Pre-Titration Stir Time				
Enter the initial mixing time prior to the start of the titration.				
10 seconds				
The range is from 0 to 180 seconds.				
Accept	Escape	Delete Digit		

Pre-titration stir time is disabled if 0 seconds is entered.

2.4.5.10. Measurement Mode


Option: Signal Stability or Timed Increment

During titration, the acquisition of the potential (mV) value of the solution can be done by using either *Signal Stability* or *Timed Increment* option.



Signal Stability

When *Signal Stability* is selected, the titrator acquires the potential (mV) only when stable conditions are reached. The principles of signal stability are plotted below:



The signal stability window (condition) represents the time interval (Δt) during which the potential measured in solution (mV) is confined inside the potential interval (ΔE).

The new signal value is acquired if the stability condition is reached after the minimum (t_{\min}) wait time.
If the stability condition is not reached and the maximum (t_{\max}) wait time has elapsed, the potential is acquired.

Signal Stability				
Enter mV variation (ΔE) in the time interval (Δt) min and max wait time period to the next sample measurement.				
<div style="text-align: center;"> 0.3 mV - ΔE </div>				
<div style="text-align: center;"> 2 seconds - Δt </div>				
<div style="text-align: center;"> 3 seconds - t_{\min} wait </div>				
<div style="text-align: center;"> 30 seconds - t_{\max} wait </div>				
Accept	Escape	Delete Digit	Next	

ΔE Maximum change in potential during *delta t*

The range is from 0.1 to 99.9 mV.

Δt The time interval during which the potential is measured.

The range is from 1 to 10 seconds.

t_{\min} wait The minimum elapsed time before a stability check. This is also the minimum elapsed time between two doses.

The range is from 2 seconds to *t max wait* time.

t_{\max} wait The maximum elapsed time between two successive doses. If the *t max wait* has elapsed, a new dose is added even if the signal stability condition is not reached.


The range is from *t min wait* time to 180 seconds.

Timed Increment

Option: 2 to 180 seconds

When *Timed Increment* is selected, the titrator acquires the potential (mV) at a fixed time interval (no signal stability check).

The time period between two acquisitions must be set according to the reaction and the response time of the electrode.



Timed Increment				
Enter the period of time to wait until the next dose.				
<div style="text-align: center;"> 5 seconds </div>				
The range is from 2 to 180 seconds.				
Accept	Escape	Delete Digit		

2.4.5.11. Electrode Type

Option: Up to 20 characters

Electrode Type

Select the highlighted letter by using the arrow keys then press "Enter".
Select the empty field for a space.
Press Accept to save the electrode type.

A	B	C	D	E	F	G	H	I	J	K	L					
M	N	O	P	Q	R	S	T	U	V	W	X	Y				
Z	a	b	c	d	e	f	g	h	i	j	k	l				
m	n	o	p	q	r	s	t	u	v	w	x	y				
z	À	Á	Â	Ã	Ä	Å	Ç	È	É	Ê	Ë	Ì	Í	Î	Ï	Ñ
Ò	Ó	Ô	Õ	Ö	Ù	Ú	Û	Ü	Ý	à	á	â	ã	ä	å	
ä	ç	è	é	ê	ë	ì	í	î	ï	ñ	ò	ó	ô	õ	ö	
ù	ú	û	ü	ý	î	*	\	_	\$	'	^	#	:			
0	1	2	3	4	5	6	7	8	9	%	.	,				
?	!	()	[]	<	>	=	/	+	-					

■ PH

Accept	Escape	Delete Letter	Cursor Left	Cursor Right
--------	--------	---------------	-------------	--------------

2.4.5.12. Blank Option

Option: Disabled, V-Blank, Blank-V

This feature allows the user to select the procedure for the blank calculations (where V is the volume of titrant dispensed during the titration and blank is the volume of titrant consumed by the blank sample).

View/Modify Method

Id: USER0002 Modified: 14:53 Jul 12, 2018

Select the option to be modified.

Method Revision:	1.0
Stirrer Configuration:	
Titration pump:	Pump 2
Dosing Type:	Dynamic
End Point Mode:	Fixed
Pre-Titration Volume:	
Pre-Titration Stir Time:	
Measurement Mode:	Sign
Electrode Type:	
Blank Option:	No Blank
Calculations:	Sample Calc. by Volume
Dilution Option:	Disabled
Titration Name:	0.1N HCl
Titration Conc.:	0.1000 N (eq/L)

Select	Escape	Print Method	Page Up	Page Down
--------	--------	--------------	---------	-----------

If one of the options (V-Blank or Blank-V) is selected in the **View / Modify Method** screen, the *blank value* will be active on the **View / Modify Method** screen and the value of the blank can be set (in liters).

Blank Value

Enter the blank volume in liters.

0.00125 L

Accept	Escape	Delete Digit	Exponent
--------	--------	--------------	----------

2.4.5.13. Calculations

The final result is calculated using the endpoint volume (titrant volume at the equivalence point or at the fixed endpoint), and a formula selected by the user.

Calculations												
Select either the calculation to be performed or modify the variables.												
<table border="1"> <tbody> <tr><td>Edit Variable Values</td></tr> <tr><td>No Formula (mL only)</td></tr> <tr><td>No Formula (L only)</td></tr> <tr><td>Sample Calc. by Weight</td></tr> <tr><td>Sample Calc. by Volume</td></tr> <tr><td>Stdz. Titrant by Weight</td></tr> <tr><td>Stdz. Titrant by Volume</td></tr> <tr><td>Generic Formula</td></tr> </tbody> </table>					Edit Variable Values	No Formula (mL only)	No Formula (L only)	Sample Calc. by Weight	Sample Calc. by Volume	Stdz. Titrant by Weight	Stdz. Titrant by Volume	Generic Formula
Edit Variable Values												
No Formula (mL only)												
No Formula (L only)												
Sample Calc. by Weight												
Sample Calc. by Volume												
Stdz. Titrant by Weight												
Stdz. Titrant by Volume												
Generic Formula												
Select	Escape											

Standard Titration Calculations

Edit Variable Values

Edit the variables in a previously selected calculation.

For each formula, selected variables can be changed.

No Formula (mL only)

Only the volume of titrant (mL) required to reach the endpoint will be displayed.

No Formula (L only)

Only the volume of titrant (L) required to reach the endpoint is displayed.

Sample Calculations by Weight

Titrant units

Option: M (mol / L), N (eq / L), g / L, mg / L

Final result units

Option: ppt (g / kg), ppm (mg / kg), ppb (μg / kg), % (g / 100 g), mg / g, mg / kg, mol / kg, mmol / g, eq / kg, meq / kg

This calculation is used when the concentration of an analyte is determined by the weight of the sample. The results are based on the initial sample weight (in grams).

The titrator will calculate the results based on the titrant and sample units selected.

A formula example is shown below using M (mol / L) as the titrant unit and ppt (g / kg) as the final result unit:

Variables can be set according to the amount of sample and titrant used.

Calculating Sample Concentration				
M (mol/L) --> ppt (g/kg)				
The calculation is:				
$\frac{U \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}} \times \frac{\text{g}}{\text{mol}}}{\frac{\text{g}}{\text{kg}} \times 1000\text{g}}$				
Select the variables to change value. U = volume dispensed in liters.				
1.000 mol/L -> titrant conc. 1.000 mol/mol -> (sample/titrant) 1.000 g/mol -> mw of sample 1.000 g -> sample weight				
Select	Escape	Save / Exit		

Sample Calculations by Volume

Titrant Units

Option: M (mol / L), N (eq / L), g / L, mg / L

Final Result Units

Option: ppt (g / L), ppm (mg / L), ppb (μg / L), M (mol / L), N (eq / L), mg / L, μg / L, mmol / L, mg / mL, mg / 100 mL, g / 100 mL, eq / L, meq / L

This calculation is used when the concentration of an analyte is determined in terms of the volume of sample. The results are based on the initial sample volume (in milliliters).

The titrator will calculate the results based on the selected units.

Calculating Sample Concentration				
N (eq/L) --> ppt (g/L)				
The calculation is:				
$\frac{U \times \frac{\text{eq}}{\text{L}} \times \frac{\text{mol}}{\text{eq}} \times \frac{\text{g}}{\text{mol}}}{\frac{\text{mL}}{\text{L}} \times 1000\text{mL}}$				
Select the variables to change value.				
1.000 eq/L -> titrant conc. 1.000 mol/eq -> (sample/titrant) 1.000 g/mol -> mw of sample 1.000 mL -> sample volume				
Select	Escape	Save / Exit		

A formula example is shown below using N (eq / L) as the titrant units and g / L as the final result units. Variables can be set according to the amount of sample and titrant used.

Calculating Sample Concentration

M (mol/L) --> mol/L

The calculation is:

$$\frac{U \times \frac{\text{mol}}{\text{L}} \times \frac{\text{mol}}{\text{mol}}}{\text{mL} \times \frac{\text{L}}{1000\text{mL}}}$$

Select the variables to change value.
U = volume dispensed in liters.

1.000 mol/L -> titrant conc.
1.000 mol/mol -> (sample/titrant)
100.00 mL -> sample volume

Select	Escape	Save / Exit	
--------	--------	----------------	--

Standardize Titrant by Weight

Option: M (mol /L), N (eq / L), g / L, mg / L

This calculation is used when the concentration of the titrant is determined using a solid standard. Determination of the titrant concentration is based on the primary standard weight (in grams).

The calculation is based on the selected titrant unit. If the titrant unit is M (mol / L), the formula used to calculate the result is displayed below.

Calculating Titrant Concentration

The titrant concentration unit is M (mol/L).

The calculation is:

$$\frac{\text{g} \times \frac{\text{mol}}{\text{g}} \times \frac{\text{mol}}{\text{mol}}}{U}$$

Select the variables to change value.
U = volume dispensed in liters.

0.200 g -> standard weight
204.23 g/mol -> mw of standard
1.000 mol/mol -> (titrant/standard)

Select	Escape	Save / Exit	
--------	--------	----------------	--

Standardize Titrant by Volume

Option: M (mol / L), N (eq / L), g / L, mg / L

This calculation is used when the concentration of the titrant is determined using a primary standard solution. Determination of the titrant concentration is based on the primary standard volume (in milliliters).

The calculation is based on the selected titrant unit. If the titrant unit is N (eq / L), the formula used to calculate the result is displayed below.

Calculating Titrant Concentration

The titrant concentration unit is N (eq/L).

The calculation is:

$$\frac{\text{mL} \times \frac{\text{L}}{1000\text{mL}} \times \frac{\text{eq}}{\text{L}}}{\text{U}}$$

Select the variables to change value.
U = volume dispensed in liters.

1.684 mL -> standard volume

2.375 eq/L -> standard conc.

Select	Escape	Save / Exit	
--------	--------	----------------	--

Generic Formula

Final results units:

Option: ppt (g / kg), ppt (g / L), ppm (mg / kg), ppm (mg / L), ppb (μg / kg), ppb (μg / L), % (g / 100 g), M (mol / L), mg / g, N (eq / L), g / L, mg / kg, mg / L, mol / kg, μg / L, mol / L, mmol / g, eq / kg, mmol / L, meq / kg, mg / mL, mg / 100 mL, g / 100 mL, eq / L, meq / L, no unit

Users can define their own calculation formula based on the final result units in a solid or liquid sample.

The titrator will calculate the results based on the selected unit.

The formula can be either for titrant standardization or sample analysis.

Calculating Sample Concentration

Final unit is mg/L.

The calculation is:

$$\frac{C \times U \times F1 \times F2 \times F3}{S}$$

Select the variables to change value.
U = volume dispensed in liters.

1.000 C -> (titrant conc.)

1.000 F1 -> (general factor)

1.000 F2 -> (general factor)

1.000 F3 -> (general factor)

Select	Escape	Save / Exit	
--------	--------	----------------	--

C the concentration of the titrant

F1, F2, F3 general factor

S sample size, in grams or milliliters

V the volume delivered, in liters, to reach the endpoint

General factors

- Weight conversion** mol / L, eq / L, g / L, mg / L
- Reaction ratio** mol / mol, mol / eq, eq / mol
- Unit conversion** L to mL, g to mg
- Weight conversion** kg, g, mg, μ g, mole, mmole

2.4.5.14. Dilution Option

Option: Enabled or Disabled

When the initial sample is diluted, a titration is made with an aliquot of the diluted sample, dilution calculations can be used. The calculations are based on the original sample weight or volume in order to express the results for the initial sample.

Dilution Parameters

Select the option.

Final Dilution Volume:	100.000 mL
Aliquot Volume:	10.000 mL
Analyte size to be diluted:	1.0000 mL

Select	Escape		
--------	--------	--	--

- Final Dilution Volume** The volume of the sample after dilution
- Aliquot Volume** Volume of sample taken from the dilution for titration
- Analyte size to be diluted** The initial sample weight or volume

2.4.5.15. Titrant Name

Option: Up to 15 characters

Titrant Name

Select the highlighted letter by using the arrow keys then press "Enter".
 Select the empty field for a space.
 Press Accept to save the entered text.

■	A	B	C	D	E	F	G	H	I	J	K	L				
	M	N	O	P	Q	R	S	T	U	V	W	X	Y			
	Z	a	b	c	d	e	f	g	h	i	j	k	l			
	m	n	o	p	q	r	s	t	u	v	w	x	y			
	z	À	Á	Â	Ã	Ä	Å	Ç	È	É	Ê	Ë	Ì	Í	Î	Ï
	Ò	Ó	Ô	Õ	Ö	Ù	Ú	Û	Ü	ß	à	á	â	ã		
	ä	å	ç	è	é	ê	ë	ì	í	î	ï	ò	ó	ô	õ	ö
	ù	ú	û	ü	ÿ	ı	*	\	_	\$	'	^	#	:		
	0	1	2	3	4	5	6	7	8	9	%	.	,			
	?	!	()	[]	<	>	=	/	+	-				

■ 0.1N NaOH

Accept	Escape	Delete Letter	Cursor Left	Cursor Right
--------	--------	---------------	-------------	--------------

2.4.5.16. Titrant Concentration

Enter the concentration of the titrant to be used. When determining the titrant concentration, only the concentration unit is displayed.

Titrant Concentration				
Enter the titrant concentration.				
0.10123 M (mol/L)				
Accept	Escape	Delete Digit	Exponent	

2.4.5.17. Analyte Size

Option: 0.001 to 250.0

Enter the size of the sample (for sample-concentration determination) or standard (for titrant-concentration determination).

Sample Volume				
Enter the initial sample volume in milliliters.				
1.0000 mL				
This volume will be used when fixed sample size is selected.				
Accept	Escape	Delete Digit	Exponent	

2.4.5.18. Analyte Entry

Option: Fixed or Manual

Analyte Entry						
Select the entry mode of analyte.						
<table border="1" style="width: 100%;"> <tr> <td style="text-align: center;">Fixed Weight or Volume</td> </tr> <tr> <td style="text-align: center;">Manual Weight or Volume</td> </tr> </table>					Fixed Weight or Volume	Manual Weight or Volume
Fixed Weight or Volume						
Manual Weight or Volume						
Verify the correct formula is being used, i.e. weight or volume analyte type.						
Select	Escape					

Fixed Weight or Volume
Manual Weight or Volume

For each titration will use a set weight or volume in the calculations.

For each titration the exact weight or volume can be entered at the beginning of each titration.

2.4.5.19. Maximum Titrant Volume

Option: 0.100 to 100.000 mL

The maximum titrant volume used in the titration must be set according to the analysis.

If the titration endpoint (fixed or equivalence point) is not reached, the titration will be terminated after the maximum titrant volume has been dispensed. The error message "Limits Exceeded" will appear on the display.

Maximum Titrant Volume				
Enter the maximum titrant volume to be dispensed.				
15.000 mL				
Recommend the total volume of the burette.				
Accept	Escape	Delete Digit		

2.4.5.20. Potential Range

Option: -2000.0 to 2000.0 mV

The input potential range can be set by the user. The titration will be terminated and an error message will appear if the potential is outside these limits.

These limits provide protection against a titration that does not generate an endpoint due to potential over-range.

Potential Range				
Enter the upper and lower potential.				
2000.0 mV - Upper Limit				
-2000.0 mV - Lower Limit				
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	

2.4.5.21. Volume / Flow Rate

The flow rate for the dosing system can be set by the user in an interval of 0.3 to 2 times the burette volume:

5 mL burette 0.3 to 10 mL/min

10 mL burette 0.3 to 20 mL/min

25 mL burette 0.3 to 50 mL/min

50 mL burette 0.3 to 100 mL/min

The flow rate is set for all burette operations.

Flow Rate				
Enter the titrant/reagent flow rate.				
50.0 mL/min				
The range is from 0.3 to twice the total volume of the burette.				
ACCEPT	Escape	Delete Digit		

Note: The titrator will automatically detect the burette size and display the correct high limit volume.

2.4.5.22. Signal Averaging

Option: 1, 2, 3, 4 readings

This option enables filtering on the mV / pH reading.


If 1 Reading is selected, the filtering is disabled. The titrator will take the last reading and place it into a "moving window" along with the last 2, 3 or 4 readings (depending on the selected option). The average of those readings is displayed and used for calculations. Averaging more readings is helpful when a noisy signal is received from the electrode.

View/Modify Method																																
Id: USER0001 Modified: 14:39 Jun 28, 2018																																
Select the option to be modified.																																
<table border="1"> <tbody> <tr> <td>Measurement Mode:</td> <td>Signal Stability</td> </tr> <tr> <td>Electrode Type:</td> <td>pH</td> </tr> <tr> <td>Blank Option:</td> <td>No Blank</td> </tr> <tr> <td>Calculations:</td> <td>Sample Calc. by Volume</td> </tr> <tr> <td>Dilution Option:</td> <td>Disabled</td> </tr> <tr> <td>Titrant Name:</td> <td>0.1N NaOH</td> </tr> <tr> <td>Titrant Conc.:</td> <td>1.0000000</td> </tr> <tr> <td>Analyte Size:</td> <td>1 Reading</td> </tr> <tr> <td>Analyte Entry:</td> <td>2 Readings</td> </tr> <tr> <td>Maximum Titrant Volume:</td> <td>3 Readings</td> </tr> <tr> <td>Potential Range:</td> <td>4 Readings</td> </tr> <tr> <td>Volume/Flow Rate:</td> <td>25 mL</td> </tr> <tr> <td>Signal Averaging:</td> <td>1 Reading</td> </tr> <tr> <td>Significant Figures:</td> <td>XXXXXX</td> </tr> </tbody> </table>					Measurement Mode:	Signal Stability	Electrode Type:	pH	Blank Option:	No Blank	Calculations:	Sample Calc. by Volume	Dilution Option:	Disabled	Titrant Name:	0.1N NaOH	Titrant Conc.:	1.0000000	Analyte Size:	1 Reading	Analyte Entry:	2 Readings	Maximum Titrant Volume:	3 Readings	Potential Range:	4 Readings	Volume/Flow Rate:	25 mL	Signal Averaging:	1 Reading	Significant Figures:	XXXXXX
Measurement Mode:	Signal Stability																															
Electrode Type:	pH																															
Blank Option:	No Blank																															
Calculations:	Sample Calc. by Volume																															
Dilution Option:	Disabled																															
Titrant Name:	0.1N NaOH																															
Titrant Conc.:	1.0000000																															
Analyte Size:	1 Reading																															
Analyte Entry:	2 Readings																															
Maximum Titrant Volume:	3 Readings																															
Potential Range:	4 Readings																															
Volume/Flow Rate:	25 mL																															
Signal Averaging:	1 Reading																															
Significant Figures:	XXXXXX																															
Select	Escape	Print Method	Page Up	Page Down																												

2.4.5.23. Significant Figures

Option: Two (XX), Three (XXX), Four (XXXX) or Five (XXXXX)

This option allows you to set the format for displaying the final titration result.



2.4.6. PRINTING

To print method parameters, press Method Options from the main screen.

Press Print Method and wait a few seconds until the printer completes the job.

If no printer is connected to the dedicated socket, or if the printer is offline, an error message will appear on the display. See [2.10.2.1. Connecting to a Printer](#) section for details on connecting a printer to the titrator.

2.5. TITRATION MODE

2.5.1. RUNNING A TITRATION

Before beginning a titration, make sure that the following conditions are met:

- At least one pump is properly installed.
- A burette is inserted in the pump and filled with titrant.
- The aspiration tube is inserted in the titrant bottle and primed. The dispensing tube is over the titration beaker.
- The standard or sample has been carefully weighed / measured into the beaker.
- The electrode(s) and the temperature probe are submersed in the beaker.
- The desired method is selected and the parameters are set to the optimal values.

2.5.1.1. Starting a Titration

To start a new analysis, press start stop from the main screen.

When an analysis begins:


- The stirrer will turn on, if enabled. See [2.4.5.3. Stirrer Configuration](#) section for more information.
- The pre-titration volume will be dispensed, if enabled. See [2.4.5.8. Pre-Titration Volume](#) section for more information.
- After the pre-titration volume is added the pre-titration stir time starts, if enabled. See [2.4.5.9. Pre-Titration Stir Time](#) section for more information.
- The titrator will start the analysis and continue to deliver titrant until the endpoint is detected or the titration is terminated.

2.5.1.2. Suspending a Titration

While a titration or analysis is in progress, you can temporarily stop it by pressing Suspend. This will stop the dosing pump if it is running.

To continue the titration or analysis press Resume.

2.5.1.3. Viewing the Titration Curve

During a titration, the potentiometric curve and the derivative curve (equivalence point only) can be displayed on the **Titration Graph** screen by pressing .

The potentiometric curve and the derivative curve are scaled to fit simultaneously inside the display.

When a titration endpoint is successfully detected, the volume is displayed on the graph and marked with an "x".

The contents of the graph as related to an endpoint type are as follows:

Equivalence endpoint (pH) The pH readings and the selected derivative vs. volume of titrant are displayed (see Figure 1).

Equivalence endpoint (mV) The mV readings and the selected derivative vs. volume of titrant are displayed (see Figure 2).

Fixed endpoint (pH) The pH readings vs. volume of titrant are displayed (see Figure 3).

Fixed endpoint (mV) The mV readings vs. volume of titrant are displayed (see Figure 4).




Figure 1. Equivalence endpoint (pH)




Figure 2. Equivalence endpoint (mV)




Figure 3. Fixed endpoint (pH)






Figure 4. Fixed endpoint (mV)

 Changes the y-axis from the pH (mV) reading to the derivate value (equivalence point titrations only).

 Saves the graph as a bitmap (available when titration is complete).

2.5.2. STOPPING A TITRATION

The titration or analysis is terminated when one of the following conditions is met:

Titration completed

This is the only mode with valid final result values. The endpoint or stable reading was successfully detected, the final results will be displayed.

Manually terminated

The current titration or analysis was terminated by the user before the endpoint was detected.

Limits exceeded

The maximum titrant volume was delivered without reaching the endpoint. An error message is displayed on the screen.

Critical error


A critical error occurred and the titration was stopped. These errors are typically related to the dosing system. An error message is displayed on the screen.

Potential out of range

The measured values from the electrode are outside the potential range. An error message is displayed on the screen.

2.6. pH MODE


By pressing Mode from the main screen, the titrator can be switched to **Titration, pH, mV** or **ISE** modes.



When one of these keys is pressed, the titrator will enter the selected mode:

- Titration Switches to **Titration** mode.
- pH Switches to **pH** mode.
- mV Switches to **mV** mode.
- ISE Switches to **ISE** mode.

2.6.1. DISPLAY



- 1. Time and Date
- 2. pH value
- 3. Status Bar
- 4. Temperature compensation status
- 5. Temperature reading
- 6. Units
- 7. Stirrer information
- 8. Logging info
- 9. mV reading
- 10. Virtual option keys

If one of these keys is pressed, the associated function is immediately performed. Some of the keys are active only in specific screens:

General Options Gives access to options that are not directly related to the measurement process. See [2.3. General Options](#) section for more information.

Save Reading Stores the current pH reading. See [2.6.4. Logging](#) section for more information.

Or

Start Log Starts the interval log. See [2.6.4. Logging](#) section for more information.

pH Calibr. Enters the pH calibration screen. See [2.6.3. pH Calibration](#) section for more information.

pH Setup Enters the pH setup screen, parameters are associated with pH measurements and calibration. See [2.6.2. pH Setup](#) section for more information.

Mode Allows the user to switch between the available measurement modes: **Titrator**, **pH**, **mV** or **ISE** mode.

2.6.2. pH SETUP

To access pH Setup, press **pH Setup** option key while in pH mode.

```

pH Setup

Select a menu option.

Buffer Entry Type:           Manual
First Cal Point:             Point
Edit Custom Buffers
Edit Buffer Group
Calibration Reminder:       Periodic
Set Reminder Period:        10d:02h:30m
Clear Calibration
pH GLP Data
Logging Interval:           0h:00m:02s
Stability Criteria:         Medium
pH Resolution:              X.XXX
Stirrer Configuration:      Stirrer 1
Stirring Speed:             1200 RPM

Select  Escape
  
```

Use **▲** and **▼** keys to highlight the desired option.

Press **Select** or **enter** to access the selected option.

2.6.2.1. Buffer Entry Type

Option: **Automatic**, **Semiautomatic**, **Manual**

```

pH Setup

Select a menu option.

Buffer Entry Type:           Manual
First Cal Point:             Point
Edit Custom Buffers
Edit Buffer Group
Calibration Reminder:       Periodic
Set Reminder Period:        10d:02h:30m
Clear Calibration
pH GLP Data
Logging Interval:           Disabled
Stability Criteria:         Medium
pH Resolution:              X.XXX
Stirrer Configuration:      Disabled

Automatic
Semiautomatic
Manual

Select  Escape
  
```

Automatic The instrument automatically selects the pH calibration point as the closest buffer from the predefined buffer group. See [2.6.2.4. Edit Buffer Group](#) section for more information.

Semiautomatic The instrument automatically selects the closest buffer from the available buffers (standard and custom buffers).

Manual The calibration buffer must be manually selected during calibration from the available buffer list (standard and custom buffers).

2.6.2.2. First-Calibration Point

Option: Point or Offset

pH Setup

Select a menu option.

Buffer Entry Type:	Manual
First Cal Point:	Point
Edit Custom Buffers	Point Offset
Edit Buffer Group	
Calibration Reminder:	
Set Reminder Period:	
Clear Calibration	
pH GLP Data	
Logging Interval:	Disabled
Stability Criteria:	Medium
pH Resolution:	X.XXX
Stirrer Configuration:	Disabled

Select	Escape		
--------	--------	--	--

Point The slope values adjacent to the calibration points will be reevaluated (normal calibration).

Offset The existing slope values will not be changed.

2.6.2.3. Edit Custom Buffers

If you wish to use buffers other than the standard ones, use the Edit Custom Buffers option to set the desired pH value. Up to five pH custom buffers can be set.

Note: Custom buffers are not temperature compensated, enter the value of the buffer at the calibration temperature.





Edit Custom Buffers


Press <Edit> to edit selected buffer.
Press <Remove Buffer> to delete the custom buffer.



Cust	Cust	Cust	Cust	Cust
6.870	9.230	12.750	----	----

Use arrows keys to select the buffer.

Remove Buffer	Escape	Edit	◀	▶
---------------	--------	------	---	---

1. Use the  and  keys to select the desired buffer.
2. Press  to delete the selected buffer.
3. Press  to edit the selected buffer.




4. Use the numeric keypad to enter the pH buffer value.
5. Press  to save the value.
6. Press  to return to pH Setup menu.

2.6.2.4. Edit Buffer Group

Option: Up to five buffers

Select up to five buffers from the available buffers (Hanna[®] or custom) to be used for automatic buffer recognition. Within the buffer group, pH values must be at least 1.5 pH apart.

If the buffer group already contains five pH buffers, at least one pH buffer has to be removed in order to add another buffer.



Use the arrow keys to select the pH buffer to be included / removed in / from the buffer group.

 or  Adds or removes the selected pH buffer to / from buffer group.

 Returns to pH Setup menu.

2.6.2.5. Calibration Reminder

Option: Daily, Periodic, Disabled

Calibration Reminder							
Select a menu option.							
<table border="1"> <tr> <td>Daily</td> </tr> <tr> <td>Periodic</td> </tr> <tr> <td>Disabled</td> </tr> </table>					Daily	Periodic	Disabled
Daily							
Periodic							
Disabled							
Select	Escape						

Daily The calibration reminder will appear daily at a specified time.

Periodic The calibration reminder will appear after the set time since the last calibration has elapsed.

Disabled The calibration reminder will not appear.

2.6.2.6. Set Reminder Period

Option: Disabled to 31 days, 23 hours and 59 minutes

If *Daily* or *Periodic* option was selected for the calibration reminder, the reminder period must also be set.

For a daily reminder period, the time of day can be set.

For a periodic reminder period, the number of days, hours and minutes can be set.

Periodic Calibration Reminder				
Enter the time period that must be passed since the last calibration, whereafter the calibration reminder appears.				
10	2	30		
days	hours	minutes		
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

Next Moves the cursor to the next field.

Accept Saves the changes or Escape to return to the previous screen.

Off Disables the calibration reminder and return to pH setup.

2.6.2.7. Clear Calibration

This option clears the existing pH calibration for the selected channel. If the calibration is cleared, the factory calibration will be used.

Clear Clears the previous calibration or **Escape** to return to the previous screen without clearing the calibration.

Clear Calibration				
Press <Clear> to clear all calibration points.				
Press <Escape> to return without clearing the calibration points.				
Clear	Escape			

2.6.2.8. pH GLP Data

Display the pH calibration data.

pH GLP Data				
Analog 1				
Last Calibration: 10:13 May 24, 2018				
Offset: -0.1 mV Average Slope: 100.7%				
1.679pH (Hanna) 316.2mV 26.3°C A				
10:10:30 May 24, 2018				
4.010pH (Hanna) 177.5mV 26.3°C A				
10:09:11 May 24, 2018				
7.010pH (Hanna) -0.6mV 26.3°C A				
10:08:40 May 24, 2018				
10.010pH (Hanna) -179.1mV 26.3°C A				
10:09:43 May 24, 2018				
12.450pH (Hanna) -325.6mV 26.3°C A				
10:13:15 May 24, 2018				
	Escape			

2.6.2.9. Logging Interval

Option: 2 seconds to 8h 59 min 59 sec


Set the logging interval to be used- for automatic logging.

Select Off to enable manual logging.

Logging Interval				
Enter the data logging interval.				
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <div style="border: 1px solid black; padding: 2px;">0</div> hours </div> <div style="text-align: center;"> <div style="border: 1px solid black; padding: 2px;">0</div> minutes </div> <div style="text-align: center;"> <div style="border: 1px solid black; padding: 2px;">2</div> seconds </div> </div>				
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

2.6.2.10. Signal Stability Criteria

Option: Fast, Medium, Accurate




Fast Quicker results, less accuracy

Medium Medium speed results, medium accuracy

Accurate Slower results, high accuracy


2.6.2.11. pH Resolution

Option: One (X.X), Two (X.XX), Three (X.XXX) decimal places




2.6.2.12. Stirrer Configuration

Option: Disabled or Stirrer 1



2.6.2.13. Stirring Speed

Option: 200 to 2500 RPM




2.6.3. pH CALIBRATION

Calibrate the instrument often, especially if high accuracy is required.

The instrument should be recalibrated:

- Whenever the pH electrode is replaced.
- At least once a week.
- After testing aggressive chemicals.
- When “No pH Calibration” or “pH Calibration Expired” message appears on the display.



Preparation

Pour small quantities of the buffer solutions into clean beakers. If possible, use plastic beakers to minimize any EMC interferences. For accurate calibration and to minimize cross-contamination, use two beakers for each buffer solution: one for rinsing the electrode and one for calibration. If you are measuring in the acidic range, use pH 7.01 or 6.86 as the first buffer and pH 4.01 / 3.00 or 1.68 as the second buffer. If you are measuring in the alkaline range, use pH 7.01 or 6.86 as the first buffer and pH 10.01 / 9.18 or 12.45 as the second buffer.

For extended range measurements (acidic and alkaline), perform a five-point calibration by selecting five buffers across the entire pH range.

Calibration procedure

During calibration, the user has a choice of 8 standard buffers: (pH 1.68, 3.00, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45) and up to 5 custom buffers.

For accurate measurements it is recommended to perform a five-point calibration. However, at least a two-point calibration is suggested. For pH titrations, the selected buffers should bracket your endpoint (e.g. if your endpoint value is at 8.5, use 7.01 or 6.86 and 9.18 or 10.01 for calibration).

To begin calibration:


1. Press pH
Calibr.. If the instrument was calibrated before, previous calibration can be cleared by pressing Clear
Cal..

Note: It is very important to clear calibration history when a new electrode is used.

2. Immerse the pH electrode and the temperature probe approximately 4 cm (1.5") into a buffer solution and stir gently.
3. If necessary, select the pH calibration buffer value with Next
Buffer or Previous
Buffer.
4. Once the reading has stabilized, press Accept to update the calibration. The calibration buffer will be added to the Calibrated Buffers section.
5. Rinse the pH electrode and the temperature probe, then immerse them into the next buffer solution and follow the above procedure or press Escape to exit the calibration.

Notes:

- The new calibration points will replace old ones if the difference between them is ± 0.2 pH.
- Buffers used in previous calibrations will not have a solid background.
- If calibrating with a standard buffer in MTC mode, the pH value and temperature can be modified by pressing Edit. The values can be adjusted using the numeric keys. Press Accept to save the new values.



- In ATC mode, the pH value for custom buffers can be modified by pressing Edit.
- If the Automatic Buffer entry type was selected for the calibration procedure, the titrator will automatically select the buffer closest to the measured pH value from the buffer group.
- If the Semiautomatic Buffer entry type was selected, use the Previous
Buffer or Next
Buffer to select the buffer. Only buffers in the buffer group will be displayed.

Calibration messages:

Wrong Buffer. Please check the buffer.

The message is displayed when the difference between the pH reading and the value of the selected calibration buffer is significant. Check if you have selected the appropriate calibration buffer.

Wrong buffer temperature.

The message is displayed if the buffer temperature is out of the defined temperature range.

Clean the electrode or check the buffer. Press Accept to update calibration.

This message alerts the user that some dirt or deposits could be on the electrode, or the buffer is contaminated.

Slope too low. Please check the buffer.

This message appears if the current slope is under 80% of the default slope. Recalibrate the instrument using fresh buffers.

Slope too high. Press to clear the old calibration.

This message appears as a result of an erroneous slope condition.


2.6.4. LOGGING

Data logging is available in pH mode. It can be logging on demand (Manual Logging) or automatically (Interval Logging) at predefined time intervals.

The the logging report can be customized. See [2.9.3.5. Setting Up pH / mV / ISE Report](#) for more information.

2.6.4.1. Interval Logging


The logging interval is set in the pH Setup screen.


Press  to start the log.

The logging interval and name of logging file will be displayed on the measure screen.

To stop the automatic logging, press Stop.


2.6.4.2. Manual Logging

To manually log pH readings, press  from the pH measurement screen.

A new record will be added to the report every time  is pressed.

2.7. mV MODE

By pressing  from the main screen, the titrator can be switched to **Titrator**, **pH**, **mV** or **ISE** modes.



When one of these keys is pressed, the titrator will enter the selected mode:

 Switches to **Titrator** mode.


 Switches to **pH** mode.

 Switches to **mV** mode.

 Switches to **ISE** mode.

2.7.1. DISPLAY

The mV screen is shown below:



If one of these keys is pressed, the associated function is immediately performed. Some of the keys are active only in specific screens:

General Options Gives access to options that are not directly related to the measurement process. See [2.3. General Options](#) section for more information.

Save Reading Stores the current pH reading. See [2.7.4. Logging](#) section for more information.

Or


Start Log Starts the interval log. See [2.7.4. Logging](#) section for more information.

mV Calibr. Enters the pH calibration screen. See [2.7.3. Relative mV Calibration](#) section for more information.

pH Setup Enters the pH setup screen, parameters are associated with pH measurements and calibration. See [2.7.2. mV Setup](#) section for more information.

Mode Allows the user to switch between measurement modes: **Titration**, **pH**, **mV** or **ISE** mode.

2.7.2. mV SETUP



2.7.2.1. Clear Relative mV Offset

Clear Clears the relative mV offset or **Escape** to return to the previous screen.

<p>Clear Relative mV Offset</p> <p>Press Clear to clear the relative mV offset.</p> <p>Press Escape to return without clearing the relative mV offset.</p>				
Clear	Escape			

2.7.2.2. Logging Interval

Option: 2 seconds to 8h 59min 59sec

Press Off to enable manual logging.

<p>Logging Interval</p> <p>Enter the data logging interval.</p> <p style="text-align: center;"> 0 0 2 hours minutes seconds </p> <p>Press Next to move to the next entry.</p>				
Accept	Escape	Delete Digit	Next	Off

2.7.2.3. Stability Criteria

Option: Fast, Medium, Accurate

<p>mV Setup</p> <p>Select a menu option.</p> <table border="1" style="width: 100%;"> <tr> <td style="width: 80%;">Clear Relative mV Offset</td> <td style="width: 20%;"></td> </tr> <tr> <td>Logging Interval:</td> <td style="text-align: right;">0h:00m:02s</td> </tr> <tr> <td>Stability Criteria:</td> <td style="text-align: right;">Fast</td> </tr> <tr> <td>Stirrer Configuration:</td> <td style="text-align: right;">1</td> </tr> <tr> <td>Stirring Speed:</td> <td style="text-align: right;">1</td> </tr> </table> <div style="border: 1px solid black; padding: 5px; margin-top: 5px; width: fit-content;"> <p style="text-align: center;"> Fast Medium Accurate </p> </div>					Clear Relative mV Offset		Logging Interval:	0h:00m:02s	Stability Criteria:	Fast	Stirrer Configuration:	1	Stirring Speed:	1
Clear Relative mV Offset														
Logging Interval:	0h:00m:02s													
Stability Criteria:	Fast													
Stirrer Configuration:	1													
Stirring Speed:	1													
Select	Escape													

Fast Quicker results, less accuracy

Medium Medium speed results, medium accuracy

Accurate Slower results, high accuracy

2.7.2.4. Stirrer Configuration

Option: Stirrer 1 or Disabled

mV Setup												
Select a menu option.												
<table border="1"> <tr> <td>Clear Relative mV Offset</td> <td>Disabled</td> </tr> <tr> <td>Logging Interval:</td> <td>Fast</td> </tr> <tr> <td>Stability Criteria:</td> <td>Stirrer 1</td> </tr> <tr> <td>Stirring Speed:</td> <td>Disabled</td> </tr> </table>					Clear Relative mV Offset	Disabled	Logging Interval:	Fast	Stability Criteria:	Stirrer 1	Stirring Speed:	Disabled
Clear Relative mV Offset	Disabled											
Logging Interval:	Fast											
Stability Criteria:	Stirrer 1											
Stirring Speed:	Disabled											
<table border="1"> <tr> <td>Disabled</td> <td>Stirrer 1</td> </tr> </table>					Disabled	Stirrer 1						
Disabled	Stirrer 1											
Select	Escape											

2.7.2.5. Stirring Speed

Option: 200 to 2500 RPM

Stirring Speed						
Enter the speed of the stirrer within below range.						
<table border="1"> <tr> <td>1100</td> <td>RPM</td> </tr> </table>					1100	RPM
1100	RPM					
The range is from 200 to 2500 RPM.						
Accept	Escape	Delete Digit				

2.7.3. RELATIVE mV CALIBRATION

Relative mV																
Analog 1																
Set the value for the relative mV offset.																
<table border="1"> <tr> <td>Absolute mV:</td> <td>2.7</td> <td>mV</td> </tr> <tr> <td></td> <td></td> <td>Stirrer 1</td> </tr> <tr> <td></td> <td></td> <td>1100/1100 RPM</td> </tr> <tr> <td>Relative mV:</td> <td>2.7</td> <td>mV</td> </tr> </table>					Absolute mV:	2.7	mV			Stirrer 1			1100/1100 RPM	Relative mV:	2.7	mV
Absolute mV:	2.7	mV														
		Stirrer 1														
		1100/1100 RPM														
Relative mV:	2.7	mV														
Low limit: -1997.3 mV																
High limit: 2002.7 mV																
Accept	Escape	Delete Digit														

Accept

 Accepts the value.

Escape

 Cancels this operation and return to the previous screen.

Delete Digit

 Deletes the last digit.


2.7.4. LOGGING

Data logging is available in mV mode. It can be logging on demand (Manual Logging) or automatically (Interval Logging) at predefined time intervals.

The the logging report can be customized. See [2.9.3.5. Setting Up pH / mV / ISE Report](#) section for more information.

2.7.4.1. Interval Logging


The logging interval is set in the mV Setup screen.


Press  to start the log.

The logging interval and name of logging file will be displayed on the measure screen.


To stop the automatic logging, press Stop.


2.7.4.2. Manual Logging

To manually log mV readings, press  from the mV measurement screen.

A new record will be added to the report every time  is pressed.

2.8. ISE MODE

By pressing  from the main screen, the Titrator can be switched to **Titration**, **pH**, **mV** or **ISE** modes.



When one of these keys is pressed, the titrator will enter the selected mode:

 Switches to **Titration** mode.


 Switches to **pH** mode.

 Switches to **mV** mode.

 Switches to **ISE** mode.

2.8.1. DISPLAY

The ISE screen is shown below.




If one of these keys is pressed, the associated function is immediately performed. Some of the keys are active only in specific screens:

- General Options Gives access to options that are not directly related to the measurement process. See [2.3. General Options](#) section for more information.
- Save Reading Stores the current concentration reading. See [2.8.4. Logging](#) section for more information.
- OR
- Start Log Starts the interval log. See [2.8.4. Logging](#) section for more information.
- ISE Calibr. Enters the ISE calibration screen. See [2.8.3. ISE Calibration](#) section for more information.
- ISE Setup Enters the ISE setup screen. Parameters are associated with ISE measurements and calibration.
- Mode Allows the user to switch between measurement modes: **Titration**, **pH**, **mV** and **ISE** mode.

2.8.2. ISE SETUP

To access the ISE Setup, press ISE Setup option key in ISE mode.



2.8.2.1. Calibration Group

Option: All Standards or Standards Group

The screenshot shows the 'ISE Setup' menu with the following options and values:

Calibration Group:	All Standards
Temperature Compensat:	Disabled
Isopotential Point:	All Standards
Edit Custom Standards:	Standards Group
Edit Standards Group:	
Calibration Reminder:	Disabled
Set Reminder Period:	Disabled
Clear Calibration	
ISE GLP Data	
Electrode Type:	Silver
Concentration Unit:	ppm
Logging Interval:	0h:00m:02s
Stability Criteria:	Medium

At the bottom of the screen, there are five buttons: Select, Escape, and three empty buttons.

All Standards Includes both standard and custom solutions.

Standards Group Includes only the standards selected by the user.

2.8.2.2. Temperature Compensation

Option: Enabled or Disabled

Note: When Temperature compensation is enabled, the isopotential point must also be set.

The screenshot shows the 'ISE Setup' menu with the following options and values:

Calibration Group:	All Standards
Temperature Compensation:	Disabled
Isopotential Point:	Standards Group
Edit Custom Standards:	Disabled
Edit Standards Group:	Enabled
Calibration Reminder:	Disabled
Set Reminder Period:	Disabled
Clear Calibration	
ISE GLP Data	
Electrode Type:	Silver
Concentration Unit:	ppm
Logging Interval:	0h:00m:02s
Stability Criteria:	Medium

At the bottom of the screen, there are five buttons: Select, Escape, and three empty buttons.

2.8.2.3. Isopotential Point (Temperature Compensation)

Option: 1.00 E^{-2} to 1.00 E^{+5} ppm

This option allows the user to set an isopotential point for the selected electrode when temperature compensation is enabled. The isopotential point is edited in ppm units only. The isopotential point will vary for different electrodes, if measurements are going to be made at several temperatures, the value should be entered if it is known.

Isopotential Point				
Enter the value for isopotential point.				
20.0 ppm				
Low limit: 1.00E-2 ppm				
High limit: 1.00E+5 ppm				
Accept	Escape	Delete Digit		Exponent

2.8.2.4. Editing Custom Standards


Option: Up to five

Edit Custom Standards				
Press <Edit> to edit selected standard.				
Press <Remove Standard> to delete the custom standard.				
Use arrows keys to select the standard.				
Remove Standard	Escape	Edit	◀	▶

1. Use the ◀ and ▶ keys to select the standard.
2. Press Remove Standard to delete the standard.
3. Press Edit to edit the selected custom standard; use the numeric keys to edit the standard.

2.8.2.5. Editing Standard Group

Option: Up to 5 standards




Use the arrow keys to select the standard to be included in / removed from the standard group.

Add or Remove ; Adds or removes the selected standard to / from standard group.

Escape ; Returns to ISE Setup menu.

2.8.2.6. Calibration Reminder

Option: Daily, Periodic, Disabled



Daily The calibration reminder will appear daily, at specified time.

Periodic The calibration reminder will appear after the set time since the last calibration has elapsed.

Disable The calibration reminder will not appear.

2.8.2.7. Setting Reminder

If Daily or Periodic option was selected for the calibration reminder, the reminder period must also be set.

For a daily reminder, the time of day can be set.

For a periodic reminder period the number of days, hours and minutes can be set.

Periodic Calibration Reminder				
Enter the time period that must be passed since the last calibration, whereafter the calibration reminder appears.				
10 days	2 hours	30 minutes		
Press Next to move to the next entry.				
Accept	Escape	Delete Digit	Next	Off

Next Moves the cursor to the next field.

Accept Saves the changes or Escape to return to the previous screen.

Off Disables the calibration reminder and return to ISE setup menu.

2.8.2.8. Clearing Calibration

This option clears the existing ISE calibration. If the calibration is cleared, a new calibration must be done in order to take measurements.

Clear Clears the previous calibration or Escape to return to the previous screen.

Clear Calibration				
Press <Clear> to clear all calibration points.				
Press <Escape> to return without clearing the calibration points.				
Clear	Escape			

2.8.2.9. ISE GLP Data

Displays the ISE calibration data

ISE GLP Data				
Analog 1				
Last Calibration:		13:42 May 24, 2018		
Slope: 100.8%		ISE: Silver		
Isopotential Point: 20.0 ppm				
1.00E-1 ppm,	0.1mV	28.1°C	A	
13:39:43 May 24, 2018				
1.00 ppm,	59.5mV	28.1°C	A	
13:40:39 May 24, 2018				
2.00 ppm,	77.6mV	28.1°C	A	
13:41:25 May 24, 2018				
10.0 ppm,	120.0mV	28.1°C	A	
13:41:45 May 24, 2018				
100 ppm,	181.0mV	28.2°C	A	
13:42:17 May 24, 2018				
Escape				

2.8.2.10. Electrode Type

Option: Ammonia, Bromide, Cadmium, Calcium, Carbon Dioxide, Chloride, Cupric, Cyanide, Fluoride, Iodide, Lead, Nitrate, Potassium, Silver, Sodium, Sulfate, Sulfide, five custom electrodes

Electrode Type				
Select a menu option.				
<div style="border: 1px solid black; padding: 5px;"> Ammonia Bromide Cadmium Calcium Carbon Dioxide Chloride Cupric Cyanide Fluoride Iodide Lead Nitrate Potassium Silver </div>				
Select	Escape	View	Page Up	Page Down

View See the ion constants (name, molar weight, electric charge / slope).


Escape Returns to the setup screen.

Ion Constants				
View Ion constants.				
<div style="border: 1px solid black; padding: 5px;"> Name: Silver Molar Weight: 107.868 g/mol Electric Charge / Slope: 1 / 59.16 </div>				
Escape				

The Ion Constants for Custom Electrodes can be modified.


2.8.2.14. ISE Significant Digits

Option: One (X), Two (XX), Three (XXX).




2.8.2.15. Stirrer Configuration

Option: Disabled, Stirrer 1



2.8.2.16. Stirring Speed

Option: 200 to 2500 RPM



2.8.3. ISE CALIBRATION

It is recommended to calibrate the instrument frequently if high accuracy is required. The instrument should also be recalibrated whenever the “Calibrate Electrode” message appears on the LCD.

Due to electrode conditioning time, the electrode must be immersed for several seconds to stabilize. The user will be guided step by step during calibration with easy-to-follow messages on the display. This will make the calibration a simple and error-free procedure.

Preparation

Pour small quantities of the standard solution into clean beakers. If possible, use plastic beakers to minimize any EMC interferences. For accurate calibration and to minimize cross-contamination, use two beakers for each standard solution: one for rinsing the electrode and one for calibration.

Note: For accurate measurements, add the appropriate ISA (Ionic Strength Adjustment) to the calibration standards.


Calibration procedure

Before calibrating, make sure that the electrode type and concentration unit has been selected in ISE Setup.

Up to a five-points calibration is possible using any combination of five standard solutions and five custom solutions.

The ISE calibration and measurement can be performed with or without temperature compensation. If the temperature compensation option is enabled, the isopotential point of the electrode must be set in ISE Setup.

1. Press from the main screen. If the instrument was calibrated before and the calibration was not cleared, the old calibration can be cleared by pressing .
2. Immerse the ISE and the temperature probe approximately 2 cm into the standard with the lowest concentration.
3. Select the standard concentration with or .
4. When the reading has stabilized, press to update the calibration. The calibration point value will be added to the Calibrated Standard list.
5. Select and repeat the procedure with all of the available standards.
6. Press to exit the calibration.




2.8.4. LOGGING

Data logging is available in ISE mode. It can be logging on demand (Manual Logging) or automatically (Interval Logging) at predefined time intervals.

The the logging report can be customized. See [2.9.3.5. Setting Up pH / mV / ISE Report](#) section for more information.

2.8.4.1. Interval Logging


The logging interval is set in the ISE Setup screen.

Press  to start the log.

The logging interval and name of logging file will be displayed on the measure screen.

To stop the automatic logging, press Stop.


2.8.4.2. Manual Logging

To manually log pH readings, press  from the ISE measurement screen.

A new record will be added to the report every time  is pressed.

2.9. AUXILIARY FUNCTIONS

2.9.1. BURETTE

To access the **Burette** screen, press  from the main titration screen.

Highlight the desired option and then press .

Burette				
Select a menu option.				
<div style="border: 1px solid black; padding: 5px;"> Prime Burette Rinse Tip Manual Dispense Purge Burette </div>				
The current pump is: Pump 1 Current burette volume is 5 mL.				
Select	Escape	Choose Pump		

 Allows you to select the desired pump for burette operations (it is only active if two pumps are connected).


Pump Setting				
Select the current pump.				
<div style="border: 1px solid black; padding: 5px;"> Pump 1 Pump 2 </div>				
Select	Escape			

2.9.1.1. Priming the Burette

Option: Up to 5

The *Prime Burette* option is used to fill the burette with titrant before starting a titration. The priming process consists of several cycles of filling and emptying the burette with titrant.

Two rinse cycles are shown in the figure below. The dispensing tube is connected on the right side and the aspiration tube on the left side.



Note: Before starting this operation, the aspiration tube must be inserted in the titrant bottle. A waste container should be placed under the dispensing tip to collect the waste solution.

To prime the burette, select Prime Burette, enter the number of rinses and press Accept.

We recommend at least three rinses to assure that the air bubbles are completely removed.

Total Burette Rinses

Enter the total number of burette rinses.

3

A minimum of three rinses is recommended.


ACCEPT	Escape	Delete Digit		
--------	--------	-----------------	--	--

2.9.1.2. Rinsing Burette Tip

A 2 mL dose of titrant will be dispensed from the burette when this operation is selected, this will eliminate any air in the dispensing tip.

2.9.1.3. Manual Dispense

Manual Dispense option allows a defined titrant volume to be dosed. Select the *Manual Dispense* option and press Select.



Use the numeric keypad to enter the volume to be dispensed.

The manual dispense volume must be between the limits shown below:

5 mL burette 0.001 to 4.750 mL

10 mL burette 0.001 to 9.500 mL

25 mL burette 0.005 to 23.750 mL


50 mL burette 0.005 to 47.500 mL

2.9.1.4. Purging the Burette

This option allows the burette to be emptied before cleaning or storing the burette. The burette is flushed twice.

Note: Before starting this operation, remove the aspiration tube from the titrant bottle.

The figures below show the steps in a purge burette operation.




2.9.2. STIRRER

The stirrer can be turned on and off by pressing stir.


During the titration process, the stirring speed can be manually adjusted using the ▲ and ▼ keys.

2.9.3. RESULTS

From the **Data Parameters** screen, you can access the following options:




2.9.3.1. Reviewing Last Analysis Report




The report contains information based on the selections made in the **Setup Titration Report** and **Setup ISE / pH / mV Report** screen.






- View Graph Review the graph.
- Print Report Print the titration report.
- Escape Return to the previous screen.
- Page Up Page Down Scroll through the pages.

2.9.3.2. Reviewing Available Reports

Up to 100 reports can be saved on the titrator. To view one of the saved reports, highlight a report and then press .




The report contains only the information selected in the **Setup Titration Report** and **Setup pH / mV / ISE Report** screens during report configuration.

-  Review the selected graph.
-  Review the selected report.
-  Print the selected report.
-  Delete the selected report.
-  Return to the previous screen.

2.9.3.3. GLP Data

Option: Up to 20 characters



- Sample Name** Allows the sample name to be recorded in each report. The sample name will increase by one, with each new titration or logging report, if the last character is a number.
- Company Name** Allows the company name to be recorded in each report.
- Operator Name** Allows the operator name to be recorded in each report.
- Electrode Name** Allows the electrode name to be recorded in each report.
- Fields 1, 2, 3** Allows any additional information to be recorded in each report.

The fields must be selected from **Setup Titration Report** screen in order to be displayed in the titration report. See [2.9.3.6. Setting Up Titration Report](#) section for more information.

2.9.3.4. Meter Information

Displays titrator configuration data.

Meter Information			
SERIAL NUMBER 931 Titrator			
Titrator Serial Number:	12133404404		
Analog Board1 Serial Number:	30134202202		
Pump 1 Serial Number:	70094513513		
Stirrer 1 Serial Number:	70091703703		
SOFTWARE VERSION			
Titrator Software Version:	v1.00		
Base Board Software Version:	v1.00		
Pump 1 Software Version:	v1.00		
Stirrer 1 Software Version:	v1.00		
Analog 1 Calibration Date:	May 22, 2018		
Escape	Print		

Titrator Serial Number	The serial number of the titrator base board.
Analog Board 1 Serial Number	The serial number of the analog board.
Pump 1 (or 2) Serial Number	The serial number of the connected pump.
Titrator Software Version	The current software version installed on the titrator.
Base Board Software Version	The current software version present on the base board of the titrator.
Pump 1 (or 2) Software Version	The current software version for the pump.
Analog 1 Calibration Date	Manufacturer calibration date of the analog board.

Note: If more than 1 year elapsed from the calibration date of the analog board 1, the message "Analog 1 Calibration Due" will appear on the main screen. The analog board need to be recalibrated.

2.9.3.5. Setting Up pH / mV / ISE Report


Customize a unique report to record the pH, mV, and ISE measurements. An asterisk means that it will be included in the report.

Setup pH/mV/ISE Report				
Select fields to be saved in the report.				
* Result and Units				
* Potential				
* Temperature and Units				
* Date and Time				
* Calibration Data				
Sample Name				
Company Name				
Operator Name				
Electrode Name				
Field 1				
Field 2				
Field 3				
Software Versions				
Serial Numbers				
Select	Escape	Save Report	Page Up	Page Down

- Select Adds the highlighted information to the report.
- Unselect Removes the highlighted information from the report.
- Escape Returns to the Data Parameter Screen. Report is not updated.
- Save Report Updates the report with the selected items. Report previously saved will not be updated.
- Page Up Page Down Scrolls through the options.

2.9.3.6. Setting Up Titration Report

Customize a unique report to record the titration results. An asterisk means that it will be included in the titration report.



- Select Adds the highlighted information to the report.
- Unselect Removes the highlighted information from the report.
- Escape Returns to the Data Parameter Screen. Report is not updated.
- Save Report Updates the report with the selected items. Report previously saved will not be updated.
- Page Up Page Down Scrolls through the options.

2.10. MAINTENANCE & PERIPHERALS

The 25-mL burette included with the titrator exceeds the ISO 8655 standard for accurate delivery of liquids by a motor-driven piston burette.


2.10.1. BURETTE MAINTENANCE

2.10.1.1. Burette Assembly

The burette is delivered with a 25-mL syringe inside and with all of the accessories mounted. See [2.1. Setup](#) section for more information. The burette assembly consists of a rigid housing which holds the glass syringe, a 3-way valve and titrant tubing.

2.10.1.2. Changing the Burette

Remove the burette from the pump assembly by sliding it forward and then slide the new burette into place.




2.10.1.3. Disassembling the Burette

The aspiration and the dispensing tubes have fittings and tube protectors. The aspiration tube is mounted on the left side and the dispensing tube is mounted on the right side of the burette.

To remove the dispensing tube and the aspiration tube follow these steps:


1. Remove the blue tube protector (A) by sliding it off the clear titrant tubing.
2. Remove the tube lock (B) from the burette holder.
3. Turn the fitting (C) counterclockwise to remove it from the burette holder.
4. Slide the clear titrant tubing through the fitting.



2.10.1.4. Assembling the Burette

To attach the dispensing tube and the aspiration tube, follow these steps:

1. Insert the flat-shaped end of the titrant tubing into the valve outlet (A) and screw the fitting clockwise to tighten. The highest of the 9 cuts should be vertical in the final position.
2. Bend the tube up into the vertical position to enter the highest cut of the fitting (C).
3. Replace the tube lock fitting (D).
4. Replace the blue tube protector (E) by sliding it over the clear titrant tubing, the protector will sit in the tube lock fitting.



2.10.1.5. Cleaning the Burette

To clean the burette, follow these steps:


1. If the burette is filled with titrant, remove the aspiration tube from the titrant bottle and purge burette. See [2.9.1.4. Purging the Burette](#) section for more information.
2. Insert the aspiration tube into cleaning solution, deionized water or titrant solvent.
3. Go through two cycles of filling and emptying the burette. See [2.9.1.4. Purging the Burette](#) section for more information.
4. During second cycle, remove the aspiration tube from the cleaning solution, deionized water or solvent and allow the air to replace the liquid in the burette. This will clean the aspiration tube.

If this simple cleaning procedure is not adequate, continue with these steps:

1. Remove the burette assembly from the pump.
2. Remove the dispensing and aspiration tubes. Clean them separately or insert new ones.
3. Remove the protective cap from the bottom of the burette assembly by using the burette removal tool.
4. Remove the syringe from the burette assembly by unscrewing it with your fingers.
5. Extract the piston from the syringe.
6. Clean both the piston and the syringe with appropriate cleaning solution. Rinse with deionized water.
7. Remove the excess liquid.

Warning: Avoid contacting the titrant with bare hands. Avoid spilling titrant. Clean the external side of the syringe and piston to remove aggressive chemicals. Do not touch the white PTFE part of the piston or internal walls of the burette with bare hands or greasy materials.

8. Reinsert the piston into the syringe.
9. Reinsert the syringe by screwing it in the valve with your fingers.
10. Reinsert the protective cap to the bottom of the burette assembly. Carefully position the cap into the burette.
11. Slide the burette into the burette stand. Notice the position of the piston shaft to the pump couple.
12. Priming the burette three times with new titrant is recommended.



2.10.1.6. Burette Preparation (Titrant Filling)

Before starting a titration, the burette must be properly filled with titrant in order to obtain an accurate and repeatable result. To fill the burette, follow the next steps and recommendations:

1. If necessary, clean the burette and make sure it is empty.
2. From the main screen press .
3. Highlight *Prime Burette* option and press .
4. Enter the number of times the burette needs to be rinsed (minimum three rinses are recommend allowing air bubbles to be evacuated).
5. Press .

To avoid the presence of the air bubbles inside the burette, make sure to have a continuous liquid flow inside the burette. A little air just above the liquid level at the first filling is normal. The next filling will evacuate all of the air; no air will be left in the valve. Sometimes during this process, slight finger tapping on the tubes helps remove any residual air bubbles.

If air bubbles are still present:

1. Remove the aspiration tube from the titrant bottle.
2. Repeat burette preparation procedure.
3. If this is not successful, clean the burette again.


2.10.2. PERIPHERALS

Warning: Connection or disconnection of POWER, PUMP ASSEMBLY, PRINTER or RS232 INTERFACE must only be done when titrator and external devices are turned off.

2.10.2.1. Connecting to a Printer

A variety of parallel printers can be connected to the parallel port of the titrator using a DB25 cable.


Warning: The titrator and the external printer must be both turned off before they are connected.



2026 Titrator image shown above.














2.10.2.2. Connecting an External PC Keyboard

This connection allows the use of an external PS/2 PC Keyboard in addition to the titrator's keypad.




2026 Titrator image shown above.

The correspondence between the titrator's keypad and the United States 101-type external keyboard is detailed below:

External PC Keyboard (United States 101)	Titrator Keypad
Function key F-1	
Function key F-2	
Function key F-3	
Function key F-4	
Function key F-5	Option key 1 (from left to right)
Function key F-6	Option key 2 (from left to right)
Function key F-7	Option key 3 (from left to right)
Function key F-8	Option key 4 (from left to right)
Function key F-9	Option key 5 (from left to right)
Function key F-10	
Arrow key: Up	
Arrow key: Down	
Arrow key: Left	
Arrow key: Right	
Page Up	
Page Down	
Numeric keys: 0 to 9	
Enter	
Alphanumeric keys	Allow alphanumeric entries

2.10.2.3. Connecting to a Computer


The titrator can be connected to a computer using a USB cable. HI900 PC application needs to be installed on the PC.



2026 Titrator image shown above.

To connect the PC to the titrator follow the steps below:

1. Connect the cable to the USB port on the rear panel of the titrator.
2. Connect the cable to the USB port on the PC.



The HI900 PC application allows the transfer of methods and reports between the titrator and PC. See [2.3.12. USB Link with PC](#) section for more information.

2.11. ACCESSORIES

2.11.1. SOLUTIONS

2.11.1.1. pH Calibration Buffers

Ordering Information	Description
HI7001M	pH 1.68 buffer solution, 230 mL
HI7001L	pH 1.68 buffer solution, 500 mL
HI7004M	pH 4.01 buffer solution, 230 mL
HI7004L	pH 4.01 buffer solution, 500 mL
HI7006M	pH 6.86 buffer solution, 230 mL
HI7006L	pH 6.86 buffer solution, 500 mL

Ordering Information	Description
HI7007M	pH 7.01 buffer solution, 230 mL
HI7007L	pH 7.01 buffer solution, 500 mL
HI7009M	pH 9.18 buffer solution, 230 mL
HI7009L	pH 9.18 buffer solution, 500 mL
HI7010M	pH 10.01 buffer solution, 230 mL
HI7010L	pH 10.01 buffer solution, 500 mL

2.11.1.2. pH Calibration Buffers in FDA Approved Bottle

Ordering Information	Description
HI8004L	pH 4.01 buffer solution, 500 mL
HI8006L	pH 6.86 buffer solution, 500 mL
HI8007L	pH 7.01 buffer solution, 500 mL

Ordering Information	Description
HI8009L	pH 9.18 buffer solution, 500 mL
HI8010L	pH 10.01 buffer solution, 500 mL

2.11.1.3. pH Technical Calibration Buffers

Ordering Information	Description
HI5016	pH 1.68 buffer solution, 500 mL
HI5003	pH 3.00 buffer solution, 500 mL
HI5004	pH 4.01 buffer solution, 500 mL
HI5068	pH 6.86 buffer solution, 500 mL

Ordering Information	Description
HI5007	pH 7.01 buffer solution, 500 mL
HI5091	pH 9.18 buffer solution, 500 mL
HI5010	pH 10.01 buffer solution, 500 mL
HI5124	pH 12.45 buffer solution, 500 mL

2.11.1.4. pH Millesimal Calibration Buffers

Ordering Information	Description
HI6016	pH 1.679 buffer solution, 500 mL
HI6003	pH 3.000 buffer solution, 500 mL
HI6004	pH 4.010 buffer solution, 500 mL
HI6004-01	pH 4.010 buffer solution, 1 L
HI6068	pH 6.862 buffer solution, 500 mL
HI6007	pH 7.010 buffer solution, 500 mL

Ordering Information	Description
HI6007-01	pH 7.010 buffer solution, 1 L
HI6091	pH 9.177 buffer solution, 500 mL
HI6010	pH 10.010 buffer solution, 500 mL
HI6010-01	pH 10.010 buffer solution, 1 L
HI6124	pH 12.450 buffer solution, 500 mL

2.11.1.5. Electrode Cleaning Solutions

Ordering Information	Description
HI7061M	General purpose cleaning solution, 230 mL
HI7061L	General purpose cleaning solution, 500 mL
HI7073M	Protein cleaning solution, 230 mL
HI7073L	Protein cleaning solution, 500 mL
HI7074M	Inorganic cleaning solution, 230 mL
HI7074L	Inorganic cleaning solution, 500 mL
HI7077M	Oil & fat cleaning solution, 230 mL
HI7077L	Oil & fat cleaning solution, 500 mL

2.11.1.6. Electrode Cleaning Solutions in FDA Approved Bottle

Ordering Information	Description
HI8061L	General purpose solution, 500 mL
HI8073L	Protein cleaning solution, 500 mL
HI8077L	Oil & fat cleaning solution, 500 mL

2.11.1.7. Electrode Storage Solutions

Ordering Information	Description
HI70300M	Storage solution, 230 mL
HI70300L	Storage solution, 500 mL

2.11.1.8. Electrode Storage Solutions in FDA Approved Bottle

Ordering Information	Description
HI80300M	Storage solution, 230 mL
HI80300L	Storage solution, 500 mL

2.11.1.9. Electrode Refill Electrolyte Solutions

Ordering Information	Description
HI7071	3.5 M KCl with AgCl reference electrolyte solution, 30 mL
HI7072	1 M Potassium nitrate electrode fill solution
HI7075	1.7 M Potassium nitrate, 0.7 M potassium chloride electrode fill solution
HI7076	1 M Sodium chloride electrode fill solution
HI7078	0.5 M Ammonium sulfate electrode fill solution
HI7082	3.5 M KCl reference electrolyte solution, 30 mL

2.11.1.10. Electrode Refill Electrolyte Solutions in FDA Approved Bottle

Ordering Information	Description
HI8071	3.5 M KCl with AgCl reference electrolyte solution, 30 mL
HI8082	3.5 M KCl reference electrolyte solution, 30 mL

2.11.1.11. ORP Pretreatment Solutions

Ordering Information	Description
HI7091L	Reducing pretreatment solution, 500 mL
HI7092M	Oxidizing pretreatment solution, 230 mL
HI7092L	Oxidizing pretreatment solution, 500 mL

2.11.1.12. Titration Reagents

Ordering Information	Description
HI70429	0.05 M Silver nitrate titration reagent, 1 L
HI70433	0.01 N Stabilized iodine titration reagent, 1 L
HI70439	0.1 M Sodium thiosulfate titration reagent, 1 L
HI70440	0.02 N Stabilized iodine titration reagent, 1 L
HI70441	0.04 N Stabilized iodine titration reagent, 1 L
HI70448	0.02 M Silver nitrate titration reagent, 1 L
HI70449	0.02 M EDTA titration reagent, 1 L
HI70455	0.01 N Sodium hydroxide titration reagent, 1 L
HI70456	0.1 N Sodium hydroxide titration reagent, 1 L
HI70457	1 N Sodium hydroxide titration reagent, 1 L
HI70458	0.01 M Sulfuric acid titration reagent, 1 L
HI70459	0.05 M Sulfuric acid titration reagent, 1 L
HI70462	0.01 N Hydrochloric acid titration reagent, 1 L
HI70463	0.1 N Hydrochloric acid titration reagent, 1 L
HI70464	1 N Hydrochloric acid titration reagent, 1 L

2.11.1.13. Ion-Selective Electrode Calibration Standards

Ordering Information	Description
HI4001-01	0.1 M Ammonia standard
HI4001-02	100 ppm Ammonia standard (as N)
HI4001-03	1000 ppm Ammonia standard (as N)
HI4002-01	0.1 M Bromide standard
HI4003-01	0.1 M Cadmium standard
HI4004-01	0.1 M Calcium standard
HI4005-01	0.1 M Carbon dioxide standard
HI4005-03	1000 ppm Carbon dioxide standard (as CaCO ₃)
HI4007-01	0.1 M Chloride standard
HI4007-02	100 ppm Chloride standard
HI4007-03	1000 ppm Chloride standard
HI4008-01	0.1 M Cupric standard
HI4010-01	0.1 M Fluoride standard
HI4010-02	100 ppm Fluoride standard
HI4010-03	1000 ppm Fluoride standard
HI4011-01	0.1 M Iodide standard
HI4012-01	0.1 M Lead standard
HI4012-21	0.1 M Sulfate standard
HI4013-01	0.1 M Nitrate standard
HI4013-02	100 ppm Nitrate standard
HI4013-03	1000 ppm Nitrate standard
HI4014-01	0.1 M Potassium standard
HI4015-01	0.1 M Silver standard

2.11.2. SENSORS

2.11.2.1. pH Electrodes

Ordering Information	Description
HI1043B	Glass-body, double junction, refillable, combination pH electrode Use: strong acid and base, paint and solvents
HI1053B	Glass-body, triple ceramic, conic shape, refillable, combination pH electrode Use: emulsions, fats and creams, soil and semi-solids samples
HI1083B	Glass-body, micro, Viscolene, nonrefillable, combination pH electrode Use: biotechnology and micro titration
HI1131B	Glass-body, double junction, refillable, combination pH electrode Use: general purpose
HI1330B	Glass-body, semimicro, single junction, refillable, combination pH electrode Use: laboratory, vials, and test tubes
HI1331B	Glass-body, semimicro, single junction, refillable, combination pH electrode Use: flasks
HI1230B	Plastic-body (PEI), double junction, gel-filled, combination pH electrode Use: general purpose
HI2031B	Glass-body, conical tip, refillable, combination pH electrode Use: dairy and semi-solid products
HI1332B	Plastic-body (PEI), double junction, refillable, combination pH electrode Use: chemicals, field applications and quality control testing
FC100B	Plastic-body (PVDF), double junction, refillable, combination pH electrode Use: cheese
FC200B	Plastic-body (PVDF), single junction, conical tip, non-refillable Viscolene electrolyte, combination pH electrode Use: milk, yogurt, dairy products, and semi-solid foods
FC210B	Glass-body, double junction, conical tip, non-refillable Viscolene electrolyte, combination pH electrode Use: milk, yogurt, and cream
FC220B	Glass-body, single junction, refillable, combination pH electrode Use: milk, yogurt, cream, sauce, and fruit juices
FC911B	Plastic-body (PVDF), double junction, refillable, combination pH electrode Use: sauce, juices, dairy products and other liquid or slurry forms of food
HI1413B	Glass-body, single junction, flat tip, non-refillable Viscolene electrolyte, combination pH electrode Use: surfaces, skin, leather, paper, and emulsions

2.11.2.2. ORP Electrodes

Ordering Information	Description
HI3131B	Glass-body, refillable, combination platinum ORP electrode Use: laboratories and general purpose
HI3230B	Plastic-body (PEI), gel-filled, combination platinum ORP electrode Use: municipal water and quality control
HI4430B	Plastic-body (PEI), gel-filled, combination gold ORP electrode Use: oxidants and ozone

2.11.2.3. Half-Cell Electrodes

Ordering Information	Description
HI5311	Glass-body, silver / silver chloride (Ag / AgCl) reference half-cell electrode, double junction, refillable with 4 mm banana plug with 1 m (3.3') cable Use: general purpose with wide temperature range
HI5315	Plastic-body (PEI), double junction, silver / silver chloride (Ag / AgCl) reference half-cell electrode, refillable with 4 mm plug with 1 m (3.3') cable. Use: Ion-Selective Electrodes
HI5412	Glass-body, single Calomel reference half-cell electrode, refillable with 4mm plug with 1 m (3.3') cable Use: general purpose with constant temperature range

















2.11.2.4. Ion-Selective Electrodes


Ordering Information	Description
HI4101	Ammonia ion selective electrode
HI4002 / HI4102	Bromide ion selective electrode
HI4003 / HI4103	Cadmium ion selective electrode
HI4004 / HI4104	Chloride ion selective electrode
HI4105	Carbon dioxide ion selective electrode
HI4007 / HI4107	Chloride ion selective electrode
HI4008 / HI4108	Cupric ion selective electrode
HI4009 / HI4109	Cyanide ion selective electrode
HI4010 / HI4110	Fluoride ion selective electrode
HI4011 / HI4111	Iodide ion selective electrode
HI4012 / HI4112	Lead ion selective electrode
HI4013 / HI4113	Nitrate ion selective electrode
HI4014 / HI4114	Potassium ion selective electrode
HI4015 / HI4115	Silver / Sulfide ion selective electrode
FC300B	Sodium electrode

2.11.2.5. Temperature Sensor


Ordering Information	Description
HI7662-TW	Temperature probe with 1 m (3.3') paneled cable

2.11.3. TITRATOR COMPONENTS


	<p>HI930100 Pump assembly</p>		<p>HI900205 5 mL Syringe</p>
	<p>Burette with: HI930105 - 5 mL syringe HI930110 - 10 mL syringe HI930125 - 25 mL syringe HI930150 - 50 mL syringe</p>		<p>HI900210 10 mL Syringe</p>
	<p>HI930191 Blank support</p>		<p>HI900225 25 mL Syringe</p>
	<p>HI930190 Blank burette support</p>		<p>HI900250 50 mL Syringe</p>
	<p>HI900942 Tool for burette cap removal</p>		<p>HI930301 Overhead stirrer & 3 propellers</p>
	<p>HI900260 3 Way valve</p>		<p>HI930302 Replacement propellers (3 pcs.)</p>
	<p>HI900270 Aspiration tube with fitting and protection tube</p>		<p>HI930303 High chemical resistance propellers (3 pcs.)</p>
	<p>HI930280 Dispensing tube with dispensing tip, fitting, protection tube and tube guide</p>		<p>HI930320 Overhead holder support rod</p>




HI930310
Overhead electrode holder




HI900947
Power adapter (european plug)




HI930330
Titrant bottle cap assembly




HI900946
Power adapter (USA plug)




HI930311
Electrode adapter for
overhead holder




HI920013
USB-C to USB-A cable




HI930315
Titrant bottle holder




HI930900U
USB Storage device



HI7662-TW
Temperature probe



HI930801
Instruction manual binder



HI900945
Shorting cap

PART 3. APPLICATIONS

HI0001EN — 0.1N SODIUM HYDROXIDE TITRANT CONCENTRATION

Description

Method for the standardization (titer determination) of 0.1N Sodium Hydroxide (NaOH) titrant solution against Potassium Hydrogen Phthalate (KHP). The results are expressed in **N (eq/L)**.

Reference

AOAC Official Methods of Analysis, Official Method 936.16

Electrode

- [HI1131B](#) Combination pH Electrode
- [HI7662-T](#) Temperature Probe

Reagents

- [HI70456](#) 0.1N Sodium Hydroxide (1 L)
- [HI70401](#) Potassium Hydrogen Phthalate (20 g)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7071](#) Electrode Fill Solution (30 mL x 4)
- [HI7004L](#) pH 4.01 Buffer Solution (500 mL)
- [HI7007L](#) pH 7.01 Buffer Solution (500 mL)
- [HI7010L](#) pH 10.01 Buffer Solution (500 mL)
- [HI740036P](#) 100 mL Plastic Beaker (10 pcs)
- Analytical Balance with 0.0001 g resolution

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.1N sodium hydroxide ([HI70456](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press from the main screen. Use the arrow keys to highlight [HI0001EN 0.1N Sodium Hydroxide](#) and press .

Electrode Preparation

- Press from the main screen, if necessary select the analog board and press .
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.



Sample Preparation

- Crush approximately 3 grams of potassium hydrogen phthalate ([HI70401](#)) and dry it for 2 hours at 120 °C. Cool to room temperature in a desiccator.
- Place a clean 100 mL plastic beaker on the analytical balance.
- Zero the balance.
- Carefully weigh approximately 0.20 grams of dried potassium hydrogen phthalate into the beaker. Ensure that all of the potassium hydrogen phthalate is on the bottom of the beaker.
- Record the exact weight of the sample once the balance has stabilized with an accuracy of 0.0001 grams.
- Remove the beaker from the balance and add deionized water to the 50 mL mark on the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.






- Press . You will be prompted to enter the weight of the analyte (weight of potassium hydrogen phthalate). Use the numeric keypad to enter the exact weight and press  to start the analysis.

Note: Ensure that the potassium hydrogen phthalate dissolves completely during the pre-titration stir time. Erroneous results may occur if the sample does not dissolve completely prior to titration. If necessary the pre-titration stir time can be increased.

- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **N (eq/L) of sodium hydroxide**.
- Remove the pH electrode, temperature probe and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

Note: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1N sodium hydroxide titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1N sodium hydroxide.
- Press  from the main screen.
- Using the arrow keys, highlight Titrant Conc. and press .
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press .
- Press  to exit the **View/Modify Method** screen. Use the arrow keys to highlight **Save Method** and press .

METHOD PARAMETERS

Name: 0.1N Sodium Hydroxide
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titration Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.030 mL
 Max Vol: 0.500 mL
 delta E: 4.500 mV
 End Point Mode: pH 1EQ point, 1st Der
 Recognition Options:
 Threshold: 500 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 5.000 mL
 Pre-Titration Stir Time: 60 sec
 Measurement Mode: Signal Stability
 delta E: 0.3 mV
 delta t: 2 sec
 Min wait: 3 sec
 Max wait: 30 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Weight
 Dilution Option: Disabled
 Titrant Name: 0.1N NaOH
 Analyte Size: 0.20000 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Stdz. Titrant by Weight
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.200 g
 mw of standard: 204.23 g/mol
 Titrant/Standard: 1.000 eq/mol

$$\frac{\text{eq}}{\text{L}} \text{NaOH} = \frac{0.200 \times 1.000}{204.23 \times V (\text{L})}$$

RESULTS

Titration Report
 Method Name: 0.1N Sodium Hydroxide
 Time & Date: 17:03 Jun 07, 2018
 Report ID: Ti_00053

Titration Results
 Method Name: 0.1N Sodium Hydroxide
 Time & Date: 17:03 Jun 07, 2018
 Analyte Size: 0.20920 g
 End Point Volume: 10.215 mL
 pH Equivalence Point: 8.394
 Result: 0.10027 N(eq/L)
 Initial & Final pH: 4.173 to 9.570
 Titration Duration: 6:25 [mm:ss]
 Titration went to Completion
 Analyst Signature: _____

HI0002EN — 0.1N HYDROCHLORIC ACID TITRANT CONCENTRATION

Description

Method for the standardization (titer determination) of 0.1N Hydrochloric Acid (HCl) titrant solution against Sodium Hydroxide (NaOH). The results are expressed in **N (eq/L)**.

Reference

AOAC Official Methods of Analysis, Official Method 936.15

Electrode

- HI1131B Combination pH Electrode
- HI7662-T Temperature Probe

Reagents

- HI70463 0.1N Hydrochloric Acid (1 L)
- HI70456 0.1N Sodium Hydroxide (1 L)
- HI70436 Deionized Water (1 gal)

Accessories

- HI70300L Storage Solution (500 mL)
- HI7071 Electrode Fill Solution (30 mL x 4)
- HI7004L pH 4.01 Buffer Solution (500 mL)
- HI7007L pH 7.01 Buffer Solution (500 mL)
- HI7010L pH 10.01 Buffer Solution (500 mL)
- HI740036P 100 mL Plastic Beaker (10 pcs)
- 10 mL Class A Volumetric Pipette

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.1N hydrochloric acid (HI70463) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press from the main screen. Use the arrow keys to highlight HI0002EN 0.1N Hydrochloric Acid and press .

Electrode Preparation

- Press from the main screen, if necessary select the analog board and press .
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.

SAMPLE PREPARATION

- Use a Class A volumetric pipette to transfer exactly 10.00 mL of 0.1N sodium hydroxide (HI70456) to a clean 100 mL beaker
- Add deionized water to the 50 mL mark on the beaker.

ANALYSIS

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press . The titrator start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **N (eq/L) of hydrochloric acid**.

- Remove the pH electrode, temperature probe and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

Note: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1N hydrochloric acid titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1N hydrochloric acid.
- Press Method Options from the main screen.
- Using the arrow keys, highlight Titrant Conc. and press Select.
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press Accept.
- Press Escape to exit the **View/Modify Method** screen. Use the arrow keys to highlight Save Method and press Select.

METHOD PARAMETERS

Name: 0.1N Hydrochloric Acid
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.030 mL
 Max Vol: 0.500 mL
 delta E: 6.000 mV
 End Point Mode: pH 1EQ point, 1st Der
 Recognition Options:
 Threshold: 500 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 5.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 3 sec
 Max wait: 15 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Volume
 Dilution Option: Disabled
 Titrant Name: 0.1N HCl
 Analyte Size: 10.0000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 15.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Stdz. Titrant by Volume
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Standard volume: 10.000 mL
 Standard conc.: 0.100 eq/L

$$\frac{\text{eq HCl}}{\text{L}} = \frac{10.000 \times 0.100}{\text{V (L)} \times 1000}$$

RESULTS

Titration Report

Method Name: 0.1N Hydrochloric Acid
 Time & Date: 14:55 July 30, 2018
 Report ID: Ti_00002

Titration Results

Method Name: 0.1N Hydrochloric Acid
 Time & Date: 14:55 July 30, 2018
 Analyte Size: 10.000 mL
 End Point Volume: 9.979 mL
 pH Equivalence Point: 5.059
 Result: 0.10020 N(eq/L)
 Initial & Final pH: 12.135 to 4.989
 Titration Duration: 2:45 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI0003EN — 0.1M SODIUM THIOSULFATE TITRANT CONCENTRATION

Description

Method for the standardization (titer determination) of 0.1M Sodium Thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$) titrant solution against Potassium Iodate (KIO_3). The results are expressed in **M (mol/L)**.

Reference

Standard Methods for the Examination of Water and Wastewater 19th Edition, Method 4500-Cl B

Electrode

- [HI3131B](#) Combination ORP Electrode

Reagents

- [HI70439](#) 0.1M Sodium Thiosulfate (1 L)
- [HI70407](#) Potassium Iodate (20 g)
- [HI70425](#) 16% Sulfuric Acid (500 mL)
- [HI70468](#) Potassium Iodide (35 g)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7071](#) Electrode Fill Solution (30 mL x 4)
- [HI740036P](#) 100 mL Plastic Beakers (10 pcs)
- Analytical Balance 0.0001 g
- 100 mL Class A Volumetric Flask
- 10 mL Class A Volumetric Pipette

Device Preparation

- Connect the ORP electrode to the titrator.
- Install a 25 mL burette filled with 0.1M sodium thiosulfate ([HI70439](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press from the main screen. Use the arrow keys to highlight [HI0003EN 0.1M Sodium Thiosulfate](#) and press .

Electrode Preparation

- Prepare the ORP electrode according to the procedure in the manual.



Sample Preparation

- Crush approximately 2 grams of potassium iodate ([HI70407](#)) and dry it for 2 hours at 120 °C. Cool to room temperature in a desiccator.
- Carefully weigh approximately 0.35 grams of dried potassium iodate.
- Record the exact weight of the sample once the balance has stabilized with an accuracy of 0.0001 grams.
- Carefully transfer the salt to a 100 mL Class A volumetric flask. Add approximately 80 mL of deionized water, and mix to dissolve. Once the salt is completely dissolved bring the flask to volume with deionized water, mix well.
- Use a Class A volumetric pipette to transfer exactly 10.00 mL of the solution to a clean 100 mL plastic beaker.
- Add deionized water to the 50 mL mark on the beaker.
- Add 5.00 mL of 16% sulfuric acid ([HI70425](#)) and 1.5 grams of potassium iodide ([HI70468](#)) to the beaker.

Analysis






- Place the beaker under the stirrer assembly and lower it to immerse the ORP electrode and stirrer. Ensure that the reference junction of the ORP electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press . You will be prompted to enter the weight of the analyte (weight of potassium iodate). Use the numeric keypad to enter the exact weight and press  to start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **M (mol/L) of sodium thiosulfate**.
- Remove the ORP electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

Note: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1M sodium thiosulfate titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1M sodium thiosulfate.
- Press  from the main screen.
- Using the arrow keys, highlight Titrant Conc. and press .
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press .
- Press  to exit the **View/Modify Method** screen. Use the arrow keys to highlight Save Method and press .

METHOD PARAMETERS

```
Name: 0.1M Sodium Thiosulfate
Method Revision: 3.0
Stirrer Configuration:
  Stirrer: Stirrer 1
  Stirring Speed: 1400 RPM
Pump Configuration:
  Titrant Pump: Pump 1
Dosing Type: Dynamic
  Min Vol: 0.030 mL
  Max Vol: 0.600 mL
  delta E: 6.500 mV
End Point Mode: mV 1EQ point, 1st Der
Recognition Options:
  Threshold: 50 mV/mL
  Range: NO
  Filtered Derivatives: NO
Pre-Titration Volume: 5.000 mL
Pre-Titration Stir Time: 0 sec
Measurement Mode: Signal Stability
  delta E: 0.3 mV
  delta t: 2 sec
  Min wait: 2 sec
  Max wait: 20 sec
Electrode Type: ORP
Blank Option: No Blank
Calculations: Stdz. Titrant by Weight
Dilution Option: Enabled
  Final Dilution Volume: 100.000 mL
  Aliquot Volume: 10.000 mL
Titrant Name: 0.1M Na2S2O3
Analyte Size: 0.35000 g
Analyte Entry: Manual
Maximum Titrant Volume: 15.000 mL
Potential Range: -2000.0 to 2000.0 mV
Volume/Flow Rate: 25 mL/50.0 mL/min
Signal Averaging: 1 Reading
Significant Figures: XXXXX
```

CALCULATIONS

```
Calculations: Stdz. Titrant by Weight
Titrant units: M (mol/L)
Titrant volume dosed: V (L)
Standard weight: 0.350 g
Dilution Factor: 0.100
  Final Dilution volume: 100.000 mL
  Aliquot Volume: 10.000 mL
mw of standard: 214.00 g/mol
Titrant/Standard: 6.000 mol/mol

$$\frac{\text{mol}}{\text{L}} \text{Na}_2\text{S}_2\text{O}_3 = \frac{0.350 \times 0.10 \times 6.0}{214.00 \times V(\text{L})}$$

```

RESULTS

```
Titration Report
Method Name: 0.1M Sodium Thiosulfate
Time & Date: 17:03 Jun 07, 2018
Report ID: Ti_00073
```

```
Titration Results
Method Name: 0.1M Sodium Thiosulfate
Time & Date: 17:03 Jun 07, 2018
Analyte Size: 0.35020 g
End Point Volume: 9.635 mL
mV Equivalence Point: 233.0
Result: 0.10191 M (mol/L)
Initial & Final mV: 361.8 to 173.4
Titration Duration: 2:51 [mm:ss]
Titration went to Completion
```

Analyst Signature: _____

HI0010EN — 0.1M FERROUS AMMONIUM SULFATE TITRANT CONCENTRATION

Description

Method for the standardization (titer determination) of 0.1M Ferrous Ammonium Sulfate (FAS) titrant solution against Potassium Dichromate ($K_2Cr_2O_7$). The results are expressed in **M (mol/L)**.

Reference

Standard Methods for the Examination of Water and Wastewater 21st Edition, Method 5220B

Electrode

- HI3131B Combination ORP Electrode

Reagents

- HI70444 25% Sulfuric Acid
- HI70436 Deionized Water (1 gal)
- Ferrous Ammonium Sulfate (ACS Grade)
- Potassium Dichromate (ACS Grade)

Accessories

- HI70300L Storage Solution (500 mL)
- HI7071 Electrode Fill Solution (30 mL x 4)
- HI740036P 100 mL Plastic Beakers (10 pcs)
- Analytical Balance with 0.0001 g resolution
- 100 mL Class A Volumetric Flask
- 500 mL Class A Volumetric Flask
- 10 mL Class A Volumetric Pipette

Titrant Preparation

- Carefully weigh 19.607 grams of ferrous ammonium sulfate.
- Carefully transfer the salt to a 500 mL Class A volumetric flask. Add approximately 300 mL of deionized water, and mix to dissolve.
- Add 40.00 mL of 25% sulfuric acid (HI70444) to the flask. Invert the solution to mix.
- Allow the flask to return to room temperature.
- Bring the flask to volume with deionized water, mix well.

Device Preparation

- Connect the ORP electrode to the titrator.
- Install a 25 mL burette filled with 0.1M ferrous ammonium sulfate on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press Select Method from the main screen. Use the arrow keys to highlight HI0010EN 0.1M FAS and press Select.

Electrode Preparation

- Prepare the ORP electrode according to the procedure in the manual.

Sample Preparation

- Carefully weigh approximately 0.49 grams of dried potassium dichromate.
- Record the exact weight of the sample once the balance has stabilized with an accuracy of 0.0001 grams.
- Carefully transfer the salt to a 100 mL Class A volumetric flask. Add approximately 80 mL of deionized water, and mix to dissolve. Once the salt is completely dissolved bring the flask to volume with deionized water, mix well.
- Use a Class A volumetric pipette to transfer exactly 10.00 mL of the solution to a clean 100 mL plastic beaker.
- Add 25.00 mL of 25% sulfuric acid (HI70444) to the beaker.
- Add deionized water to the 50 mL mark on the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the electrode and stirrer. Ensure that the reference junction of the ORP electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press . You will be prompted to enter the weight of the analyte (weight of potassium dichromate). Use the numeric keypad to enter the exact weight and press to start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **M (mol/L) of ferrous ammonium sulfate**.
- Remove the ORP electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

Note: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.1M ferrous ammonium sulfate titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.1M ferrous ammonium sulfate.
- Press from the main screen.
- Using the arrow keys, highlight Titrant Conc. and press .
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press .
- Press to exit the **View/Modify Method** screen. Use the arrow keys to highlight Save Method and press .

METHOD PARAMETERS

Name: 0.1M FAS
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.030 mL
 Max Vol: 0.500 mL
 delta E: 4.500 mV
 End Point Mode: mV 1EQ point, 1st Der
 Recognition Options:
 Threshold: 35 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 5.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 0.5 mV
 delta t: 3 sec
 Min wait: 2 sec
 Max wait: 20 sec
 Electrode Type: ORP
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Weight
 Dilution Option: Enabled
 Final Dilution Volume: 100.000 mL
 Aliquot Volume: 10.000 mL
 Titrant Name: 0.1M FAS
 Analyte Size: 0.49000 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL

Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Stdz. Titrant by Weight
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.490 g
 Dilution Factor: 0.100
 Final Dilution volume: 100.000 mL
 Aliquot Volume: 10.000 mL
 mw of standard: 294.18 g/mol
 Titrant/Standard: 6.000 mol/mol

$$\frac{\text{mol}}{\text{L}} \text{ FAS} = \frac{0.490 \times 0.10 \times 6.0}{294.18 \times V(\text{L})}$$

RESULTS

Titration Report
 Method Name: 0.1M FAS
 Time & Date: 15:59 August 1, 2018
 Report ID: Ti_00015

Titration Results
 Method Name: 0.1M FAS
 Time & Date: 15:59 August 1, 2018
 Analyte Size: 0.491 g
 End Point Volume: 9.879 mL
 mV Equivalence Point: 667.4
 Result: 0.10137 M (mol/L)
 Initial & Final mV: 791.3 to 598.0
 Titration Duration: 3:05 [mm:ss]
 Titration went to Completion
 Analyst Signature: _____

HI0200EN – 0.02M SILVER NITRATE TITRANT CONCENTRATION

Description

Method for the standardization (titer determination) of 0.02M Silver Nitrate (AgNO_3) titrant solution against Sodium Chloride (NaCl). The results are expressed in **M (mol/L)**.

Reference

AOAC Official Methods of Analysis, Official Method 941.18

Electrode

- [HI4115](#) Silver/Sulfide Combination ISE



Reagents

- [HI70448](#) 0.02M Silver Nitrate (1 L)
- [HI70406](#) Sodium Chloride (20 g)
- [HI70427](#) 1.5M Nitric Acid Solution (500 mL)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI7072](#) Electrode Fill Solution (4 x 30 mL)
- Analytical Balance with 0.0001 g resolution
- 150 mL Glass Beaker
- 100 mL Class A Volumetric Flask
- 5 mL Class A Volumetric Pipette

Device Preparation

- Connect the Silver/Sulfide electrode to the titrator.
- Install a 25 mL burette filled with 0.02M silver nitrate ([HI70448](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press  from the main screen. Use the arrow keys to highlight [HI0200EN 0.02M Silver Nitrate](#) and press .

Electrode Preparation

- Prepare the Silver/Sulfide electrode according to the procedure in the manual.



Sample Preparation

- Crush approximately 2 grams of sodium chloride ([HI70406](#)) and dry it for 2 hours at 140 °C. Cool to room temperature in a desiccator.
- Weigh 0.20 g of dried sodium chloride with an accuracy of 0.0001 g. Transfer the salt to a 100 mL volumetric flask. Add approximately 80 mL of distilled water and mix. Dissolve completely before bringing to volume.
- Use a Class A volumetric pipette to transfer exactly 5.00 mL of prepared standard solution to a 150 mL glass beaker and add distilled water to the 100 mL mark on the beaker.
- Add 10.00 mL of 1.5M nitric acid ([HI70427](#)) to the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the Silver/Sulfide electrode and stirrer. Ensure that the reference junction of the electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press . You will be prompted to enter the weight of the analyte (weight of sodium chloride). Use the numeric keypad to enter the exact weight and press  to start the analysis.

- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **M (mol/L) of silver nitrate**.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

Note: For improved accuracy, repeat this procedure a minimum of three times and calculate the average value.

For methods utilizing 0.02M silver nitrate titrant solution, follow the steps below to enter the titer/standardized value.

- Select the method utilizing 0.02M silver nitrate.
- Press Method Options from the main screen.
- Using the arrow keys, highlight Titrant Conc. and press Select.
- Use the numeric keypad to enter the standardized (titer) value of the titrant then press Accept.
- Press Escape to exit the **View/Modify Method** screen. Use the arrow keys to highlight Save Method and press Select.

METHOD PARAMETERS

Name: 0.02M Silver Nitrate
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.030 mL
 Max Vol: 0.500 mL
 delta E: 8.000 mV
 End Point Mode: mV 1EQ point, 1st Der
 Recognition Options:
 Threshold: 100 mV/mL
 Range: NO
 Filtered Derivatives: YES
 Pre-Titration Volume: 6.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 20 sec
 Electrode Type: Silver/Sulfide
 Blank Option: No Blank
 Calculations: Stdz. Titrant by Weight
 Dilution Option: Enabled
 Final Dilution Volume: 100.000 mL
 Aliquot Volume: 5.000 mL
 Titrant Name: 0.02M AgNO₃
 Analyte Size: 0.20000 g
 Analyte Entry: Manual
 Maximum Titrant Volume: 15.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Stdz. Titrant by Weight
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Standard weight: 0.200 g
 Dilution Factor: 0.05
 Final Dilution volume: 100.000 mL
 Aliquot Volume: 5.000 mL
 mw of standard: 58.440 g/mol
 Titrant/Standard: 1.000 mol/mol

$$\frac{\text{mol}}{\text{L}} \text{AgNO}_3 = \frac{0.200 \times 0.05 \times 1.0}{58.440 \times V(\text{L})}$$

RESULTS

Titration Report
 Method Name: 0.02M Silver Nitrate
 Time & Date: 15:52 August 1, 2018
 Report ID: Ti_00037

Titration Results
 Method Name: 0.02M Silver Nitrate
 Time & Date: 15:52 August 1, 2018
 Analyte Size: 0.1923 g
 End Point Volume: 9.065 mL
 mV Equivalence Point: 273.1
 Result: 0.01815 M (mol/L)
 Initial & Final mV: 146.9 to 291.0
 Titration Duration: 2:21 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1004EN – ALKALINITY OF WATER

0 to 2500 mg/L CaCO₃, pH 4.5 Endpoint

Description

Method for the determination of total (methyl red) alkalinity in water by titration of a sample to pH 4.5. The results are expressed in **mg/L (ppm) as calcium carbonate**.

For the determination of phenolphthalein alkalinity, set the endpoint to pH 8.3.

Reference

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 2320B

Electrode

- [HI1131B](#) Combination pH Electrode
- [HI7662-T](#) Temperature Probe

Reagents

- [HI70463](#) 0.1N Hydrochloric Acid (1 L)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7082](#) Electrode Fill Solution (4 x 30 mL)
- [HI7004L](#) pH 4.01 Buffer Solution (500 mL)
- [HI7007L](#) pH 7.01 Buffer Solution (500 mL)
- [HI7010L](#) pH 10.01 Buffer Solution (500 mL)
- [HI740036P](#) 100 mL Plastic Beaker (10 pcs)
- 50 mL Class A Volumetric Pipette

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.1N hydrochloric acid ([HI70463](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- For the determination of the exact concentration of the 0.1N hydrochloric acid, follow [HI0002EN 0.1N Hydrochloric Acid Titrant Concentration](#).
- Press Select Method from the main screen. Use the arrow keys to highlight [HI1004EN Alkalinity of Water](#) and press Select.

Electrode Preparation

- Press Mode from the main screen, if necessary select the analog board and press pH.
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.


Sample Preparation

- Use a Class A volumetric pipette to transfer exactly 50.00 mL of sample to a clean 100 mL plastic beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature sensor and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press . The titrator will start the analysis.
- At the end of the titration, when pH 4.50 is reached, "Titration Completed" will appear with the result. The result is expressed in mg/L as calcium carbonate.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

METHOD PARAMETERS

Name: Alkalinity of Water
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.050 mL
 Max Vol: 0.500 mL
 delta E: 5.000 mV
 End Point Mode: Fixed 4.500 pH
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 20 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.1N HCl
 Titrant Conc.: 0.1000 N(eq/L)
 Analyte Size: 50.000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 25.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Sample Calc. by Volume
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Final result units: mg/L
 Titrant Conc.: 0.1000 N(eq/L)
 Sample/Titrant: 0.500 mol/eq
 mw of standard: 100.09 g/mol
 Sample Volume: 50.000 mL

$$\frac{\text{mg}}{\text{L}} \text{CaCO}_3 = \frac{V(\text{L}) \times 1000 \times 0.10 \times 0.5 \times 100.09 \times 1000}{50.00}$$

RESULTS

Titration Report

Method Name: Alkalinity of Water
 Time & Date: 14:36 August 1, 2018
 Report ID: Ti_00036

Titration Results

Method Name: Alkalinity of Water
 Time & Date: 14:36 August 1, 2018
 Analyte Size: 50.000 mL
 End Point Volume: 9.336 mL
 pH Fixed End Point: 4.500
 Result: 934.44 mg/L
 Initial & Final pH: 10.232 to 4.419
 Titration Duration: 3:23 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1005EN – ACIDITY OF WATER

0 to 2500 mg/L, pH 8.3 Endpoint

Description

Method for the determination of total (phenolphthalein) acidity in water by titration of a sample to pH 8.3. The results are expressed in **mg/L (ppm) as calcium carbonate**.

For the determination of methyl orange acidity, set the endpoint to pH 3.7.

Reference

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 2310B

Electrode

- [HI1131B](#) Combination pH Electrode
- [HI7662-T](#) Temperature Probe

Reagents

- [HI70456](#) 0.1N Sodium Hydroxide (1 L)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7082](#) Electrode Fill Solution (4 x 30 mL)
- [HI7004L](#) pH 4.01 Buffer Solution (500 mL)
- [HI7007L](#) pH 7.01 Buffer Solution (500 mL)
- [HI7010L](#) pH 10.01 Buffer Solution (500 mL)
- [HI740036P](#) 100 mL Plastic Beaker (10 pcs)
- 50 mL Class A Volumetric Pipette

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.1N sodium hydroxide ([HI70456](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- For the determination of the exact concentration of the 0.1N sodium hydroxide, follow [HI0001EN 0.1N Sodium Hydroxide Titrant Concentration](#).
- Press from the main screen. Use the arrow keys to highlight [HI1005EN Acidity in Water](#) and press .

Electrode Preparation

- Press from the main screen, if necessary select the analog board and press .
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.


Sample Preparation

- Use a Class A volumetric pipette to transfer exactly 50.00 mL of sample to a clean 100 mL plastic beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature sensor and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: The dispensing tip should be slightly submerged in the sample.

- Press , the titrator will start the analysis.
- At the end of the titration, when pH 8.30 is reached, "Titration Completed" will appear with the result. The result is expressed in mg/L as calcium carbonate.
- Remove the pH electrode, temperature probe and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

METHOD PARAMETERS

Name: Acidity of Water
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.050 mL
 Max Vol: 0.500 mL
 delta E: 5.000 mV
 End Point Mode: Fixed 8.300 pH
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 20 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.1N NaOH
 Titrant Conc.: 0.1000 N(eq/L)
 Analyte Size: 50.000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 25.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Sample Calc. by Volume
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Final result units: (mg/L)
 Titrant Conc.: 0.1000 N(eq/L)
 Sample/Titrant: 0.500 mol/eq
 mw of standard: 100.09 g/mol
 Sample Volume: 50.000 mL

$$\frac{\text{mg}}{\text{L}} \text{CaCO}_3 = \frac{V(\text{L}) \times 1000 \times 0.10 \times 0.5 \times 100.09 \times 1000}{50.0}$$

RESULTS

Titration Report

Method Name: Acidity of Water
 Time & Date: 14:54 August 1, 2018
 Report ID: Ti_00023

Titration Results

Method Name: Acidity of Water
 Time & Date: 14:54 August 1, 2018
 Analyte Size: 50.000 mL
 End Point Volume: 5.879 mL
 pH Fixed End Point: 8.300
 Result: 588.43 (mg/L)
 Initial & Final pH: 2.465 to 8.398
 Titration Duration: 3:42 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1007EN – CHLORIDE IN WATER

0 to 150 ppm (mg/L)

Description

Method for the determination of chloride in water. The results are expressed as **ppm (mg/L) as Chloride**.

REFERENCE

Standard Methods for the Examination of Water and Wastewater 21st edition, Method 4500-Cl

Electrode

- [HI4115](#) Silver/Sulfide Combination ISE



Reagents

- [HI70448](#) 0.02M Silver Nitrate (1 L)
- [HI70427](#) 1.5M Nitric Acid Solution (500 mL)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI7072](#) Electrode Fill Solution (4 x 30 mL)
- 150 mL Glass Beaker
- 100 mL Class A Volumetric Pipette
- 10 mL Class A Volumetric Pipette

DEVICE PREPARATION

- Connect the Silver/Sulfide electrode to the titrator.
- Install a 25 mL burette filled with 0.02M silver nitrate ([HI70448](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- For the determination of the exact concentration of the 0.02M Silver Nitrate, follow [HI0200EN 0.02M Silver Nitrate Titrant Concentration](#).
- Press  from the main screen. Use the arrow keys to highlight [HI1007EN Chloride in Water](#) and press .

Electrode Preparation

- Prepare the Silver/Sulfide electrode according to the procedure in the manual.


Sample Preparation

- Use a class A volumetric pipette to transfer exactly 100.00 mL of sample to a clean 150 mL beaker.
- Add 10.00 mL of 1.5M nitric acid ([HI70427](#)) to the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the electrode and stirrer. Ensure that the reference junction of the electrode is 5 to 6 mm below the surface. If necessary add extra deionized water.

Note: *The dispensing tip should be slightly submerged in the sample.*

- Press , the titrator will start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **ppm (mg/L) of chloride**.
- Remove the electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

METHOD PARAMETERS

Name: Chloride in Water
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.030 mL
 Max Vol: 0.500 mL
 delta E: 5.000 mV
 End Point Mode: mv 1EQ point, 1st Der
 Recognition Options:
 Threshold: 100 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 20 sec
 Electrode Type: Silver/Sulfide
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.02M AgNO3
 Titrant Conc.: 2.0000E-2 M (mol/L)
 Analyte Size: 100.000 mL
 Analyte Entry: Manual
 Maximum Titrant Volume: 25.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Sample Calc. by Volume
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Final result units: (mg/L)
 Titrant Conc.: 2.0000E-2 M (mol/L)
 Sample/Titrant: 1.000 mol/mol
 mw of sample: 35.453 g/mol
 Sample Volume: 100.000 mL

$$\frac{\text{mg}}{\text{L}} = \frac{V(\text{L}) \times 1000 \times 0.02 \times 1.0 \times 35.45 \times 1000}{100.0}$$

RESULTS

Titration Report

Method Name: Chloride in Water
 Time & Date: 15:11 August 1, 2018
 Report ID: Ti_00052

Titration Results

Method Name: Chloride in Water
 Time & Date: 15:11 August 1, 2018
 Analyte Size: 100.000 mL
 End Point Volume: 4.781 mL
 mV Fixed End Point: 280.3
 Result: 33.897 ppm (mg/L)
 Initial & Final mV: 94.8 to 298.5
 Titration Duration: 1:24 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1008EN – NEUTRALIZATION WITH SULFURIC ACID

0 to 200 meq/L

Description

Method for the determination of strong or weak base concentration by titration of a sample to the equivalence point with sulfuric acid. The results are expressed as **meq/L**.

Electrode3

- [HI1131B](#) Combination pH Electrode
- [HI7662-T](#) Temperature Probe

Reagents

- [HI70459](#) 0.05M Sulfuric Acid (1 L)
- [HI70436](#) Deionized Water (1 gal)

Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7082](#) Electrode Fill Solution (4 x 30 mL)
- [HI7004L](#) pH 4.01 Buffer Solution
- [HI7007L](#) pH 7.01 Buffer Solution
- [HI7010L](#) pH 10.01 Buffer Solution
- [HI740036P](#) 100 mL Plastic Beaker (10 pcs)
- 10 mL Class A Volumetric Pipette

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.05M sulfuric acid ([HI70459](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- For the determination of the exact concentration of the 0.05M sulfuric acid, follow [HI0103EN 0.05M Sulfuric Acid Titrant Concentration](#).
- Press from the main screen. Use the arrow keys to highlight [HI1008EN Neutralization w/H2SO4](#) and press .

Electrode Preparation

- Press from the main screen, if necessary select the analog board and press .
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.

Sample Preparation

- Use a class A volumetric pipette to transfer exactly 10.00 mL of sample to a clean 100 mL plastic beaker.
- Add deionized water to the 50 mL mark on the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface.

Note: *The dispensing tip should be slightly submerged in the sample.*

- Press , the titrator will start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **meq/L**.
- Remove the pH electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

METHOD PARAMETERS

Name: Neutralization w/ H2SO4
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.050 mL
 Max Vol: 0.500 mL
 delta E: 20.000 mV
 End Point Mode: pH 1EQ point, 1st Der
 Recognition Options:
 Threshold: 50 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 15 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.05M H2SO4
 Titrant Conc.: 5.0000E-2 M (mol/L)
 Analyte Size: 10.000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 20.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Sample Calc. by Volume
 Titrant units: M (mol/L)
 Titrant volume dosed: V (L)
 Final result units: meq/L
 Titrant Conc.: 5.0000E-2 M (mol/L)
 Sample/Titrant: 2.000 eq/mol
 Sample Volume: 10.000 mL

$$\frac{\text{meq}}{\text{L}} = \frac{V(\text{L}) \times 1000 \times 0.05 \times 2.0 \times 1000}{10.0}$$

RESULTS

Titration Report

Method Name: Neutralization w/ H2SO4
 Time & Date: 09:46 August 1, 2018
 Report ID: Ti_00027

Titration Results

Method Name: Neutralization w/ H2SO4
 Time & Date: 09:46 August 1, 2018
 Analyte Size: 10.000 mL
 End Point Volume: 9.562 mL
 mV Equivalence Point: 7.966
 Result: 95.620 meq/L
 Initial & Final pH: 11.655 to 6.248
 Titration Duration: 1:24 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1009EN – NEUTRALIZATION WITH SODIUM HYDROXIDE

0 to 200 meq/L

Description

Method for the determination of strong or weak acid concentration by titration of a sample to the equivalence point with sodium hydroxide. The results are expressed as **meq/L**.

Electrode

- [HI1131B](#) Combination pH Electrode
- [HI7662-T](#) Temperature Probe



Reagents

- [HI70456](#) 0.1N Sodium Hydroxide (1 L)
- [HI70436](#) Deionized Water (1 gal)



Accessories

- [HI70300L](#) Storage Solution (500 mL)
- [HI7082](#) Electrode Fill Solution (4 x 30 mL)
- [HI7004L](#) pH 4.01 Buffer Solution
- [HI7007L](#) pH 7.01 Buffer Solution
- [HI7010L](#) pH 10.01 Buffer Solution
- [HI740036P](#) 100 mL Plastic Beaker (10 pcs)
- 10 mL Class A Volumetric Pipette

Device Preparation

- Connect the pH electrode and temperature probe to the titrator.
- Install a 25 mL burette filled with 0.1N sodium hydroxide ([HI70456](#)) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- For the determination of the exact concentration of the 0.1N sodium hydroxide, follow [HI0001EN 0.1N Sodium Hydroxide Titrant Concentration](#)
- Press  from the main screen. Use the arrow keys to highlight [HI1009EN Neutralization w/NaOH](#) and press .

Electrode Preparation

- Press  from the main screen, if necessary select the analog board and press .
- Calibrate the electrode using pH 4.01, 7.01 and 10.01 buffers. Refer to the instruction manual for calibration procedure.


Sample Preparation

- Use a class A volumetric pipette to transfer exactly 10.00 mL of sample to a clean 100 mL plastic beaker.
- Add deionized water to the 50 mL mark on the beaker.

Analysis

- Place the beaker under the stirrer assembly and lower it to immerse the pH electrode, temperature probe and stirrer. Ensure that the reference junction of the pH electrode is 5 to 6 mm below the surface.

Note: *The dispensing tip should be slightly submerged in the sample.*

- Press , the titrator will start the analysis.
- At the end of the titration, after detection of the equivalence point, "Titration Completed" will appear with the result. The result is expressed in **meq/L**.
- Remove the pH electrode and stirrer from the sample and rinse them thoroughly with deionized water.
- Record the result.

METHOD PARAMETERS

Name: Neutralization w/ NaOH
 Method Revision: 3.0
 Stirrer Configuration:
 Stirrer: Stirrer 1
 Stirring Speed: 1400 RPM
 Pump Configuration:
 Titrant Pump: Pump 1
 Dosing Type: Dynamic
 Min Vol: 0.050 mL
 Max Vol: 0.500 mL
 delta E: 20.000 mV
 End Point Mode: pH 1EQ point, 1st Der
 Recognition Options:
 Threshold: 50 mV/mL
 Range: NO
 Filtered Derivatives: NO
 Pre-Titration Volume: 0.000 mL
 Pre-Titration Stir Time: 0 sec
 Measurement Mode: Signal Stability
 delta E: 1.0 mV
 delta t: 2 sec
 Min wait: 2 sec
 Max wait: 15 sec
 Electrode Type: pH
 Blank Option: No Blank
 Calculations: Sample Calc. by Volume
 Dilution Option: Disabled
 Titrant Name: 0.1N NaOH
 Titrant Conc.: 0.1000 N(eq/L)
 Analyte Size: 10.000 mL
 Analyte Entry: Fixed
 Maximum Titrant Volume: 20.000 mL
 Potential Range: -2000.0 to 2000.0 mV
 Volume/Flow Rate: 25 mL/50.0 mL/min
 Signal Averaging: 1 Reading
 Significant Figures: XXXXX

CALCULATIONS

Calculations: Sample Calc. by Volume
 Titrant units: N (eq/L)
 Titrant volume dosed: V (L)
 Final result units: meq/L
 Titrant Conc.: 5.0000E-2 M (mol/L)
 Sample/Titrant: 0.1000 N(eq/L)
 Sample Volume: 10.000 mL

$$\frac{\text{meq}}{\text{L}} = \frac{V(\text{L}) \times 1000 \times 0.1 \times 1.0 \times 1000}{10.0}$$

RESULTS

Titration Report

Method Name: Neutralization w/ NaOH
 Time & Date: 10:29 August 2, 2018
 Report ID: Ti_00017

Titration Results

Method Name: Neutralization w/ NaOH
 Time & Date: 10:29 August 2, 2018
 Analyte Size: 10.000 mL
 End Point Volume: 15.970 mL
 pH Equivalence Point: 8.431
 Result: 159.70 meq/L
 Initial & Final pH: 2.675 to 10.316
 Titration Duration: 3:20 [mm:ss]
 Titration went to Completion

Analyst Signature: _____

HI1011EN – TROUBLESHOOTING 1

Description

Method for verifying the dosing and potentiometric signal accuracy of the titrator. This method should be used to troubleshoot a titrator equipped with a 25 mL burette. The titrator dispenses a 20.00 mL pre-titration volume, waits 20 seconds and dispenses an additional 20.00 mL dose, bringing the total volume to 40.00 mL. This procedure can also be used to check the stability of the mV and temperature channels. The dosing accuracy of the 25 mL burette is ± 0.025 mL ($\pm 0.1\%$ of the full volume). If the results are not correct, check all fittings for leakage, and burette and tubing for air bubbles. Repeat the measurement.

Reference

ISO/TC 48/SC1N 380E and 383E: "Piston and/or Plunger Operated Volumetric Apparatus"

Accessories

- HI762000C 0 °C Temperature Key
- HI762070C 70 °C Temperature Key
- HI70436 Deionized Water (1 gal)
- HI7662-T Temperature Probe
- Shorting Cap
- Narrow Neck Beaker
- Analytical Balance with 0.0001g resolution

Device Preparation

- Connect the shorting cap to the BNC socket on Analog Board 1
- Install a 25 mL burette filled with room temperature deionized water (HI70436) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press Select Method from the main screen. Use the arrow keys to highlight HI1011EN Troubleshooting 1 and press Select.

Large Dose Dispensing Procedure

- Add a small amount of deionized water to a narrow neck beaker.
- Place the narrow neck beaker on an analytical balance and zero the balance.
- Place the dosing tip through the neck of the beaker. Take care not to immerse it in the liquid during dispensing and not to touch the beaker walls.
- Press start stop.
- Write down the exact weight displaced on the balance after each dose.
- This procedure can be repeated on pump 2.

Other burette sizes can be checked using the following settings, see instruction manual for accuracy:

Burette Volume	Pre-titration Volume	Max. Titrant Volume
5 mL	4.000 mL	8.000 mL
10 mL	8.000 mL	16.000 mL

METHOD PARAMETERS

Name:	Troubleshooting 1	Measurement Mode:	Timed Increment
Method Revision:	3.0	Time interval:	20 sec
Stirrer Configuration:		Electrode Type:	Shorting Cap
Stirrer:	Stirrer 1	Blank Option:	No Blank
Stirring Speed:	0 RPM	Calculations:	No Formula (mL only)
Pump Configuration:		Titrant Name:	DI Water
Titrant Pump:	Pump 1	Maximum Titrant Volume:	40.000 mL
Dosing Type:	Linear - 20.000 mL	Potential Range:	-2000.0 to 2000.0 mV
End Point Mode:	Fixed 10.0 mV	Volume/Flow Rate:	25 mL/50.0 mL/min
Pre-Titration Volume:	20.000 mL	Signal Averaging:	1 Reading
Pre-Titration Stir Time:	0 sec	Significant Figures:	XXXXXX

Calculations

$$V = m * \frac{1}{\rho} * \left(1 + \frac{\rho_{\text{air}}}{\rho_{\text{L}}} - \frac{\rho_{\text{air}}}{\rho_{\text{std}}} \right)$$

V	Volume of measure mass of water (mL)
m	Measure mass of water (g)
ρ_{L}	Density of dispensed water (g/mL)
ρ_{air}	Density of ambient air (g/mL)
ρ_{std}	Density of calibration standard weight (g/mL)

Alternative Calculations

If the actual values of the above parameters are not accessible the following equation can be used:

$$V = M * F$$

V	Volume of measured mass of water (mL)
F	Transformation factor

The values from the table have been calculated by correcting the air and water density with temperature, assuming the density of dry air $\rho_{\text{air}} = 0.0012$ g/mL and density of calibration steel standard weigh $\rho_{\text{STD}} = 8$ g/mL.

Temperature (°C)	Factor
17.0	1.002290
18.0	1.002467
19.0	1.002654
20.0	1.002853
21.0	1.003061
22.0	1.003282
23.0	1.003512

Temperature (°C)	Factor
24.0	1.003752
25.0	1.004002
26.0	1.004261
27.0	1.004531
28.0	1.004809
29.0	1.005097
30.0	1.005395

Temperature Channel Fast Check Procedure

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Connect the HI762000C 0 °C temperature key to the RCA socket (temperature sensor input) on Analog Board 1.
- On the main screen select Mode, if necessary select the analog board and press mV.
- The titrator should display ATC 0.0 ± 0.4 °C with no fluctuations or drift.
- Connect the HI762070C 70 °C temperature key to the RCA socket (temperature sensor input) on Analog Board 1.
- The titrator should display ATC 70.0 ± 0.4 °C with no fluctuations or drift.
- This procedure can be repeated on analog board 2.

Temperature & mV Channel Logging Procedure

- Connect the shorting cap to the BNC socket on Analog Board 1.
- Connect the HI762000C 0°C temperature key to the RCA socket (temperature sensor input) on Analog Board 1.
- On the main screen select Mode, if necessary select the analog board and press mV.
- Press mV Setup and use the arrow keys to highlight Logging Interval. Set the logging interval to 15 seconds and press Accept. Press Escape to return to the main screen.
- Press the results key and use the arrow keys to highlight Setup pH/mV/ISE Report, press Select.
- Select Potential and Temperature and Units. All other fields should be unselected.
- Press Save Report to return to the Data Parameters screen.
- Press Escape to return to the main screen.
- Once on the main screen press Start Log to start the automatic log.
- Let the log run for about 10 minutes. Press Stop Log to stop the automatic log.
- Press results, use the arrow keys to highlight Review Last Analysis Report, and press Select.
- The mV column should display 0.0 ± 0.1 mV and the temperature column should display 0.0 °C ± 0.4 °C.
- This procedure can be repeated using the HI762070C 70 °C temperature key and on analog board 2.

HI1012EN – TROUBLESHOOTING 2

Description

Method for verifying the dosing of the titrator. This method should be used to troubleshoot a titrator equipped with a 25 mL burette. The titrator dispenses a 10.00 mL pre-titration volume, waits 20 seconds and dispenses an additional 0.5 mL dose twenty times, waiting 20 seconds between each dose, bringing the total volume to 20 mL. This procedure can also be used to check the stirrer functionality. The dosing accuracy of the 25 mL burette is ± 0.025 mL ($\pm 0.1\%$ of the full volume).

If the results are not correct, check all fittings for leakage, and burette and tubing for air bubbles. Repeat the measurement.

Reference

ISO/TC 48/SC1N 380E and 383E: "Piston and/or Plunger Operated Volumetric Apparatus"

Accessories

- HI762000C 0 °C Temperature Key
- HI70436 Deionized Water (1 gal)
- HI7662-T Temperature Probe
- Shorting Cap
- Narrow Neck Beaker
- Analytical Balance with 0.0001g resolution

Device Preparation

- Connect the shorting cap to the BNC socket on Analog Board 1
- Install a 25 mL burette filled with room temperature deionized water (HI70436) on pump one and verify that no air bubbles are present in the burette or tubing. If necessary prime the burette until all the air has been removed completely.
- Press Select Method from the main screen. Use the arrow keys to highlight HI1012EN Troubleshooting 2 and press Select.

Small Dose Dispensing Procedure

- Add a small amount of deionized water to a narrow neck beaker. By doing this the air space in the beaker will be vapor-saturated minimizing evaporation.
- Place the narrow neck beaker on an analytical balance and zero the balance.
- Place the dosing tip through the neck of the beaker. Take care not to immerse it in the liquid during dispensing and not to touch the beaker walls.
- Press start stop.
- Write down the exact weight displaced on the balance after each dose.
- This procedure can be repeated on pump 2.

Other burette sizes can be checked using the following settings, see instruction manual for accuracy:

Burette Volume	Pre-titration Volume	Max. Titrant Volume
5 mL	4.000 mL	8.000 mL
10 mL	8.000 mL	16.000 mL

METHOD PARAMETERS

Name:	Troubleshooting 2	Measurement Mode:	Timed Increment
Method Revision:	3.0	Time interval:	10 sec
Stirrer Configuration:		Electrode Type:	Shorting Cap
Stirrer:	Stirrer 1	Blank Option:	No Blank
Stirring Speed:	0 RPM	Calculations:	No Formula (mL only)
Pump Configuration:		Titrant Name:	DI Water
Titrant Pump:	Pump 1	Maximum Titrant Volume:	20.000 mL
Dosing Type:	Linear - 0.500 mL	Potential Range:	-2000.0 to 2000.0 mV
End Point Mode:	Fixed 10.0 mV	Volume/Flow Rate:	25 mL/50.0 mL/min
Pre-Titration Volume:	10.000 mL	Signal Averaging:	1 Reading
Pre-Titration Stir Time:	0 sec	Significant Figures:	XXXXXX

Calculations

$$V = m * \frac{1}{\rho} * \left(1 + \frac{\rho_{\text{air}}}{\rho_{\text{L}}} - \frac{\rho_{\text{air}}}{\rho_{\text{std}}} \right)$$

V	Volume of measure mass of water (mL)
m	Measure mass of water (g)
ρ_{L}	Density of dispensed water (g/mL)
ρ_{air}	Density of ambient air (g/mL)
ρ_{std}	Density of calibration standard weight (g/mL)

Alternative Calculations

If the actual values of the above parameters are not accessible the following equation can be used:

$$V = M * F$$

V	Volume of measured mass of water (mL)
F	Transformation factor

The values from the table have been calculated by correcting the air and water density with temperature, assuming the density of dry air $\rho_{\text{air}} = 0.0012$ g/mL and density of calibration steel standard weigh $\rho_{\text{STD}} = 8$ g/mL.

Temperature (°C)	Factor
17.0	1.002290
18.0	1.002467
19.0	1.002654
20.0	1.002853
21.0	1.003061
22.0	1.003282
23.0	1.003512
24.0	1.003752
25.0	1.004002
26.0	1.004261
27.0	1.004531
28.0	1.004809
29.0	1.005097
30.0	1.005395

Stirring Speed Fast Check Procedure

- On the main screen select Mode, if necessary select the analog board and press mV.
- Press mV Setup and use the arrow keys to highlight Stirrer Configuration. Use the arrow keys to highlight Stirrer 1. Press Accept.
- Use the arrow keys to highlight String Speed. Use the numeric keypad to enter 200 rpms then press Accept.
- Press Escape to exit the mV Setup screen.
- From the main screen, press stir, use the up arrow key to increase the stir speed slowly to 2500 rpms.
- Check that the propeller continues to increase speed, following the commands.
- This procedure can be repeated on stirrer 2.

PART 4. TITRATION THEORY

4.1. TITRATION THEORY

4.1.1. INTRODUCTION

A titration is a quantitative, volumetric procedure used in analytical chemistry to determine the concentration of an analyte (the species being measured) in solution. The concentration of the analyte is determined by slowly adding a titrant to the solution. As the titrant is added, a chemical reaction occurs between the titrant and the analyte.

Titration reactions are relatively fast, simple reactions that can be expressed using a chemical equation. The titration reaction continues as the titrant is added until all of the analyte is consumed and the analyte reacts completely and quantitatively with the titrant.

The point at which all of the analyte has been reacted is called the equivalence point, also known as the theoretical or stoichiometric endpoint. This point is accompanied by an abrupt physical change in the solution, which sharply defines the endpoint of the reaction. The physical change associated with the titration endpoint can be produced by the titrant, or an indicator, and can be detected visually or by physical measurements.

Titration cannot be used to determine the quantity of all analytes.

The chemical reaction between the titrant and analyte must fulfill four requirements:

- Must be fast and occur within approximately one second after the titrant is added
- Must go to completion
- Must have well-known stoichiometry (reaction ratios)
- A convenient endpoint or inflection point

Titration provides many advantages over alternative methods; they are highly precise, quickly performed and require relatively simple apparatus and instrumentation.

4.1.2. USES OF TITRATIONS

- Acid content of plant effluents, food (e.g.: cheese and wine), plating and etching baths, petroleum products, drugs
- Base content of fertilizer (containing ammonia), bleach, minerals
- Hardness in water
- Metal content of alloys, minerals, ores, clays, waters, plating baths, paints, paper, plant materials, biological fluids, petroleum products
- Moisture content in foodstuffs, petrochemicals, plastics, pharmaceutical products
- Redox reagent concentrations such as available chlorine in potable water, peroxide, traces of oxidants and reductants in food, reductants in high temperature or high pressure boiler water, vitamin analysis

4.1.3. ADVANTAGES & DISADVANTAGES

Advantages of titration as an analytical technique:

- More precise results than many instrumental methods, such as measurement by electrode, the accuracy of the measurement is up to 0.1%
- Simple methods, reasonable capital costs and easy training
- Suitability to measure major components of a mixture or product
- Automation can reduce time and labor spent on each analysis

Disadvantages of titration as an analytical technique:

- The time it takes to prepare standards and titrants
- Good technique is required to achieve precise results (training and practice required)
- Not suitable for determining trace or minor components of a mixture or product
- Limited dynamic range, it may require additional sample preparation (dilution) and repeat analyses

4.2. TYPES OF TITRATIONS

4.2.1. TITRATIONS ACCORDING TO THE MEASUREMENT METHOD

4.2.1.1. Amperometric Titrations

An amperometric titration is performed by placing two electrodes (typically a metal ion-selective electrode and a reference electrode) into the sample solution and keeping the potential of the metal electrode at a selected voltage. The current that flows, due to the oxidation or reduction of a reactant or product, is plotted vs. volume of titrant to provide the titration curve and locate the equivalence point. Changes in the current are due to changes in the concentration of a particular species (being oxidized or reduced at the electrode).

Generally, the reaction between the analyte and titrant forms a new species. Depending on the titration, the reactants are electroactive and the products are not or vice-versa. Amperometric titration curves look like two straight lines intersecting at the equivalence point, this is due to the change in the electroactivity of the solution. **Figure 1A**, Amperometric titrations, shows an active analyte and non-reactive titrant. **Figure 1B** and **1D**, Amperometric titrations, shows a nonreactive analyte and a reactive titrant. **Figure 1C**, Amperometric titrations, shows a reactive analyte and titrant.

Many metal ions can be amperometrically titrated using a precipitation, complexation or redox reaction. Some metal ions and species that can be determined in this manner include silver, barium, halides, potassium, magnesium, palladium, molybdate, sulfate, tungstate, zinc, bismuth, cadmium, fluoride, indium, thallium, iodine and gold.




Figure 1: Amperometric titrations

4.2.1.2. Potentiometric Titrations

Potentiometric titrations are done by measuring the voltage across the solution using an electrode system. An electrode system consists of an indicator electrode and a reference electrode. As titrant is added, the variations in the potential of the indicator electrode, with respect to the reference electrode, are monitored to show the progress of the titration.

Potentiometry is the measurement of a potential under conditions of zero current flow. The measured potential can then be used to determine the analytical quantity of interest, generally a component concentration of the analyte solution. The potential that develops in the electrochemical cell is the result of the free energy change that would occur if the chemical phenomena were to proceed until the equilibrium condition has been satisfied.

There are many types of titrations where potentiometry can be used, e.g. pH electrodes for acid-base titrations, platinum ORP electrodes in redox titrations, ion-selective electrodes, such as chloride or fluoride for a specific ion titration, and silver electrodes for argentometric (silver-based) titrations.

In **Figure 2A**, Potentiometric titrations, the pH of the solution is plotted against the volume of titrant. In **Figure 2B**, Potentiometric titrations, the electrode potential is plotted against the volume of titrant.




Figure 2: Potentiometric titrations

4.2.1.3. Spectrophotometric Titrations

The name comes from the method used to detect the endpoint of the titration, not its chemistry. Highly colored indicators that change color during the course of the titration are available for many titrations. More accurate data on the titration curve can be obtained if the light absorption is monitored instrumentally using a light source, a simple monochromator and a photodetector, rather than visually determining the color or light absorption change. Light absorption by either an indicator or by one of the reactants or products can be used to monitor the titration.

In **Figure 3A**, Spectrophotometric titrations, the absorption of a metal-indicator complex is being monitored. The absorption is constant while the metal is complexed by the ethylenediaminetetraacetic acid (EDTA) titrant. The metal indicator complex was stripped, causing a sharp break in the titration curve. The point where all the metal is complexed and stripped from the indicator is the equivalence point. This point is marked by "e.p." on the graph.

In **Figure 3B**, Spectrophotometric titrations, the metal complex is being measured while being titrated with EDTA. The new complex being formed is not colored and does not absorb light. The extrapolated intersection of the two lines determines the equivalence point.




Figure 3: Spectrophotometric titrations

4.2.2. TITRATIONS ACCORDING TO THE REACTION TYPE

4.2.2.1. Acid-Base Titrations

Acid-base titrations are the most common type of titrations. They are based upon a reaction between an acid and a base, a stoichiometric neutralization or the exchange of protons. Virtually all acid-base titrations are carried out using a strong acid or a strong base as the titrant. The endpoint of a titration carried out with a weak acid or a weak base would be difficult to detect due to a small change in pH at the equivalence point.

Chemical indicators can be used to determine the endpoint. The indicator will change color to signify that the end of the titration has been reached. The color of the indicator is dependent upon the concentration of ions in the solution. An acid-base indicator is composed of a conjugate weak acid-weak base pair, where the two forms exhibit different colors depending on the pH of the solution. For an indicator, the acid ionization constant (K_a) is usually written as:

$$K_a = \frac{[H_3O^+][In^-]}{[HIn]}$$

HIn is the acid form of the indicator and In^- is the base form. At the center of the change region, the ratio of $[In^-]$ to $[HIn]$ is one, $[H_3O^+] = K_a$ and $pH = pK_a$. The color change region is usually ± 1 pH unit around this point.

Table 1, Aqueous acid-base chemical indicators, contains a list of aqueous acid-base chemical indicators, the pH range, the pK_a and the expected color (acid and base form). It is generally recommended to select a chemical indicator that has a pK_a as close to the endpoint of the titration as possible.

When chemical indicators are not suitable, a potentiometric pH titration can also be used. The pH of the solution is plotted versus the volume of titrant added.

Table 1: Aqueous acid-base chemical indicators

pH Range	Indicator	pK_a	Acid Form	Base Form
0.0 - 1.6	Methyl Violet		Yellow	Blue
1.2 - 2.8	Thymol Blue	1.65	Red	Yellow
3.2 - 4.4	Methyl Orange	3.46	Red	Yellow
3.8 - 5.4	Bromocresol Green	4.90	Yellow	Blue
4.8 - 6.0	Methyl Red	5.00	Red	Yellow
5.2 - 6.8	Chlorophenol Blue	6.25	Yellow	Red
6.0 - 7.6	Bromothymol Blue	7.30	Yellow	Blue
6.6 - 8.0	Phenol Red	8.00	Yellow	Red
7.4 - 9.0	Metacresol Purple	8.30	Yellow	Purple
8.0 - 9.6	Thymol Blue	9.20	Yellow	Blue
8.2 - 10.0	Phenolphthalein	9.50	Clear	Pink
9.4 - 10.6	Thymolphthalein		Clear	Blue
10.1 - 12.0	Alizarin Yellow R		Yellow	Red
11.4 - 12.6	Indigo Carmine		Blue	Yellow

Figure 4, Acid-base titration, shows a traditional strong acid-strong base titration curve, the volume of sodium hydroxide (NaOH) added to the solution is plotted against the pH of the solution. Note the abrupt change in the pH at the equivalence point.




Figure 4: Acid-base titration

4.2.2.2. Argentometric Titrations

Argentometric titrations use silver (nitrate) as the titrant and are generally precipitation titrations, as many silver salts are insoluble. These titrations are commonly used to titrate and determine the concentration of bromide, chloride, cyanide, iodide and sulfide. Argentometric titrations can be done with Mohr's indicator. After all of the chloride has reacted, a red silver chromate precipitate is formed or the titration can be easily followed with a silver ISE (or chloride ISE for chloride titrations) and a reference electrode.

Figure 5, Argentometric titration, shows the titration of a sodium chloride solution with silver nitrate (AgNO_3). The volume of AgNO_3 is plotted against the potentiometric signal from a chloride ISE.




Figure 5: Argentometric titration

4.2.2.3. Complexometric Titrations

A complex is a species where a central metal ion is covalently bonded to one or more electron donating groups called ligands. In a complexometric titration, metal ions are titrated using a titrant that binds strongly to it. Often these titrants contain EDTA or CDTA, polydentate ligands that form very stable coordination compounds with metal ions. The complexation reaction must be fast in order to be useful for direct titration. Some metal ions react too slowly with EDTA for a direct titration.

An indicator electrode that responds to the metal ion can be used to monitor the titration progress. The titration curve will appear similar to a potentiometric titration. Complexation indicators change color at the endpoint as all metal ions are "consumed" or complexed by the titrant.

Figure 6, Complexometric titration, shows a typically complexometric titration curve when using an indicator electrode that responds to the metal ion.




Figure 6: Complexometric titration

4.2.2.4. Ion-Selective Titrations

The most popular ion-selective titration is an acid-base titration. The hydrogen ion concentration is specifically measured and monitored during the titration process to locate the equivalence point. Using an ion-selective electrode (ISE) as the indicator electrode, the potentiometric signal (in mV) is used to directly follow a specific ion's concentration (or activity).

Examples of ISE titrations include titrating fluoride with an aluminum titrant using a fluoride ISE, chloride with silver nitrate using a chloride ISE, sodium with a sodium ISE, etc. The equivalence point can be determined by plotting the mV value vs. the amount of titrant added.

4.2.2.5. Non-Aqueous Solvent Acid-Base Titrations

Non-aqueous solvents must be used to titrate very weak acids and bases due to the inherent leveling effect water has on all acids and bases dissolved in it. A wide variety of weak acids and bases can be titrated using non-aqueous solvents. Mixtures of acids or bases can often be individually analyzed in a single sequential titration.

Titration of Acids

Weak acids with pK_a 's up to about 11 can be titrated in non-aqueous solvents. These include carboxylic acids, enols, phenols, imides, sulfonic acids and inorganic acids. Water or lower alcohols are suitable for titrating medium to strong acids (pK_a less than 5). Titrating a weaker acid with a strong base titrant requires a solvent less acidic than water or ethanol/methanol. Solvents such as acetone, acetonitrile, *t*-butyl alcohol, dimethylformamide, isopropanol and pyridine have been found to work well for acid-base titrations of strong, medium and weak acids/bases. Titrants include alcoholic potassium hydroxide and various sodium or potassium alkoxides in a 10:1 mixture of benzene/methanol. The best titrants are quaternary ammonium hydroxides (such as tetrabutylammonium hydroxide) due to good solubility of tetraalkylammonium salts of the titrated acids and the clean potentiometric titration curve obtained. **Figure 7**, Non-aqueous titration, shows an example of titration with tributylmethylammonium hydroxide titrant.




Figure 7: Non-aqueous titration

Titration of Bases

Weak bases with pK_b 's up to about 11, which do not ionize with water, can be titrated in non-aqueous solvents. These bases include aliphatic and aromatic amines, basic nitrogen heterocycles, alkali metal and amine salts of acids, and many other organic basic compounds. Titrating a weak base with a strong acid titrant requires a basic solvent that is as weak as possible. Water and alcohols allow the titration of medium strength bases, such as aliphatic amines ($pK_b = 4$ to 5), but not the titration of weaker bases such as pyridine ($pK_b = 8.8$). Glacial acetic acid works well for weak bases and has been used extensively. Less basic solvents such as acetone, acetonitrile and nitromethane extend the range of titratable compounds.

The endpoint for non-aqueous titrations are usually determined potentiometrically using a pH glass electrode, a modified calomel or double junction reference electrode with a low-flow rate reference junction. Good potentiometric titration curves are obtained in most solvents, except those with very low dielectric constants such as benzene or chloroform, when high electrical resistance of the solvent causes unstable potentials.

4.2.2.6. Precipitation Titrations

Precipitation titrations allow for faster analysis when compared to gravimetric analysis, where a precipitate is formed, filtered, dried and weighed to analyze a compound. Typically silver halides, silver thiocyanate and a few mercury, lead, and zinc salts are titrated using this method. The chemical reactions must form an insoluble salt and precipitate out quickly in order to be analyzed by this method. When the reaction is not quick, a back titration can be used. A measured excess of the precipitating reagent (titrant) is added to force the reaction to occur, and then unreacted titrant is titrated with a standard solution of another reagent.

4.2.2.7. Redox Titrations

There are a number of oxidation-reduction reactions that can be used to determine unknown concentration by titration. If the reaction goes to completion, is fast and has an analytical signal available to follow it, a titration can be performed. The term "fast" means that each addition of titrant is reacted completely and the sensing electrode is able to detect the change in solution in less than one second. **Figure 8**, Redox titration, shows an example of a redox titration using Cerium (IV) as a titrant.

Redox titrations are potentiometric titrations where the mV signal from a combination ORP (redox) electrode (usually with a platinum indicator electrode) is used to follow the reaction of oxidant/reductant. The electrode potential is determined by the Nernst equation and is controlled by the oxidant/reductant ratio.




Figure 8: Redox titration

Visual indicators, such as Ferroin, are also available. The oxidized and reduced form of the indicator will have different colors and can be used to determine the endpoint.

Various reductants can be determined by titrants with oxidants such as potassium permanganate, potassium chromate or iodine. Commonly used reductants that are used as titrants include sodium thiosulfate and ferrous ammonium sulfate.

As with acid-base titrations, the potential changes dramatically at the equivalence point.

4.2.2.8. Karl Fischer Titrations

This method is based on a well-defined chemical reaction between water and the Karl Fischer reagent. The chemistry provides excellent specificity for water determination. The method can be used to determine free and bound water in a sample matrix. The Karl Fischer method is widely considered to produce the most rapid, accurate and reproducible results and has the largest detectable concentration range spanning 1 ppm to 100%.

The determination of water content is one of the most commonly practiced methods in laboratories around the world. Knowledge of water content is critical to understanding chemical and physical properties of materials and ascertaining product quality. Water content determination is conducted on many sample types including pharmaceuticals and cosmetics, foods and natural products, organic and inorganic compounds, chemicals, solvents and gases, petroleum and plastic products as well as paints and adhesives. The KF method is verifiable and can be fully documented. As a result, Karl Fischer titration is the standard

method for analysis of water in a multitude of samples as specified by numerous organizations including the Association of Official Analytical Chemists, the United States and European Pharmacopoeia, ASTM, American Petroleum Institute, British Standards and DIN.

4.2.3. TITRATIONS ACCORDING TO THE TITRATION SEQUENCE

4.2.3.1. Back Titrations

Back titrations are generally used when a reaction is too slow to be directly accomplished during a “direct” titration, where the reaction goes to completion within a few seconds. In a back titration, a large excess of a reagent is added to the sample solution, helping a slow reaction to go to completion. The unreacted, excess reagent is then titrated. The difference in the total volume of the first reagent added and amount determined from the second titration is the quantity of reagent required to complete the first reaction.

4.2.3.2. Multiple Endpoint Titrations

Under certain conditions, some titrations can exhibit more than one equivalence point and be titratable to the individual endpoints to determine the concentration of each individual component. Examples of these types of titrations include acid-base (different strengths acids or bases are in a mixture), redox (each species has a different reduction potential), complexometric (different species are separately titratable), and acid-base, using polyprotic acids (the pK_a of the different protons varies enough to separate them). In **Figure 9A**, Multiple endpoint titrations, a titration of a polyprotic acid is shown, the different acid strengths of the first and second proton can be determined. **Figure 9B**, Multiple endpoint titrations, shows a titration with two different metal redox species, the different redox potentials allow the species to be separated. In **Figure 9C**, Multiple endpoint titrations, the solution being titrated contains a mixture of strong, weak, and very weak acids, the different pK_a 's allow the species to be separated.




Figure 9: Multiple endpoint titrations

4.3. TITRATION PROCEDURE

4.3.1. MANUAL TITRATION

Apparatus required for manual titration include:

- Volumetric burette, for precisely controlled delivery of titrant to the reaction vessel
- Erlenmeyer, or similar flask, that facilitates constant mixing or swirling required to ensure solution homogeneity
- Volumetric pipettes for the precise addition of samples and indicator solutions
- Titrant solutions of known concentration
- A visual or instrumental indicator for detecting the completion of the reaction

A typical manual titration consists of the following steps:

1. A volumetric pipette is used to add a known volume of sample to the flask.
2. An indicator solution or instrument probe is added to the flask.
3. A burette is used to measure the addition of titrant to the flask and dispense titrant in a controlled manner.
4. Titrant is added via the burette until the method indication signals the reaction endpoint.
5. Analyte concentration is calculated based on the concentration and volume of titrant required to reach the endpoint.



4.3.2. AUTOMATIC TITRATION

Automatic titrators are high-precision analytical instruments that deliver the titrant, monitor the physical change associated with the titration reaction, automatically stop at the endpoint and calculates the concentration of the analyte. Automatic titrators are best for repetitive titrations and high-accuracy analyses.

An automatic titrator must have an accurate liquid dispensing system. In high-accuracy systems like the [HI900-series](#) titrators, the liquid dispensing system consists of three main components: motor-driven syringe burette capable of accurately and precisely dispensing very small volumes of titrant, valve system capable of switching between titrant intake and outlet and a dispensing tip. These three main subsystem components must be as accurate as possible, with very low gear backlash in the burette pump, minimal piston seal flexing, precision ground inner diameter of the glass syringe, a low dead-volume valve, minimal evaporation/permeation, and chemically resistant tubing.

Apparatus required for automatic titration include:

- An automatic titrator, equipped with a burette
- A beaker
- An electronic stirring system, either a propeller stirrer or a magnetic stir bar and stir plate
- Volumetric pipettes for the precise addition of samples
- Standard titrant solutions of known concentration
- An electrode system that can be used to determine the endpoint of the titration

A typical automatic titration consists of the following steps:

1. Set up the automatic titrator according to the manufacturer's instructions.
2. Use a volumetric pipette to add a known volume of sample to the beaker.
3. Submerge the propeller stirrer or add the stir bar to the beaker and turn on.
4. Start the titration.
5. The titrator will automatically stop at the endpoint and determine the concentration of the analyte.

4.4. TITRATION RESULTS

4.4.1. ACCURACY

The factors most critical to achieving accurate results with the HI900 titration systems are the concentration of the sample, size of the sample and having an optimized set of method parameters.

4.4.2. REPEATABILITY

Repeatability or the agreement between replicate determinations, is expressed quantitatively as the relative standard deviation (RSD).

4.4.3. SOURCES OF ERROR

One of the advantages of volumetric analysis is excellent accuracy and precision. The sources of error can be grouped into sampling, titrant and standards, chemical reactions, endpoint determination and calculations.

4.4.3.1. Sampling Errors

- Selection of a non-homogeneous or non-representative sample
- Sample changed or was contaminated during collection, storage or transfers
- Poor technique when transferring sample to beaker or flask
- Errors in the balance (calibrate and check balance regularly)

4.4.3.2. Preparation Errors

Incorrect preparation due to:

- Poor technique in weighing the salt or when transferring to volumetric glassware
- Low-purity salts or water used to make titrant and standard
- Dirty or wet glassware
- Improper storage of titrant or standard which allows water gain, evaporation or deterioration
- Failure to standardize frequently, to adjust for change in titrant
- Failure to flush titrator tubing with a volume of titrant before standardizing
- Volume errors from pipettes and volumetric flasks (grade A glassware is required)
- Balance errors when weighing out salts (calibrate and check balance regularly)

4.4.3.3. Dispensing Errors

Incorrect dispensing due to:

- Dead valve volume and leaking valve
- Inaccuracy in motor drive and gear lash or backlash
- Poor burette or piston seal
- Non-uniform diameter of burette glass cylinder
- Chemical incompatibility with tubing or bubble generation
- Density or temperature changes in titrant
- Inadequate volume to cover electrode

4.4.3.4. Chemical Reaction Errors

- Inappropriate solvent or sample, resulting in side reactions
- Poor mixing in the titration vessel
- Reaction between titrant and sample is not rapid
- Reaction does not go to completion
- Reaction has side reactions

4.4.3.5. Endpoint Determination Errors

Most manual titrations use a visual indicator to indicate when the endpoint is reached and the titration should be stopped. Automatic titrators use instrumental methods to determine the end of a titration and the equivalence point. There are two predominant methods used to determine the equivalence point, first derivative and second derivative.

The first derivative is often used to determine the inflection point. The inflection point of the titration curve (mV vs. volume) is normally assumed to be the equivalence point. The maximum value of the first derivative (ΔmV vs. ΔV) corresponds to the theoretical equivalence point. During a titration it is rare to have a data point exactly at the first derivative maximum, the maximum value is determined by interpolating the first derivative data points.

The second derivative (ΔmV^2 vs. ΔV^2) can also be used to determine the equivalence point, and can offer advantages over the first derivative method. Second derivatives have increased sensitivity to smaller inflection points and easier numerical evaluation of the actual equivalence point. The value where the second derivative is equal to zero is the equivalence point. The second derivative requires fewer points located near the equivalence point, where data is often not obtained or not as reliable.

Errors in determining the endpoint can result from:

- Incorrect signals from the sensor
- Sensor drift
- Sensor or instrument has slow response (it is recommended to keep the sensors in good condition)
- Inappropriate setting on the titrator

4.5. CALCULATIONS

The main variables used in calculating a result from a titration are the sample volume, the concentration of the titrant, and the volume of titrant required to reach the equivalence point. At the equivalence point, an equal number of equivalents of the analyte and titrant has been added.

4.5.1. SAMPLE CALCULATION BY MASS

$$C_{\text{sample}} = \frac{V_{\text{titrant}} \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte}}}{m_{\text{sample}}} \times 100$$

C_{sample}	Sample Concentration (g/100g)
V_{titrant}	Volume of Titrant
C_{titrant}	Titrant Concentration (eq/L)
Ratio	Equivalence Ratio of analyte / titrant (mol analyte / eq titrant)
$\text{FW}_{\text{analyte}}$	Formula Weight of the Analyte (g/mol)
m_{sample}	Mass of Sample (g)

4.5.2. SAMPLE CALCULATION BY VOLUME

$$C_{\text{sample}} = \frac{V_{\text{titrant}} \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte}}}{V_{\text{sample}}} \times 100$$

C_{sample}	Sample Concentration (g/100mL)
V_{titrant}	Volume of Titrant
C_{titrant}	Titrant Concentration (eq/L)
Ratio	Equivalence Ratio of analyte / titrant (mol analyte / eq titrant)
$\text{FW}_{\text{analyte}}$	Formula Weight of the Analyte (g/mol)
V_{sample}	Volume of Sample (mL)

4.5.3. STANDARDIZE TITRANT BY MASS

Titrant standardization is the second most important calculation in titrations. A primary standard is titrated in order to determine the concentration of the titrant. This is essentially a typical titration calculated in “reverse”, where the concentration of the solution is known and the titrant is unknown.

$$C_{\text{titrant}} = \frac{m_{\text{standard}} \times \text{Ratio}}{FW_{\text{standard}} \times V_{\text{titrant}}}$$

C_{titrant}	Titrant Concentration (N)
m_{standard}	Mass of Standard (g)
Ratio	Equivalence Ratio of titrant / standard (eq titrant / mol standard)
FW_{standard}	Formula Weight of the Standard (g/mol)
V_{titrant}	Volume of Titrant (L)

4.5.4. STANDARDIZE TITRANT BY VOLUME

Titrant standardization is the second most important calculation in titrations. A primary standard is titrated in order to determine the concentration of the titrant. This is essentially a typical titration calculated in “reverse”, where the concentration of the solution is known and the titrant is unknown.

$$C_{\text{titrant}} = \frac{V_{\text{standard}} \times (1 \text{ L}/1000 \text{ mL}) \times C_{\text{standard}}}{V_{\text{titrant}}}$$

C_{titrant}	Titrant Concentration (N)
V_{standard}	Volume of Standard (mL)
C_{standard}	Concentration of Standard (eq/L)
V_{titrant}	Volume of Titrant (L)

4.5.5. BLANK TITRATION

In a blank titration a pre-titration is performed, often times on the solvent to be used for the sample titration, and the titrant volume required to reach the endpoint is noted. This blank value nullifies error due to titrant required to react with the components of the titration solution matrix. The basic titration equation can be used for a blank titration, with the single modification that the volume of titrant used in the blank titration should be subtracted from the regular titration titrant volume.

$$C_{\text{sample}} = \frac{C_{\text{titrant}} \times (V_{\text{sample}} - V_{\text{blank}}) \times \text{Ratio} \times FW_{\text{analyte}}}{m_{\text{sample}}} \times 100$$

C_{sample}	Sample Concentration (g/100 g)
C_{titrant}	Titrant Concentration (eq/L)
V_{sample}	Volume of Titrant required for the sample (L)
V_{blank}	Volume of Titrant required for the blank (L)
Ratio	Equivalence ratio of analyte / titrant (mol analyte / eq titrant)
FW_{analyte}	Formula Weight of the Analyte (g/mol)
m_{sample}	Mass of Sample (g)

4.5.6. MULTIPLE ENDPOINT TITRATION

Some titrations have two or more endpoints, each corresponding to the equivalence point for a specific reaction. Multiple endpoint titrations are similar to a blank titration in that the volume of titrant required to reach the first endpoint is subtracted from the titrant volume used to reach the next sequential endpoint.

$$C_{\text{sample1}} = \frac{V_{\text{titrant1}} \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte1}}}{m_{\text{sample}}} \times 100$$

$$C_{\text{sample2}} = \frac{(V_{\text{titrant2}} - V_{\text{titrant1}}) \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte2}}}{m_{\text{sample}}} \times 100$$

$$C_{\text{sample3}} = \frac{(V_{\text{titrant3}} - V_{\text{titrant2}}) \times C_{\text{titrant}} \times \text{Ratio} \times \text{FW}_{\text{analyte3}}}{m_{\text{sample}}} \times 100$$

C_{sample1}	Sample 1 Concentration (g/100g)
C_{sample2}	Sample 2 Concentration (g/100g)
C_{sample3}	Sample 3 Concentration (g/100g)
V_{titrant1}	Volume of titrant required to reach the first endpoint (L)
V_{titrant2}	Volume of titrant required to reach the second endpoint (L)
V_{titrant3}	Volume of titrant required to reach the third endpoint (L)
C_{titrant}	Concentration of Titrant (N)
Ratio	Equivalence Ratio of analyte / titrant (mol analyte / eq titrant)
$\text{FW}_{\text{analyte1}}$	Formula Weight of the analyte 1 (g/mol)
$\text{FW}_{\text{analyte2}}$	Formula Weight of the analyte 2 (g/mol)
$\text{FW}_{\text{analyte3}}$	Formula Weight of the analyte 3 (g/mol)
m_{sample}	Mass of Sample (g)

4.5.6.1. Back Titration

The equation used in back titration calculations is also similar to the equation for a blank titration. Instead of subtracting the initial amount of titrant needed to react with the blank, the amount of second titrant needed to react with the excess titrant added in the first titration is subtracted from the amount of the first titrant added. The difference between the two amounts is the amount of titrant necessary to reach the first equivalence point.

$$C_{\text{sample}} = \frac{(C_{\text{titrant1}} \times V_{\text{titrant1}} - C_{\text{titrant2}} \times V_{\text{titrant2}}) \times \text{Ratio} \times \text{FW}_{\text{analyte}}}{V_{\text{sample}}} \times 100$$

C_{sample}	Sample Concentration (g/100mL)
C_{titrant1}	Concentration of Titrant 1 (N)
V_{titrant1}	Volume of Titrant 1 (L)
C_{titrant2}	Concentration of Titrant 2 (N)
V_{titrant2}	Volume of Titrant 2 (L)
Ratio	Equivalence Ratio of analyte / titrant (mol analyte / eq titrant)
$\text{FW}_{\text{analyte}}$	Formula Weight of the analyte (g/mol)
V_{sample}	Volume of Sample (mL)

4.6. GLOSSARY

Acid	A chemical species that can donate one or more protons (hydrogen ions).
Acid-Base Titration	Stoichiometric neutralization titrations, based upon the reaction that occurs between an acid and a base.
Activity	A physical property corresponding to the concentration of all ions in a solution. Electrodes respond to activity.
Amperometric Titration	Titrations where the current flow between two electrodes (often a metal electrode and a reference electrode) are used to monitor the titration progress.
Analyte	The chemical species being measured in a titration.
Argentometric Titration	Titrations that use silver (nitrate) as the titrant. These titrations are typically precipitation titrations.
Automatic Titrator	An instrument designed to automatically carry out a titration. It will add the appropriate amount of titrant, determine the endpoint and calculate the results.
Back Titration	A type of titration where an excess amount of titrant is added to a sample, forcing a sluggish reaction to go to completion. The excess reagent is then "back" titrated with a second titrant.
Base	A chemical species that can accept one or more protons (hydrogen ions).
Biamperometric Indication	Uses a dual platinum pin electrode to measure the current flow through a titration solution.
Bivoltametric Indication	Uses a dual platinum pin electrode to measure the voltage required to maintain a constant current flow through a titration solution while constant voltage is applied across the platinum elements of the electrode.
Burette	A graduated cylindrical piece of laboratory glassware that is used to dispense precise amounts of solution.
Complex Ion	A species where a central metal ion is covalently bonded to one or more electron donating groups called ligands.
Complexometric Titrations	Metal ions are titrated using a titrant that binds strongly to it. The titrants often contain Ethylenediaminetetraacetic Acid (EDTA) or Cyclohexylenedinitrilotetraacetic Acid (CDTA).
Endpoint	The point where a titration is stopped because a physical change in the solution has indicated a completed titration. Titration endpoints typically coincide with the equivalence point. A fixed value endpoint (pH or mV) can be used as well. The titration will stop at the desired point regardless if the titration is complete.
Equivalence Point	The point where the quantity of titrant is stoichiometrically equal to the quantity of analyte.
Formal	The theoretical number of equivalents per liter of the solution. It is used in solutions where the exact concentration of a species may be affected by the other ions present, therefore the stated concentration may not be exactly correct.
Gravimetric Analysis	A quantitative determination of an analyte based on the mass of the solid.

Indicator Electrode	An electrode that responds to the species of interest. The electrode potential is proportional to the concentration or activity of that ion in the solution being measured.
Indicators	Chemical indicators are typically organic dyes that change form under different physical conditions, causing a color change that can be seen by an analyst. Typically used in manual titrations, chemical indicators have been replaced with electrometric indicators, which are used with automatic titrators.
Inflection Point	The point on a titration curve where the second derivative curve changes signs.
Ion-Selective Electrode (ISE)	An electrode that responds to a specific ion. The electrode potential is proportional to the concentration or activity of that ion in the solution being measured.
Karl Fischer Titration	A titration that uses a chemical reaction that is specific for determining water.
Manual Titration	A titration that is carried out by hand. The analyst must add the appropriate amount of titrant, determine the endpoint and calculate the results.
Molar	The concentration of a solute in a solution.
Mole (mol)	A quantity of a chemical species. The molecular weight of a substance in grams is equal to the mass of one mole of the substance. One mole is equal to 6.022×10^{23} atoms or molecules.
Monochromator	A device that allows only a narrow range of wavelengths to pass through it by separating the light into different wavelengths.
Multiple Endpoint Titration	A titration that reacts multiple species in solution, sequentially, using the same titrant. The concentration of each analyte can be determined from their respective endpoints.
Nernst Equation	The fundamental equation relating cell voltage to the concentration of a solution.
Neutralization	A chemical reaction where an acid and a base react to form a neutral salt and water.
Non-aqueous	A solution that does not contain water.
Non-aqueous Titration	A titration that is performed in non-aqueous solutions, typically used to titrate very weak acids and bases to eliminate the leveling effect water has on all acids and bases dissolved in it.
Normal	The concentration of a solution which accounts for any stoichiometric difference between the various species in a solution.
Oxidation / Reduction Potential (ORP)	The measurement describing whether a species wants to donate or accept electrons from other species in a redox reaction. If a solution's reduction potential is higher than the species it is reacting with, it will typically gain electrons or be reduced. If the potential is lower than the species it is reacting with, it will typically lose electrons or be oxidized.
Oxidant	The species that is accepting electrons in a redox reaction.
Pipette	Scientific apparatus that is used to deliver precise volumes of liquids.
Polyprotic Acid	Acids that are capable of donating more than one proton per acid molecule.
Potentiometric Titration	A titration in which the endpoint is determined by monitoring the voltage of the solution using an electrode.

Precipitation Titration	A titration in which the analyte reacts with the titrant to form an insoluble compound. The endpoint is typically detected with an ISE sensitive to either the analyte or titrant.
Reagent	The chemical added in a titration that causes the given reaction to occur.
Reduction-Oxidation Reaction (redox)	A chemical reaction in which the atoms involved in the reaction have their oxidation numbers changed. Reduction is the gain of electrons, which decreases the oxidation number. Oxidation is the loss of electrons, which increases the oxidation number.
Reductants	The electron donor in a redox reaction.
Reference Electrode	An electrode that supplies a constant electrode potential. It is used in combination with an "indicator" electrode, allowing for the "indicator" electrode potential to be measured.
Relative Standard Deviation (RSD)	A measure of the amount of relative variation in a set of data. It is calculated by dividing the standard deviation by the mean: $RSD = (\text{Standard Deviation of } X) * 100 / (\text{Mean of } X)$
Repeatability	The variation in sample measurements taken by a single person or instrument under the same conditions.
Spectrophotometric Titration	A titration in which the endpoint is marked by a change in color and/or color intensity.
Stoichiometry	The quantitative relationship of the reactants and products in a chemical reaction.
Titant	The chemical added in a titration that causes the given reaction to occur.
Titration	A quantitative, volumetric procedure used in analytical chemistry to determine the concentration of an analyte in solution. The concentration of the analyte is determined by slowly adding a titrant to the solution. As the titrant is added, a chemical reaction between the titrant and the analyte occurs.
Titration Curve	A graph containing the physical data obtained for a titration. The data plotted is often an independent variable (volume of titrant) vs. a dependent variable (pH of the solution). From the titration curve, the equivalence point or endpoint can be determined.

4.7. LIST OF FIGURES

Figure 1: Amperometric titrations	4-2
Figure 2: Potentiometric titrations	4-3
Figure 3: Spectrophotometric titrations	4-3
Figure 4: Acid-base titration	4-4
Figure 5: Argentometric titration	4-5
Figure 6: Complexometric titration.....	4-5
Figure 7: Non-aqueous titration	4-6
Figure 8: Redox titration.....	4-7
Figure 9: Multiple endpoint titrations	4-8

CERTIFICATION

All Hanna[®] instruments conform to the **CE European Directives**.



RoHS
compliant



Disposal of Electrical & Electronic Equipment. The product should not be treated as household waste. Instead, hand it over to the appropriate collection point for the recycling of electrical and electronic equipment, which will conserve natural resources. Ensuring proper product disposal prevents potential negative consequences for the environment and human health. For more information, contact your city, your local household waste disposal service, or the place of purchase.

RECOMMENDATIONS FOR USERS

Before using this product, make sure it is entirely suitable for your specific application and for the environment in which it is used. Any variation introduced by the user to the supplied equipment may degrade the meter's performance. For your and the meter's safety do not use or store the meter in hazardous environments.

WARRANTY

The **HI931** is warranted for two years against defects in workmanship and materials when used for its intended purpose and maintained according to instructions. This warranty is limited to repair or replacement free of charge. Damage due to accidents, misuse, tampering, or lack of prescribed maintenance is not covered.

If service is required, contact your local Hanna Instruments[®] office. If under warranty, report the model number, date of purchase, serial number (engraved on the bottom of the meter), and the nature of the problem. If the repair is not covered by the warranty, you will be notified of the charges incurred. If the instrument is to be returned to Hanna Instruments, first obtain a Returned Goods Authorization (RGA) number from the Technical Service department and then send it with shipping costs prepaid. When shipping any instrument, make sure it is properly packed for complete protection.