**Compact Dry ETB medium for enterobacteriaceae / Enterobacteriaceae / milieu pour nombre enterobacteriaca selectifs / medio para enterobacterias / medio per conta enterobacterias / meio para contagem enterobacterias**

100 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9431

40 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9432

**English**

**Deutsch**

**Français**

*Compact Dry ETB is a ready to use, selective plate contained glucose for the enumeration of Enterobacteriaceae*

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**Specimen pretreatment**

**Viable count in water or liquid foodstuff**

Drop 1 ml of specimen (dilute if necessary) on the middle of the Compact Dry plate.

**Viable count in solid foodstuff**

Add buffer solution (e.g. 0.9% saline solution, PBS, buffered peptone water, etc.) to the specimen and homogenize by stomacher®. Drop 1 ml of specimen (dilute if necessary) on the middle of the dry sheet of the Compact Dry plate.

**Viable count in swab test specimen**

Use a steril wet swab to wipe the surface which should be analyzed and put it back into the device with wiping solution. After shaking the device (approx. 10 seconds) drop 1 ml of wiping solution (dilute if necessary) on the middle of the Compact Dry plate. Use the Compact Dry Swab (Art. No. ZCS1002953, content 40 pcs) as single application (1 ml) or Promedia ST-25 (Art. No. Z0302, content 10 pcs) as multiple application (10 ml) for this procedure.

**Test instructions**

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator.
5. After incubation count the number of red / purple colored colonies underneath the plate. White paper placed under the plate helps to count the colonies.

**Incubation time** 24 ± 2 hours  
**Incubation temperature** 37 ± 1 °C

**Interpretation of the results**

Colonies grown are almost all red / purple. Red / purple colored colonies are *Enterobacteriaceae*. Red / purple and otherwise colored colonies together are the total gram negative bacteria.

**Storage and shelf life**

Keep at room temperature (+ 1 to + 30 °C).  
Total shelf life 16 months after manufacturing.

**Notes**

- Some colonies might not be clearly red / purple colored.
- High concentrations (> 250 cfu) on plates will cause the entire growth area to become red / pink. In this case dilute the specimen.
- After use please follow the current disposal regulations.
- The growth area is 20 cm<sup>2</sup>. The back of the plate has a grid carved of 1 x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies (> 300 cfu) grown, the bacteria count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids. However, dilutions are recommended.
- Compact Dry plates are produced at an ISO 9001 / ISO 13485:2003 certified site.

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