Compact Dry X-BC medium for *Bacillus cereus*

100 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9721

40 plates/Platten/plaques/placas/lastre/placas

Art-No. HS9722

English**Deutsch****Français***Compact Dry X-BC is a ready to use, chromogenic plate for presumptive detection of Bacillus cereus**Compact Dry X-BC ist eine gebrauchsfertige, chromogene Platte zum Nachweis von Bacillus cereus**Compact Dry X-BC est une plaque chromogène prête à l'utilisation pour détecter Bacillus cereus***Sample pretreatment****Viable count in water or liquid foodstuff**

Drop 1 ml of specimen (dilute if necessary) on the middle of the Compact Dry plate.

Viable count in solid foodstuff

Add buffer solution to the sample and homogenize by stomacher®. Drop 1 ml of specimen (dilute if necessary) in the middle of the dry sheet of the Compact Dry plate.

Viable count in swab test specimen

Use a steril wet swab to wipe the surface which should be analyzed and put it back into the device with wiping solution. Drop 1 ml of wiping solution (dilute if necessary) on the middle of the Compact Dry plate. Use the Compact Dry Swab (Art. No. ZCS1002953, content 40 pcs) as single application (1 ml) or Promedia ST-25 (Art. No. Z0302, content 10 pcs) as multiple application (10 ml) for this procedure.

Test instructions

1. Open the cap and drop 1 ml of specimen in the middle of the Compact Dry plate.
2. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
3. Put the cap again on the plate and write the information needed on the memorandum section.
4. Turn over the capped plate and put in the incubator.
5. After incubation mark the number of light blue / blue colonies. White paper placed under the plate helps to count the colonies.

Incubation time 24 ± 2 hours
plus a further 24 ± 2 hours

Incubation temperature 30 ± 1 °C

Interpretation of the results

Bacillus cereus forms light blue / blue colonies of by chromogenes contained in a medium. Though some bacteria other than *Bacillus cereus* may also grow and form white colonies in this plate, only light blue/blue colonies should be counted. *Bacillus thuringiensis* which indicates the biochemical properties of related species of *Bacillus cereus* may grow and form light blue / blue colonies.

Storage and shelf life

Keep at room temperature (+ 1 to + 30 °C).
Total shelf life 18 months after manufacturing.

Notes

- High cell concentrations on plates (> 300 cfu) will cause the entire growth area to become bluish / blue green. In this case dilute the sample.
- After use please follow the current disposal regulations.
- The growth area is 20 cm². The back of the plate has a grid carved of 1 cm x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates are produced at an ISO 9001 certified site.

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