

PREG-TONE II[®]
plus
“Ultrasound for Ultra Profits”[®]

fast

Testing can be completed in
less than 30 seconds.

safe

All testing done externally —
no harm to sow or litter.

accurate

Accuracy proven by many
thousands of producers.

DISTRIBUTED BY:

Renco Corporation
116 Third Ave. North
Minneapolis, Minnesota 55401, USA

Tel:

612-338-6124

Toll Free:

1-800-359-8181

Fax:

612-333-9026

Website:

www.renecorp.com

Email:

sales@renecorp.com



FORM 1000

©2000 RENCO CORPORATION. ALL RIGHTS RESERVED.

“Ultrasound for Ultra Profits”[®]



PREG-TONE II[®]
plus

ANIMAL PREGNANCY TESTER



DETECTS PREGNANCY IN SWINE

✓ fast ✓ accurate ✓ efficient

PREG-TONE II[®] plus

AVAILABLE WITH OPTIONAL EARPHONE

Easy one-hand operation



The Preg-Tone II Plus is very simple to use.

- Step one:** Apply contact fluid, such as mineral or cooking oil, to animal in area shown.
- Step two:** Turn the instrument on by pressing the patented case with your thumb.
- Step three:** Press instrument to the skin and scan interior of the animal.

DURABLE

- *High-impact plastic*
- *Acid resistant*
- *Will last for years*

INEXPENSIVE

- *9 volt replaceable battery*
- *Battery will last for months*

EASY TO USE

- *One hand operation*
- *Use any oil for skin contact*
- *No lights to watch*

INDICATORS

- *Loud Sounder*
- *Intermittent tone indicates good skin contact*
- *Continuous tone indicates pregnancy*

OPTIONAL

- *Earphone attachment available for extra noisy conditions*

QUALITY

- *Renco has over 30 years experience testing and manufacturing ultrasound equipment.*

