


# NANOCOLOR® gesamt-Phosphat 5

## total Phosphate / Phosphate total / Fosfato total













**Test 0-81**  
REF 985081

Methode / Method / Méthode / Método

690 / 710 nm


0811 0.20–5.00 mg/L P  
0812 0.5–15.0 mg/L PO<sub>4</sub><sup>3-</sup>

Min	20–80 %		Max
0.20	1.16	4.04	5.00 mg/L P
0.5	3.4	12.1	15.0 mg/L PO <sub>4</sub> <sup>3-</sup>



# NANOCOLOR® ortho-Phosphat 5

ortho-Phosphate / ortho Phosphates / orto Fosfato



**Test 0-81**  
REF 985081

Methode / Method / Méthode / Método

690 / 710 nm

0815 0.20–5.00 mg/L PO<sub>4</sub>-P

0816 0.5–15.0 mg/L PO<sub>4</sub><sup>3-</sup>

Min	20–80 %		Max
0.20	1.16	4.04	5.00 mg/L PO <sub>4</sub> -P
0.5	3.4	12.1	15.0 mg/L PO <sub>4</sub> <sup>3-</sup>

