

visicolor[®] ECO Eisen 2


Iron / Fer / Hierro

530 nm

Method(e) / Método

5261


0.04–2.00 mg/L Fe




Probe / Sample
Echantillon
Muestra



1
Küvette spülen
Rinse test tube
Rincer l'éprouvette
Lavar el tubo




2
5 mL Probe
Sample
Echantillon
Muestra



3
4 ⚡ Fe-1



4
Schütteln
Shake
Agiter
Agitar




5
1 Fe-2



6
Schütteln
Shake
Agiter
Agitar



7
Säubern
Clean
Nettoyer
Limpiar



8
7'00 min



9
Messung
Measurement
Mesure
Medición

Messtoleranzen:
Measuring tolerances:
Tolérance de mesure :
Tolerancia de la medición:
± 0.04 mg/L Fe

Meerwasser
Sea water
Eau de mer
Agua de mar

