


# NANOCOLOR<sup>®</sup> gesamt-Phosphat 45

## total Phosphate/Phosphate total/Fosfato total













**Test 0-55**  
**REF 985055**

Methode / Method / Méthode / Método

690 / 710 nm


0551 5.0–50.0 mg/L P  
0552 15–150 mg/L PO<sub>4</sub><sup>3-</sup>  
0553 10–115 mg/L P<sub>2</sub>O<sub>5</sub>

Min	20–80 %		Max
5.0	14.0	41.0	50.0 mg/L P
15	42	123	150 mg/L PO <sub>4</sub> <sup>3-</sup>



# NANOCOLOR<sup>®</sup> ortho-Phosphat 45

## ortho-Phosphate / ortho Phosphates / orto Fosfato



**Test 0-55**  
REF 985055

Methode / Method / Méthode / Método

0554 5.0–50.0 mg/L PO<sub>4</sub>-P  
0555 15–150 mg/L PO<sub>4</sub><sup>3-</sup>  
0556 10–115 mg/L P<sub>2</sub>O<sub>5</sub>

690 / 710 nm

Min	20–80 %		Max
5.0	14.0	41.0	50.0 mg/L PO <sub>4</sub> -P
15	42	123	150 mg/L PO <sub>4</sub> <sup>3-</sup>

