

Medi-Test ordering for urine information

Diagnostic test strips for urine analysis

REF	REF	Type	with test fields for determination of												
50 strips per pack	100 strips per pack		Glucose	Ketones	Ascorbic acid	Protein/Albumin	Blood	Leukocytes	Nitrite	Density	pH value	Bilirubin	Urobilinogen	Creatinine	
93001	93024	Glucose ¹⁾	●												
93025	93020	Glucose / Keton ¹⁾	●	●											
93004	93027	Protein 2 ¹⁾				●					●				
93005	93028	Keton ¹⁾		●											
93006	93029	Nitrit ¹⁾							●						
93012	–	Urbi										●	●		
93015	93037	Combi 2 ¹⁾	●			●									
93007	93030	Combi 3 A [®]	●		●	●					●				
93009	93032	Combi 5	●		●	●	●				●				
93035	93036	Combi 5 N [®]	●		●	●	●		●		●				
93055	–	Combi 5 S	●	●		●	●				●				
93018	93078	Combi 6	●			●	●	●	●				●		
93013	93034	Combi 6 A	●	●	●	●	●				●	●			
93010	93022	Combi 7	●	●	●	●	●		●		●				
–	93021	Combi 8 L	●		●	●	●	●	●	●	●				
93011	93023	Combi 9 [®]	●	●	●	●	●		●		●	●	●		
–	93056	Combi 10 [®]	●	●	●	●	●		●	●	●	●	●		
93079	93058	Combi 10 [®] L	●	●	●	●	●	●	●		●	●	●		
93077	93067	Combi 10 [®] SGL	●	●		●	●	●	●	●	●	●	●		
–	93068	URYXXON [®] Stick 10 ^{2) 3)}	●	●		●	●	●	●	●	●	●	●		
930874 ^{3) 5)}	–	Microalbumin				●								●	

Further Medi-Test urine test strips

Test strips for veterinary applications

–	930870	Combi 10 [®] VET ⁴⁾	●	●		●	●	●	●	●	●	●	●	●
---	--------	---	---	---	--	---	---	---	---	---	---	---	---	---

Test strips for detection of urine adulteration

93019	–	Medi-Test Adulteration Stick ⁴⁾ for the detection of urine adulteration prior to drugs-of-abuse tests with test fields for Creatinine, Glutaraldehyde, Nitrite, pH, Oxidants and Density
-------	---	---

All products (except 4) are CE-marked according to the directive 98/79/EC

¹⁾ suitable for patient self-testing ²⁾ for evaluation with reflectometer URYXXON[®] 500 ³⁾ for evaluation with reflectometer URYXXON[®] Relax ⁴⁾ not an IVD product (no CE-mark) ⁵⁾ pack of 24 test strips

Instruments for urine test strips

93088	URYXXON [®] Relax, automated urine analysis at the point of care
930080	URYXXON [®] 500, high performance automated urine analysis

Accessories and spare parts

93038	Medi-Test Control, solution to check Medi-Test urine test strips and URYXXON [®] instruments
93065	Printer paper for URYXXON [®] Relax, pack of 5 rolls
930890	Strip adaptor for URYXXON [®] Relax
93071	Printer paper for URYXXON [®] 300 and URYXXON [®] 500, pack of 5 rolls
93074	Barcode scanner for URYXXON [®] instruments

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

