


Refractometer MASTER-M



Can be widely used for beverages, fruit juices, soups, beef bowl sauce, water-soluble cutting oil, etc. MASTER-M is made of metal, and **MASTER-PM** is plastic. For salty and acidic samples, choose the plastic body, **MASTER-PM**. This product is the cheapest type with no temperature compensation function. Since the concentration changes as the measurement temperature changes, **MASTER-α** and **MASTER-T** are recommended if you have trouble managing the temperature at each measurement.

Specifications

Model	MASTER-M
Cat.No.	2313
Range	Brix : 0.0 to 33.0%
Minimum Scale	0.2%
Accuracy	±0.2%

Accessories

- SUPER EXTRACTOR : RE-29401
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- Daylight plate for small volume samplings : RE-2311-67M
- Daylight plate : RE-2315-60M
- * Calibration Certificate : Contact an ATAGO representative for details.