

Supplied complete with meter, probe, calibration solutions, and cleaning solutions.

Specifications HI981036 0.00 to 12.00 pH Range Resolution 0.01 pH рΗ Accuracy ±0.05 pH (@25°C/77°F) Calibration automatic, one or two-point Temperature automatic, 0 to 50°C (32 to 122°F) Compensation Glass Type LT (low temperature) CR2032 Li-ion / approx. 800 Battery Type / Life Additional hours of continuous use Specifications Auto-off 8 minutes, 60 minutes, or can be disabled Environment 0 to 50°C (32 to 122°F); RH 95% max Dimensions 51 x 148 x 21 mm (2 x 5.8 x 0.9") Weight 43 q (1.5 oz.) HI981036 is supplied with pH 4.01 buffer solution sachets (2), pH 7.01 buffer solution sachets (2), electrode acid cleaning Ordering solution sachets for meat, grease and fats (2), gelled bridge

Foodcare

HI981036

Meat pH Tester

with specialized probe

The HI981036 Foodcare Meat pH tester is an application specific designed pH tester for the measurement of pH during the meat processing process. This meter offers many advanced features including a pH electrode designed specifically for meat.

pH electrode with replaceable bridge electrolyte

 The pH electrode has an outer junction sleeve that can be removed and cleaned. Once cleaned a small amount of supplied gel electrolyte is added and the junction is refreshed improving the pH measurement and extending the life of the meter.

Low temperature (LT) glass

 The pH glass tip uses a special low temperature (LT) glass formulation with a lower resistance of approximately 50 Megaohms compared to general purpose (GP) with a resistance of about 100 Megaohms. This is beneficial when measuring food products at lower temperatures in order to have the ideal resistance for the measuring circuit.

Conical tip

 Allows for easy penetration into solids and semisolids, which is needed when wanting to take a direct measurement in meat.

PVDF body

 Polyvinylidene Fluoride (PVDF) is a food grade plastic that is resistant to most chemicals and solvents, including sodium hypochlorite. It has high abrasion resistance, mechanical strength, and resistance to ultraviolet light. PVDF is also resistant to fungal growth.

Large LCD

 An enhanced LCD displays the measurement reading, stability indicator, low battery indicator, and calibration tags.

Automatic calibration

 The Foodcare Meat pH Tester is calibrated automatically to one or two points. The calibration buffers are automatically recognized and after calibration the buffer values used are shown on the display.

• Automatic temperature compensation

Stability indicator

 An hourglass indicator is displayed on the LCD until a stable reading is obtained. Once a reading stabilizes, the indicator disappears and a reading can be taken.

Automatic shut-off

 The meter can be set to automatically turn off after 8 minutes or 60 minutes to conserve battery life in the event that the meter is left on. The auto-off feature can also be disabled.

Probe diagnostic

 During calibration the meter will display an error (Err) message as an indicator that the probe needs to be cleaned.

Long battery life

 The Foodcare Meat pH Tester has an exceptional battery life of approximately 700 hours. When the battery power is running low, a low battery indicator is displayed.

Economica

 The Foodcare Meat pH Tester is a fully featured meter at a price that is affordable for both the hobbyist and professional.



electrolyte, electrode storage solution, CR2032 3V Li-ion battery, quality certificate, and instruction manual.

Information