









## Compact baby scale with special Safe-Guard tray and approval for professional mobile or stationary medical use in medical diagnostics, verification optional

## **Features**

- Verification class III (verification is optional)
- Approved as a medical device according to 93/42/EEC or regulation (EU) 2017/745
- Modern design
- Safe-Guard tray: the special, particularly recessed tray with raised side bolsters offers optimum conditions for the safe weighing of babies
- Large, insensitive plastic tray, without corners and edges, therefore particularly hygienic and easy to clean
- Lightweight, compact design, therefore easy to transport and store
- Secure and non-slip positioning with heightadjustable rubber feet
- Easy to use
- Hold function: While weighing patients that are unable to lie still, a mean average weight value is determined
- Integrated AUTO-OFF function to preserve the hatteries

## Technical data

- · LCD display, digit height 25 mm
- Dimensions of tray W×D 600×260 mm
- Overall dimensions W×D×H 600×360×141 mm
- Batteries included, 6×1.5 V AA, operating time up to 20 h (up to 50 h without using WiFi)
- Net weight approx. 4,6 kg

## Accessories

- Mechanical height rod to fix to the scale, measuring range 30 cm - 80 cm, increments of 0,1 cm, KERN MBA-A01
- WiFi data interface to transfer data from balances to PC or tablet, must be ordered at purchase, KERN YMI-A01
- Alcohol-free cleaning cloths for disinfectant wiping, quick acting, based on modern quaternary ammonium compounds and effective against papova viruses. Particularly gentle on materials. Fulfill the legal requirements for occupational safety in accordance with TRGS 525/540. Packaging contents 100 pcs., KERN MYC-01
- External mains adapter, 100 V 240 V,
   Standard EU, CH, UK, US, AUS, KERN YKA-43

STANDARD OPTION FACTORY

KCP MOVE BATT 1 DAY 230 V WIFI +3 DAYS

**C€** 0297

\*Within the EU, official verification (conformity assessment in accordance with NAWI 2014/31/EU) is mandatory by law for scales that are intended for use as a medical device. Please order this where you state the location of installation and postcode.

